

**PURPOSE:** To provide a workflow and identification scheme for Yeasts from clinical specimens.

### **INTERPRETATION AND REPORTING RESULTS:**

Perform wet prep on suspect colonies. If yeast cells are seen, proceed as follows:

**1. Sterile sites, blood cultures, body fluid cultures, biopsy specimens, invasively collected urines (straight catheter, suprapubic aspirate, cystoscopy specimen), corneal specimens, immunocompromised patients, neonates:**

- Preliminary report as Yeast species.
- Perform Vitek 2 YST identification card.
- If yeast other than *Candida* spp. is identified, report as “Probable” and send to DynaLIFE for confirmation and susceptibility testing.
- If *Candida* spp. is identified, report with quantification and refer to DynaLIFE for susceptibility testing.

**2. Respiratory isolates including sputa and specimens obtained by bronchoscopy:**

**NOTE:** *Candida* organisms are not a cause of pneumonia and are most often contaminants of the procedure, except possibly in oncology or lung transplant patients or in neonates. Even in those cases, growth of *Candida* spp. in lower respiratory specimens, regardless of species, does not correlate with disease. Yeasts are normal inhabitants of the mouth.

- Perform Vitek 2 YST identification card.
- If *Cryptococcus* spp. is identified, report as “Probable” and refer to DynaLIFE for confirmation and susceptibility testing.
- If *Cryptococcus* spp. is not identified, report as part of normal flora without specifically mentioning the presence of yeast (report as usual oropharyngeal flora) unless present in pure growth. Refer to MIC32300 – Respiratory Culture.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

Document Name: Identification of Yeasts	Document Number: MIC40200	
	Version No: 1.0	Page: 2 of 3
	Effective: <b>DRAFT</b>	

### 3. Voided urines, superficial sites, wounds and drainage fluids:

- Perform Vitek 2 YST identification card.
- If *Candida* spp. is identified, report identification if indicated in procedure for specimen site.
- If *Candida* spp. is not identified, report as Yeast species if indicated in procedure for specimen site.

### 4. Isolates from any other sites:

- Perform Vitek 2 YST identification card.
- If *Candida albicans* is identified, report identification if indicated in procedure for specimen site.
- If *Candida albicans* is not identified, report as Yeast, not *Candida albicans* if indicated in procedure for specimen site.

### **PROCEDURE NOTES:**

- Vitek 2 YST identification card: If results are not satisfactory (<90% confidence), perform urease test to rule out *Cryptococcus* and refer isolate to DynaLIFE for identification if clinically significant.
- *Candida albicans*: budding yeast cells in smear and feet in <48 hours.
- *Candida dublinensis* is difficult to distinguish from *Candida albicans*. It is germ tube positive, but growth on Sabouraud agar at 42°C to 45°C at 48 hours is absent or poor, whereas *Candida albicans* grows at 42°C to 45°C in 48 hours. It is currently not recommended that laboratories routinely perform additional tests to differentiate these organisms. However, a history of clinical failure or persistence of *Candida albicans* despite therapy with fluconazole may alert the clinician that this may be *C. dublinensis*, as this organism can develop resistance to fluconazole during therapy.
- *Candida glabrata*: small yeasts in smear with no hyphae, better growth on Chocolate agar than Blood agar.
- Increasing resistance to fluconazole has been reported in *Candida* species.
- *Cryptococcus*: spherical pleomorphic budding yeast with no hyphae, typically mucoid due to presence of capsular material, becoming dryer and duller with age, urea positive. *Cryptococcus gattii* is difficult to distinguish from *C. neoformans*. *C. neoformans* affects mostly immunocompromised hosts (malignancy, HIV, etc.) but *Cryptococcus gattii* may cause disease in immunocompetent hosts.

<b>NOTE:</b> This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

**REFERENCES:**

- Clinical Microbiology Procedures Handbook, 4<sup>th</sup> edition, ASM Press, 2016
- Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. 2015. Manual of Clinical Microbiology, 11<sup>th</sup> edition, ASM Press, Washington, D.C.
- Vitek 2 Systems product information
- CLSI. *Abbreviated Identification of Bacteria and Yeast; Approved Guideline—Second Edition*. CLSI document M35-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2008

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0		Initial Release	L. Steven

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: