



Stanton Territorial Hospital

P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Number: MIC31200

Version No: 1.0

Page: 1 of 6

Distribution:
Microbiology Culture Manual

Effective:
Date Reviewed:
Next Review:

Document Name: GBS Screen

Approved By:

Status: **DRAFT**

PURPOSE: To screen for Group B *Streptococcus* (GBS) in vaginal/rectal specimens.

SAMPLE INFORMATION:

Type	Swab <ul style="list-style-type: none"> Amie’s with or without charcoal
Source	<ul style="list-style-type: none"> Combined introital (vaginal and anorectal area) swab. Vaginal swabs are not the specimen of choice but will be processed.
Stability	<p>If the sample is received in the laboratory and processed greater than 48 hours from collection:</p> <ul style="list-style-type: none"> Add specimen quality comment: “Delayed transport may adversely affect pathogen recovery”.
Storage Requirements	Room temperature
Criteria for rejection and follow up action	<ol style="list-style-type: none"> Unlabeled/mislabeled swab. Specimen container label does not match patient identification on requisition. Duplicate specimens obtained with same collection method within 24 hours. Dry swab.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

Document Name: GBS Screen	Document Number: MIC31200	
	Version No: 1.0	Page: 2 of 6
	Effective: DRAFT	

REAGENTS and/or MEDIA:

- LIM Broth (LIM), StrepB *Select* agar (GBS) and Blood agar (BA)
- Identification reagents: catalase and Strep latex test

SUPPLIES:

- Disposable inoculation needles
- Biosafety cabinet
- 35° CO₂ incubator
- 35° ambient air incubator
- Glass test tubes
- Wooden sticks
- Sterile pipettes

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:


- Refer to MIC60040 – Culture Media Quality Control procedure.
- Refer to Test Manual for reagent quality control procedures.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

PROCEDURE INSTRUCTIONS:

Step	Action
Processing swabs for GBS screening	
1	Place swab in LIM broth, break off leaving the swab in the broth medium. Loosely recap. Incubate broth at 35° in CO ₂ incubator.
2	At 14:00, after 18 to 24 hours incubation, subculture the broth to StrepB <i>Select</i> agar: <ul style="list-style-type: none"> • Remove the required number of StrepB <i>Select</i> plates from the refrigerator and bring to room temperature. • Saturate a sterile swab in the broth and rotate against the wall of the tube above the liquid to remove excess inoculum. • Swab the first quadrant of the agar. • Streak for isolated growth using a disposable inoculation needle. • Streak out to cover the whole plate.
3	Incubate plate in O ₂ incubator at 35° for 24 hours in the “GBS 24 hours” tray.

INTERPRETATION OF RESULTS:

Step	Action				
1	Remove culture plates after 24 hours incubation.				
2	Observe plates for blue colonies: <div style="text-align: center;">  </div>				
If blue colonies are not seen:					
3	<table border="1" style="width: 100%;"> <tr> <td style="width: 30%;">No blue colonies seen at 24 hours</td> <td> <ul style="list-style-type: none"> • Record observations in the LIS. • Re-incubate plates in O₂ incubator on the “Old urine culture” shelf. </td> </tr> <tr> <td>No blue colonies seen at 48 hours</td> <td> <ul style="list-style-type: none"> • Record observations in the LIS. • Workup complete. GBS not isolated. </td> </tr> </table>	No blue colonies seen at 24 hours	<ul style="list-style-type: none"> • Record observations in the LIS. • Re-incubate plates in O₂ incubator on the “Old urine culture” shelf. 	No blue colonies seen at 48 hours	<ul style="list-style-type: none"> • Record observations in the LIS. • Workup complete. GBS not isolated.
No blue colonies seen at 24 hours	<ul style="list-style-type: none"> • Record observations in the LIS. • Re-incubate plates in O₂ incubator on the “Old urine culture” shelf. 				
No blue colonies seen at 48 hours	<ul style="list-style-type: none"> • Record observations in the LIS. • Workup complete. GBS not isolated. 				

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

If blue colonies are seen:	
IF	THEN
Blue colonies are not isolated	<ul style="list-style-type: none"> Subculture blue colonies to blood agar plate. Perform catalase test from subculture plate.
Blue colonies are isolated	<ul style="list-style-type: none"> Perform catalase test.
IF	THEN
Catalase: POSITIVE	<ul style="list-style-type: none"> Record observations in LIS. Workup complete. GBS not isolated.
Catalase: NEGATIVE	<ul style="list-style-type: none"> Record observations in LIS. Perform Strep latex test for Group B.
IF	THEN
Strep B latex test: NEGATIVE	<ul style="list-style-type: none"> Record observations in LIS. Workup complete. GBS not isolated.
Strep B latex test: POSITIVE	<ul style="list-style-type: none"> Record observations in LIS. GBS isolated.
IF	THEN
Group B screen is positive and clinical history does not state penicillin allergy	<ul style="list-style-type: none"> Susceptibility testing not performed.
Group B screen is positive and clinical history states penicillin allergy	<ul style="list-style-type: none"> Perform susceptibility testing as per ASTM.

Please Note:

Each *Streptococcus* grouping latex test should be tested with at least one extra grouping latex suspension as a negative control.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

REPORTING RESULTS:

IF	REPORT
GBS not isolated	<ul style="list-style-type: none"> Report: “No Streptococcus agalactiae (Group B) isolated” A copy is automatically sent to the STH Obstetrics Unit.
GBS isolated and no penicillin allergy indicated in clinical history	<ul style="list-style-type: none"> Choose key 6 on STRB keypad to add isolate: “Streptococcus agalactiae (Group B)” List quantitation as “Isolated” The following canned culture comment will be added: “Current guidelines suggest intrapartum GBS prophylaxis should be with penicillin (or ampicillin) IV, cefazolin if nonsevere penicillin allergy, clindamycin in severe beta-lactam allergy, or vancomycin in severe beta-lactam allergy where clindamycin tests resistant or susceptibility is unknown. A severe allergy is defined as anaphylaxis, angioedema, respiratory distress or urticaria following administration of penicillin or a cephalosporin. If patient has severe allergy and susceptibility testing to clindamycin is needed for intrapartum antibiotic prophylaxis, please contact the laboratory IMMEDIATELY (within 5 days) for susceptibility testing at 867-669-4162.” A copy is automatically sent to the STH Obstetrics Unit.
GBS isolated and clinical history indicates penicillin allergy	<ul style="list-style-type: none"> Choose key 7 on the STRB keypad to add isolate: “Streptococcus agalactiae (Group B)” List quantitation as “Isolated” KB susceptibility panel for GBS screen is ordered. Report susceptibility results as per ASTM. A copy is automatically sent to the STH Obstetrics Unit.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

LIMITATIONS:

1. The selective mixture of antibiotics/antifungal agents inhibits the growth of the majority of microorganisms, with the exception of *Streptococci*, some *enterococci* and *lactobacilli*.
2. As with any chromogenic medium, it is important to streak at closely spaced intervals in order to obtain well-isolated colonies: the morphology and colour of the colonies will then be more typical.
3. If an inoculum contains a high density of *Streptococcus agalactiae*, the medium around the deposit may be coloured.
4. The colonies of some species, other than *Streptococcus agalactiae* (for example, *Streptococcus pyogenes*, *Streptococcus porcinus* and *Streptococcus gallolyticus*) may appear blue.
5. The intrinsic demands of some *Streptococci* can lead to absence or partial inhibition of their growth.

REFERENCES:

- Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016.
- Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. 2015. Manual of Clinical Microbiology, 11th edition, ASM Press, Washington, D.C.
- StrepB Select agar package insert. 08/2009.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0		Initial Release	L. Steven

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: