

Document Name: Water Testing – HPC Unit Dose SimPlate Method

Approved By:

Jennifer G. Daley Bernier, A/Manager, Laboratory Services

Status: **APPROVED**

PURPOSE: The presence of heterotrophic bacteria in dialysate waters can lead to the development of a Gram negative toxin mediated pyrogenic reaction, bacteremia, and chronic inflammatory response syndrome. Hot tub waters are screened for the presence of *Pseudomonas aeruginosa*, the causative agent of a superficial skin infection known as “Hot tub folliculitis”. The HPC SimPlate method utilizes enzyme technology to target the most common enzymes of waterborne bacteria. The by-product of the enzymatic reactions can be seen as fluorescence using a UV light. Additionally, hot tub waters have a MacConkey plate inoculated to aid in the detection of *Pseudomonas aeruginosa*.

SAMPLE INFORMATION:

Type	Hot tub water, Dialysate water
Volume	10mL +/- 0.2mL
Storage Requirements	Refrigerated
Criteria for rejection and follow up action	<ul style="list-style-type: none"> <10 mL water received Received >48 hours after collection

REAGENTS and/or MEDIA:

Type	SimPlate for HPC kit from IDEXX
Storage Requirements	<ul style="list-style-type: none"> Store at 2-30°C away from light Expiry date is printed on the box of media tubes

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

Document Name: Water Testing – HPC Unit Dose SimPlate Method	Document Number: MIC10315	
	Version No: 2.0	Page: 2 of 6
	Effective: 12 May, 2017	

SUPPLIES:

- 10mL syringe
- UV light
- 35°C incubator

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

Performed per shipment or lot number:

- Refer to QC Procedure MIC60500 – Water quality control.
- A TQC order is automatically generated when new kit is entered into TQC to record the QC results.


NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC10315HPCUnitDoseSimPlateMethodPRO	Print Date: 20/03/2019 3:02:00 PM

PROCEDURE INSTRUCTIONS:

Step	Action
Performing HPC Water Testing	
1	Accession waters and select source. If source is HOT TUB or SPA WATER – HPC will automatically order and a MacConkey plate label will generate. Labels for total coliform testing will also print out for these waters.
2	Label requisition and sample containers.
3	Label sterile media tube (green top), SimPlate and MAC plate (where required)
4	Mix water. Use a 10 mL syringe to remove 10 mL of water from sample container. Add to green top tube and shake. Allow powder in tube to dissolve
5	Remove SimPlate lid and pour contents of green top tube onto the center of the plate base
6	Replace the lid and gently swirl to distribute the sample. Note: air bubbles do not interfere with test
7	Tip the plate at a 90° angle so the excess water will drain into the absorbent pad at the bottom of the plate
8	Invert the plate onto the plastic lid. On lid write R: (Date + 2 day) and time incubated
9	If water source is HOT TUB, perform total coliform testing as well using the Colilert-18 kit and processed as above.
10	If water is HOT TUB, aliquot off an additional 1mL of water and flood a labeled MAC plate - incubate along with the SimPlate (35°C for 48hrs)

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

INTERPRETATION AND REPORTING:

1	While viewing under the UV lamp - count the number of wells showing fluorescence.
2	<p>Using the HPC keypad, select the number of wells that show fluorescence. This will automatically calculate the MPN (most probably number) of heterotrophic plate count bacteria in the original sample.</p>  <p>Report this number as the HPC.</p>
3	<p>HOT TUB/SPA WATERS:</p> <ul style="list-style-type: none"> • Look for growth of <i>Pseudomonas aeruginosa</i> – Non-lactose fermenter, oxidase positive – send to VITEK for GNI • Report as: “Isolated: <i>Pseudomonas aeruginosa</i>” • If MAC plate is negative- Report as: “No <i>Pseudomonas aeruginosa</i> isolated”

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

APPENDIX:

**Unit-Dose
SimPlate For HPC
Most Probable Number (MPN) Table**

# Positive Wells	MPN	95% confidence limits		# Positive Wells	MPN	95% confidence limits	
		lower	upper			lower	upper
0	<0.2	<0.03	<1.4	43	12.0	8.8	16.2
1	0.2	0.0	1.4	44	12.4	9.1	16.7
2	0.4	0.1	1.6	45	12.8	9.5	17.3
3	0.6	0.2	1.9	46	13.2	9.8	17.8
4	0.8	0.3	2.2	47	13.7	10.2	18.3
5	1.0	0.4	2.5	48	14.1	10.6	18.9
6	1.2	0.6	2.7	49	14.6	10.9	19.5
7	1.5	0.7	3.0	50	15.1	11.3	20.1
8	1.7	0.8	3.3	51	15.6	11.7	20.7
9	1.9	1.0	3.6	52	16.1	12.1	21.3
10	2.1	1.1	3.9	53	16.6	12.5	22.0
11	2.3	1.3	4.2	54	17.1	13.0	22.7
12	2.6	1.5	4.5	55	17.7	13.4	23.4
13	2.8	1.6	4.8	56	18.3	13.9	24.1
14	3.0	1.8	5.1	57	18.9	14.4	24.9
15	3.3	2.0	5.4	58	19.5	14.9	25.7
16	3.5	2.2	5.8	59	20.2	15.4	26.5
17	3.8	2.3	6.1	60	20.9	15.9	27.3
18	4.0	2.5	6.4	61	21.6	16.5	28.2
19	4.3	2.7	6.7	62	22.3	17.1	29.2
20	4.5	2.9	7.0	63	23.1	17.7	30.2
21	4.8	3.1	7.4	64	23.9	18.3	31.2
22	5.1	3.3	7.7	65	24.8	19.0	32.3
23	5.3	3.5	8.0	66	25.7	19.7	33.5
24	5.6	3.8	8.4	67	26.6	20.4	34.7
25	5.9	4.0	8.7	68	27.6	21.2	36.1
26	6.2	4.2	9.1	69	28.7	22.0	37.5
27	6.5	4.4	9.4	70	29.9	22.9	39.0
28	6.8	4.7	9.8	71	31.1	23.8	40.7
29	7.1	4.9	10.2	72	32.4	24.8	42.5
30	7.4	5.1	10.6	73	33.9	25.8	44.4
31	7.7	5.4	10.9	74	35.5	27.0	46.6
32	8.0	5.6	11.3	75	37.2	28.2	49.1
33	8.3	5.9	11.7	76	39.2	29.6	51.9
34	8.6	6.2	12.1	77	41.4	31.1	55.1
35	9.0	6.4	12.6	78	44.0	32.8	58.9
36	9.3	6.7	13.0	79	47.0	34.8	63.6
37	9.7	7.0	13.4	80	50.7	37.1	69.5
38	10.0	7.3	13.9	81	55.5	39.8	77.5
39	10.4	7.6	14.3	82	62.3	43.2	89.9
40	10.8	7.9	14.8	83	73.8	47.6	114.6
41	11.2	8.2	15.2	84	>73.8	>47.6	>114.6
42	11.6	8.5	15.7				

MPN is per ml of the 10 ml sample added to the media tube (pour-off is accounted for).

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

REFERENCES:

- *SimPlate for HPC Water Testing Product Information.* (n.d.). Retrieved from www.idexx.ca/view/xhtml/en_ca/water/simplate.jsf.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Dec2013	Initial Release	A.Darrach
1.1	31March2016	Reviewed – No Changes	C. Russell
2.0	12-May-2017	Reviewed and revised. Safety precautions and reagent storage requirements added; Updated format; New document number (old number MIC52615)	L. Steven

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.