



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Number: MIC60060

Version No: 1.0

Page: 1 of 3

Distribution:

Microbiology Quality Control Manual

Effective:

Date Reviewed:

Next Review:

Document Name: Microbiology Stain Quality Control

Approved By:

Status: **DRAFT**

PURPOSE: To standardize the quality control procedures for microbiology stains; Gram-stain and acridine orange.

REAGENTS and/or MEDIA:

- BD Crystal Violet
- BD Gram Iodine (Stabilized)
- BD Decolorizer
- BD Safranin
- Dalynn Biologicals Acridine Orange

SUPPLIES:

- Remel™ QC-Slide Gram Stain Control

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

PROCEDURE INSTRUCTIONS:

Step	Action
Performing Gram-stain quality control	
1	Gram-stain quality control is performed daily with the first run of patient slides. During the morning start-up, a QC slide will be placed on the staining rack for this purpose.
2	A QC order is auto-generated in TQC daily.
3	Stain the QC slide with the patient slide(s). Refer to MIC20115 – Gram Stain Procedure.
4	Read the QC slide before the patient’s slide(s) to ensure quality of the staining.
5	Enter the results into TQC. Refer to MIC60110.

EXPECTED RESULTS:

Positive	Gram-positive cocci
Negative	Gram-negative bacilli

Step	Action
Performing acridine orange quality control	
1	Acridine orange quality control is performed as-tested with patient slide.
2	A QC order is auto-generated in TQC when the stain is ordered for a specimen.
3	Stain a QC slide with the patient slide. Refer to MIC20100-Acridine Orange stain.
4	Read the QC slide before the patient’s slide to ensure quality of the staining.
5	Enter the results into TQC. Refer to MIC60110.

EXPECTED RESULTS:

Positive	Fluorescent cocci
Negative	Fluorescent bacilli

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

REFERENCES:

- Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016.
- Dalynn Biologicals Acridine Orange Stain package insert, 2014.
- BD Gram Stain Kits and Reagents package insert, 2017.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0		Initial Release	L. Steven

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: