

**Stanton Territorial Hospital**P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Number: MIC53400

Version No: 1.0

Page: 1 of 4

Distribution:
Microbiology Test Manual**Effective:**
Date Reviewed:
Next Review:**Document Name: Aerotolerance Test****Approved By:****Status: DRAFT**

PURPOSE: The aerotolerance test is used to designate an organism as an anaerobe.

SAMPLE INFORMATION:

Type	Organisms growing on anaerobic culture media and not present on aerobic culture media.
-------------	--

REAGENTS and/or MEDIA:

- Chocolate agar
- Brucella agar

SUPPLIES:

- Disposable inoculation needles
- Glass microscope slides
- Anaerobic trays or jars
- Anaerobic packs and indicators
- 35° ambient air and 37° CO₂ incubators

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

Document Name: Aerotolerance Test	Document Number: MIC53400	
	Version No: 1.0	Page: 2 of 4
	Effective:	

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to MIC60010 – Microbiology Quality Control for anaerobic tray and jar quality control.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

PROCEDURE INSTRUCTIONS:

Step	Action
Performing the aerotolerance test	
1	Subculture a single, well-isolated colony of each morphological type suspected to be an anaerobe. Using one colony ensures that the same organism goes onto both agar plates and the slide.
2	Touch each colony with a loop or sterile stick and subculture to Chocolate agar, Brucella agar and a labelled, glass microscope slide.
3	The Chocolate agar should be inoculated first, so that if only the Brucella agar plate grows there is no question of not having enough organisms to initiate growth.
4	Incubate the Chocolate agar in CO ₂ incubator for 24 hours. If no growth is seen, re-incubate for an additional 24 hours.
5	Incubate the Brucella agar in anaerobic jar or tray for 48 hours along with original anaerobic media.

INTERPRETATION OF RESULTS:

IF	THEN
Isolate grows on Chocolate agar at 24 or 48 hours AND Isolate grows on Brucella agar at 24 or 48 hours.	Organism is not a strict anaerobe and can be identified using routine methods.
Isolate does not grow on Chocolate agar at 24 or 48 hours AND Isolate does grow on Brucella agar at 24 or 48 hours.	Organism is a strict anaerobe.
Isolate does not grow on Chocolate agar at 24 or 48 hours AND Isolate does not grow on Brucella agar at 24 or 48 hours.	Test is inconclusive: <ul style="list-style-type: none"> • Reasons: organism died or media not inoculated. • Repeat aerotolerance testing with colony on original anaerobic media.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

LIMITATIONS:

1. Chocolate agar should be used for aerotolerance testing. *Haemophilus* spp. will grow anaerobically on Brucella agar and therefore will be mistaken for anaerobic Gram-negative rods if Chocolate agar is not used.
2. Some authorities have suggested performing aerotolerance testing in multiple environments (ambient air, CO₂, a microaerophilic environment) on problem isolates. This is not necessary for identification of the most commonly isolated anaerobes, but it may be considered for isolates when their exact atmospheric requirement is difficult to determine using the above procedure.
3. If anaerobic organism processing occurs on the open bench, all plates should be promptly incubated anaerobically as some clinical isolates may die after relatively short exposure to oxygen.

REFERENCES:

- Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016.
- Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. 2015. Manual of Clinical Microbiology, 11th edition, ASM Press, Washington, D.C.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	23 MAR 19	Initial Release	L. Steven

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: