

Document Name: Beta - Lactamase Test

Approved By:

Status: **DRAFT**

PURPOSE: The beta-lactamase test is used to rapidly detect the production of Beta-lactamase in colonies of *Neisseria gonorrhoeae*, *Moraxella catarrhalis*, *Staphylococcus* spp., *Enterococcus* spp., *Haemophilus influenzae* and anaerobic bacteria.

SAMPLE INFORMATION:

Type	Few, well isolated colonies
-------------	-----------------------------

REAGENTS and/or MEDIA:

Type	BD BBL Cefinase Disks
Stability and Storage Requirements	<ul style="list-style-type: none"> • Store unopened package at -20°C to +8°C. • After opened, store in air tight container containing desiccant at -20°C to +8°C. • Discard remaining disks 60 days after opening blister packaging. • Do not use the cartridge if the disks appear orange or red in color.

SUPPLIES:

- Glass microscope slides
- Forceps
- Sterile saline
- Sterile pipettes
- Wooden sticks
- Disposable loops

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

Document Name: Beta - Lactamase Test	Document Number: MIC50300	
	Version No: 1.0	Page: 2 of 4
	Effective: DRAFT	

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Quality control is performed weekly:
 - Positive: *Staphylococcus aureus* ATCC 29213
 - Negative: *Haemophilus influenzae* ATCC 10211
- A TQC order is automatically generated on Wednesdays to record the QC results.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

PROCEDURE INSTRUCTIONS:

Step	Action																				
Performing the beta-lactamase test																					
1	Using forceps, remove the required number of disks from the dispenser and place on a microscope slide. Use 1 disk per organism.																				
2	Using a sterile pipette, moisten each disk with a drop of sterile saline.																				
3	With a disposable loop or wooden stick, pick several similar colonies from the agar plate and smear onto the surface of the disk.																				
4	Observe the disk for up to 5 minutes for a colour change. For <i>Staphylococci</i> , observe the disk for up to 60 minutes.																				
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="background-color: #00b050; color: white;">Organism</th> <th style="background-color: #00b050; color: white;">Result</th> <th style="background-color: #00b050; color: white;">Time</th> <th style="background-color: #00b050; color: white;">Interpretation</th> </tr> </thead> <tbody> <tr> <td><i>Staphylococcus aureus</i></td> <td>Positive</td> <td>1 hour</td> <td>Resistant to penicillin, ampicillin, carbenicillin. Probably susceptible to cephalothin, methicillin, oxacillin, nafcillin and other penicillinase-resistant penicillins.</td> </tr> <tr> <td><i>Enterococcus faecalis</i></td> <td>Positive</td> <td>5 minutes</td> <td>Resistant to penicillin and ampicillin.</td> </tr> <tr> <td><i>Haemophilus influenzae</i></td> <td>Positive</td> <td>1 minute</td> <td>Resistant to ampicillin. Susceptible to cephalosporins.</td> </tr> <tr> <td><i>Neisseria gonorrhoeae</i> <i>Moraxella catarrhalis</i></td> <td>Positive</td> <td>1 minute</td> <td>Resistant to penicillin.</td> </tr> </tbody> </table>	Organism	Result	Time	Interpretation	<i>Staphylococcus aureus</i>	Positive	1 hour	Resistant to penicillin, ampicillin, carbenicillin. Probably susceptible to cephalothin, methicillin, oxacillin, nafcillin and other penicillinase-resistant penicillins.	<i>Enterococcus faecalis</i>	Positive	5 minutes	Resistant to penicillin and ampicillin.	<i>Haemophilus influenzae</i>	Positive	1 minute	Resistant to ampicillin. Susceptible to cephalosporins.	<i>Neisseria gonorrhoeae</i> <i>Moraxella catarrhalis</i>	Positive	1 minute	Resistant to penicillin.
	Organism	Result	Time	Interpretation																	
	<i>Staphylococcus aureus</i>	Positive	1 hour	Resistant to penicillin, ampicillin, carbenicillin. Probably susceptible to cephalothin, methicillin, oxacillin, nafcillin and other penicillinase-resistant penicillins.																	
	<i>Enterococcus faecalis</i>	Positive	5 minutes	Resistant to penicillin and ampicillin.																	
<i>Haemophilus influenzae</i>	Positive	1 minute	Resistant to ampicillin. Susceptible to cephalosporins.																		
<i>Neisseria gonorrhoeae</i> <i>Moraxella catarrhalis</i>	Positive	1 minute	Resistant to penicillin.																		

INTERPRETATION OF RESULTS:

IF	THEN
Yellow to red colour at site of inoculation	Beta-lactamase = Positive
No colour change on the disk	Beta-lactamase = Negative



NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

LIMITATIONS/PRECAUTIONS:

1. The efficacy of this test in predicting the β -lactam resistance of microorganisms other than *Neisseria gonorrhoeae*, *Haemophilus influenzae*, *Moraxella catarrhalis*, *Staphylococci*, *Enterococci* and certain anaerobic bacteria is unproven.
2. Resistance to β -lactam antibiotics has been, on rare occasions, reported in some of the above organisms without the production of β -lactamases. In these cases, resistance mechanisms such as permeability barriers have been postulated. Therefore, the β -lactamase test should be used as a rapid supplement and not a replacement for conventional susceptibility testing.
3. For some strains of *Staphylococci*, particularly *S.epidermidis*, an inducible β -lactamase has been described that might result in a false-negative β -lactamase reaction with a strain which is resistant to penicillin or ampicillin.

REFERENCES:

- BD BBL, Paper Disks for the Detection of β -Lactamase Enzymes package insert, 2018-09
- Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	4 APR 19	Initial Release	L. Steven

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: