

**Stanton Territorial Hospital**P.O. Box 10, 550 Byrne Road  
YELLOWKNIFE NT X1A 2N1

Document Number: MIC51800

Version No: 1.0

Page: 1 of 5

**Distribution:**

Microbiology Test Manual

**Effective:**

Date Reviewed:

Next Review:

Document Name: CAMP and Reverse CAMP Test

Approved By:

Status: **DRAFT**

**PURPOSE:** The CAMP test is used in the identification of *Streptococcus agalactiae* and many Gram-positive rods, including *Listeria monocytogenes*. The reverse CAMP test is used for the identification of *Arcanobacterium haemolyticum*.

**SAMPLE INFORMATION:****Type**

One, well isolated colony

**REAGENTS and/or MEDIA:**

- Blood agar (BA)

**SUPPLIES:**

- QC organism *Staphylococcus aureus* ATCC 25923
- QC organism *Streptococcus agalactiae* ATCC 12386
- QC organism *Streptococcus pyogenes* ATCC 19615
- Disposable loops
- 35° ambient air and 37° CO<sub>2</sub> incubators

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

<b>Document Name: CAMP and Reverse CAMP Test</b>	<b>Document Number: MIC51800</b>	
	<b>Version No: 1.0</b>	Page: 2 of 5
	<b>Effective: DRAFT</b>	

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.


All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

- Quality control is performed as tested:
  - Positive: *Streptococcus agalactiae* ATCC 12386
  - Negative: *Streptococcus pyogenes* ATCC 19615
- A TQC order is automatically generated when test is ordered to record the QC results.

<b>NOTE:</b> This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME:	Print Date:

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Performing the CAMP test</b>	
1	Streak <i>Staphylococcus aureus</i> ATCC 25923 in a straight line across the center of a Blood agar plate.
2	Streak the unknown organism in the same manner perpendicular to but NOT touching the <i>Staphylococcus</i> streak.
3	Streak the positive control organism parallel to and approximately 2.5 cm from the unknown organism.
4	Streak the negative control organism in the same way on the opposite side of the <i>Staphylococcus</i> streak.
5	<p>Label the identification of each streak on the back of the agar plate:</p> 
6	<p>Incubate the plate overnight at 35° in the CO<sub>2</sub> incubator.</p> <p><b>NOTE:</b> Refrigeration may enhance the reaction after incubation.</p>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

**INTERPRETATION OF RESULTS:**

IF	THEN
Distinct arrowhead of hemolysis at the intersection of the <i>Staphylococcus</i>	<b>CAMP Test = Positive</b>
No enhanced hemolysis at the intersection of the <i>Staphylococcus</i>	<b>CAMP Test = Negative</b>
Distinct arrow of no hemolysis at the intersection of the two haemolytic organisms.	<b>Reverse CAMP Test = Positive</b>

**PROCEDURE NOTES:**

Gram-positive bacilli that are CAMP test positive:

- *Rhodococcus equi*
- *Listeria monocytogenes*
- *Propriosebacterium avidum/granulosum*
- *Actinomyces neuii*
- *Turicella otitidis*
- *Corynebacterium glucuronolyticum*
- *Corynebacterium colyaeae*
- *Corynebacterium imitans*
- Some strains of *Corynebacterium striatum* and *Corynebacterium afermentans* group

Organisms that are reverse CAMP test positive:

- *Corynebacterium pseudotuberculosis*
- *Corynebacterium ulcerans*
- *Arcanobacterium haemolyticum*
- *Clostridium perfringens*

**LIMITATIONS/PRECAUTIONS:**

1. The test is 98% sensitive in detecting *Streptococcus agalactiae*. Isolates with a negative CAMP test may still be *Streptococcus agalactiae* and require further testing.
2. Increased nonspecific hemolysis at the intersections (a *matchstick* effect) may be seen with other Streptococci, but only *Streptococcus agalactiae* produces a definite arrowhead.
3. *Streptococcus pyogenes* can give a reaction that may be interpreted as positive, but it is PYR positive.
4. The CAMP Test separates *Listeria monocytogenes* from most other *Listeria* species.
5. If the agar is too thin or hemolysed, the reaction may be very weak.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date:

**REFERENCES:**

- Clinical Microbiology Procedures Handbook, 4<sup>th</sup> edition, ASM Press, 2016

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	5 APR 19	Initial Release	L. Steven

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

Print Date: