

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC31200 – GBS Screen	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director of Health Services	Date Approved:
Accreditation Canada Applicable Standard:	

**GUIDING PRINCIPLE:**

10-35% of women are asymptomatic carriers of Group B *Streptococcus* (*Streptococcus agalactiae*, GBS) in the genital and gastrointestinal tracts. GBS may be transmitted to the neonate at birth, and remains a leading cause of serious illness and death in newborn populations. A selective enrichment broth must be used for initial processing of prenatal GBS specimens as it significantly improves detection of GBS.

StrepBSelect is a selective, chromogenic, agar-based medium for the isolation and presumptive identification of *Streptococcus agalactiae*. This medium is particularly suited to screening for GBS in vaginal or vaginal-rectal specimens from pregnant women.

**PURPOSE/RATIONALE:**

To screen for Group B Streptococcus (GBS) in vaginal/rectal specimens.

**SCOPE/APPLICABILITY:**

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for GBS screen.

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

### SAMPLE INFORMATION:

<b>Type</b>	Swab <ul style="list-style-type: none"><li>• Amie's with or without charcoal</li></ul>
<b>Source</b>	<ul style="list-style-type: none"><li>• Combined introital (vaginal and anorectal area) swab</li><li>• Vaginal swabs are not the specimen of choice but will be processed</li></ul>
<b>Stability</b>	If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none"><li>• Add specimen quality comment: "Delayed transport may adversely affect pathogen recovery"</li></ul>
<b>Storage Requirements</b>	Room temperature
<b>Criteria for rejection and follow up action</b>	<ol style="list-style-type: none"><li>1. Unlabeled/mislabeled swab</li><li>2. Specimen container label does not match patient identification on requisition</li><li>3. Duplicate specimens obtained with same collection method within 24 hours</li></ol>

### REAGENTS and/or MEDIA:

- LIM Broth (LIM), StrepB *Select* agar (GBS) and Blood agar (BA)
- Identification reagents: catalase and Strep latex test

### SUPPLIES:

- Disposable inoculation needles
- Wooden sticks
- Glass test tubes
- Sterile pipettes

### EQUIPMENT

- Biosafety cabinet
- 35° CO<sub>2</sub> incubator
- 35° ambient air incubator

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

- Refer to MIC60040 – Culture Media Quality Control procedure
- Refer to Test Manual for reagent quality control procedures

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Processing swabs for GBS screen</b>	
<b>1</b>	Place swab in LIM broth, break off leaving the swab in the broth medium and loosely recap. Leave broth in the rack in the BSC.
<b>2</b>	The evening technologist will incubate the rack of LIM broths in the CO <sub>2</sub> incubator before leaving for the evening.
<b>3</b>	At 14:00, after 18 to 24 hours incubation, subculture the broth to StrepB <i>Select</i> agar: <ul style="list-style-type: none"> <li>• Remove the required number of StrepB <i>Select</i> plates from the refrigerator and bring to room temperature</li> <li>• Saturate a sterile swab in the broth and rotate against the wall of the tube above the liquid to remove excess inoculum</li> <li>• Inoculate GBS with the swab</li> <li>• Ensure all surfaces of swab make contact with the agar</li> <li>• Streak for isolated growth using a disposable inoculation needle</li> </ul>
<b>4</b>	Incubate the media: <ul style="list-style-type: none"> <li>• Place GBS in the O<sub>2</sub> incubator in the 1:00 pm new urine rack</li> </ul>

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

**INTERPRETATION OF RESULTS:**

Step	Action	
<b>1</b>	<ul style="list-style-type: none"> <li>Observe GBS plate at 18 to 24 hours</li> <li>Examine for blue colonies</li> </ul>	
<b>2</b>	<b>IF</b>	<b>THEN</b>
	No blue colonies seen at 24 hours	<ul style="list-style-type: none"> <li>Record observations in the LIS</li> <li>Re-incubate plates in O<sub>2</sub> incubator on the "Old urine culture" shelf</li> </ul>
	No blue colonies seen at 38-48 hours	<ul style="list-style-type: none"> <li>Record observations in the LIS</li> <li>Workup complete.</li> <li>GBS not isolated</li> </ul>
	Blue colonies seen at 24 or 38 hours	<ul style="list-style-type: none"> <li>Record observations in the LIS</li> <li>Subculture colonies to BA plate if no isolated colonies are present</li> <li>From BA sub-plate or original GBS plate, perform:                             <ol style="list-style-type: none"> <li>Catalase</li> </ol> </li> </ul>
	<b>IF</b>	<b>THEN</b>
	Catalase POSITIVE	<ul style="list-style-type: none"> <li>Record observations in LIS</li> <li>Workup complete</li> <li>GBS not isolated</li> </ul>
	Catalase NEGATIVE	<ul style="list-style-type: none"> <li>Record observations in LIS</li> <li>Perform Strep latex test for Group B</li> </ul>
	<b>IF</b>	<b>THEN</b>
	Strep B latex test NEGATIVE	<ul style="list-style-type: none"> <li>Record observations in LIS</li> <li>Workup complete</li> <li>GBS not isolated</li> </ul>
	Strep B latex test POSITIVE	<ul style="list-style-type: none"> <li>Record observations in LIS</li> <li>GBS isolated</li> </ul>
	<b>IF</b>	<b>THEN</b>
	Group B screen is positive  Penicillin allergy NOT present	<ul style="list-style-type: none"> <li>Susceptibility testing not performed</li> </ul>
	Group B screen is positive  Penicillin allergy present	<ul style="list-style-type: none"> <li>Perform susceptibility testing as per ASTM</li> </ul>

**NOTE:** Each Streptococcus grouping latex test should be tested with at least one extra grouping latex suspension as a negative control

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

**REPORTING INSTRUCTIONS:**

IF	REPORT
GBS not isolated	<ul style="list-style-type: none"> <li>Report: <b>"No Streptococcus agalactiae (Group B) isolated"</b></li> <li>A copy will automatically print to STH OBS (SOBS)</li> </ul>
GBS isolated and no penicillin allergy indicated in clinical history	<ul style="list-style-type: none"> <li>Choose key 7 on STRB keypad to add isolate: <b>"Streptococcus agalactiae (Group B)"</b></li> <li>List quantitation as <b>"Isolated"</b></li> <li>The following isolate comments will be added: <b>&amp;A336</b> and <b>&amp;IAPO</b></li> <li>A copy will automatically print to STH OBS (SOBS)</li> </ul>
GBS isolated and clinical history indicates penicillin allergy	<ul style="list-style-type: none"> <li>Choose key 8 on the STRB keypad to add isolate: <b>"Streptococcus agalactiae (Group B)"</b></li> <li>List quantitation as <b>"Isolated"</b></li> <li>The following isolate comment will be added: <b>&amp;A336</b></li> <li>KB susceptibility panel for GBS screen is ordered</li> <li>Report susceptibility results as per ASTM</li> <li>A copy will automatically print to STH OBS (SOBS)</li> </ul>

**LIMITATIONS:**

- As with any chromogenic medium, it is important to streak at closely spaced intervals in order to obtain well-isolated colonies: the morphology and colour of the colonies will then be more typical.
- If an inoculum contains a high density of *Streptococcus agalactiae*, the medium around the deposit may be coloured.
- The colonies of some species, other than *Streptococcus agalactiae* (for example, *Streptococcus pyogenes*, *Streptococcus porcinus* and *Streptococcus gallolyticus*) may appear blue.
- The intrinsic demands of some *Streptococci* can lead to absence or partial inhibition of their growth.
- Prenatal GBS screening should be done between 35-37 weeks gestation as GBS colonization can be transient and colonization early in pregnancy is not predictive of early-onset GBS disease. Late third trimester colonization status has been used as a proxy of intrapartum colonization. The negative predictive value of GBS cultures performed  $\leq 5$  weeks before delivery is 95%-98%; however, the clinical utility decreases when a prenatal culture is performed more than five weeks before delivery because the negative predictive value declines.

**CROSS-REFERENCES:**

- MIC60040 Culture Media Quality Control

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

**REFERENCES:**

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4<sup>th</sup>ed.) Washington, D.C.: ASM Press
2. Jorgensen, J. H., Pfaller, M. A., Carroll, K. C., Funke, G., Landry, M. L., Richter, S. S., Warnock, D. W. (2015). *Manual of clinical microbiology* (11<sup>th</sup>ed.). ASM Press
3. BioRad Laboratories. (2009/08). *StrepB Select* package insert

**APPROVAL:**

\_\_\_\_\_  
Date

\_\_\_\_\_

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	15 May17	Initial Release	L. Steven
2.0	8 Feb 19	Updated to include new chromogenic agar <i>StrepB Select</i>	L. Steven
3.0	30 Dec 20	Procedure reviewed and added to NTHSSA policy template	L. Steven

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.