

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC51100 – Beta-Lactamase Test	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director of Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

GUIDING PRINCIPLE:

The beta-lactamase test is used to rapidly detect the production of beta-lactamase in colonies of *Neisseria gonorrhoeae*, *Moraxella catarrhalis*, *Staphylococcus* spp., *Enterococcus* spp., *Haemophilus influenzae* and anaerobic bacteria.

PURPOSE/RATIONALE:

This standard operating procedure describes how to perform the beta-lactamase test.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) performing the beta-lactamase test.

SAMPLE INFORMATION:

Type	Few, well isolated colonies
-------------	-----------------------------

REAGENTS and/or MEDIA:

Type	BD BBL Cefinase Disks
Stability and Storage Requirements	<ul style="list-style-type: none"> • Store unopened package at -20°C to 8°C • After opened, store in air tight container containing desiccant at -20°C to 8°C • Discard remaining disks 60 days after opening blister packaging • Do not use the cartridge if the disks appear orange or red in color

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

SUPPLIES:

- Glass microscope slides
- Forceps
- Sterile saline
- Sterile pipettes
- Wooden sticks
- Disposable loops

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Quality control is performed weekly:
 - Positive: *Staphylococcus aureus* ATCC 29213
 - Negative: *Haemophilus influenzae* ATCC 10211
- A TQC order is automatically generated on Wednesdays to record the QC results

PROCEDURE INSTRUCTIONS:

Step	Action
Performing the beta-lactamase test	
1	Remove disks from the fridge and bring to room temperature.
2	Using forceps remove the required number of disks from the dispenser and place on a microscope slide. Use 1 disk per organism.
3	Using a sterile pipette, moisten each disk with a drop of sterile saline.
4	With a disposable loop or wooden stick, pick several similar colonies from the agar plate and smear onto the surface of the disk.
5	Observe the disk for up to 5 minutes for a colour change. For <i>Staphylococci</i> , observe the disk for up to 60 minutes.

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

	Organism	Result	Time	Interpretation
6	<i>Staphylococcus aureus</i>	Positive	1 hour	Resistant to penicillin, ampicillin, carbenicillin Probably susceptible to cephalothin, methicillin, oxacillin, nafcillin and other penicillinase-resistant penicillins
	<i>Enterococcus faecalis</i>	Positive	5 minutes	Resistant to penicillin and ampicillin
	<i>Haemophilus influenzae</i>	Positive	1 minute	Resistant to ampicillin Susceptible to cephalosporins
	<i>Neisseria gonorrhoeae</i> <i>Moraxella catarrhalis</i>	Positive	1 minute	Resistant to penicillin

INTERPRETATION OF RESULTS:

IF	THEN
Yellow to red colour at site of inoculation	Beta-lactamase = Positive
No colour change on the disk	Beta-lactamase = Negative
	

LIMITATIONS:

1. The efficacy of this test in predicting the β -lactam resistance of microorganisms other than *Neisseria gonorrhoeae*, *Haemophilus influenzae*, *Moraxella catarrhalis*, *Staphylococci*, *Enterococci* and certain anaerobic bacteria is unproven.
2. Resistance to β -lactam antibiotics has been, on rare occasions, reported in some of the above organisms without the production of β -lactamases. In these cases, resistance mechanisms such as permeability barriers have been postulated. Therefore, the β -lactamase test should be used as a rapid supplement and not a replacement for conventional susceptibility testing.
3. For some strains of *Staphylococci*, particularly *S.epidermidis*, an inducible β -lactamase has been described that might result in a false-negative β -lactamase reaction with a strain which is resistant to penicillin or ampicillin.

REFERENCES:

1. BD BBL. (2018-09). *Paper Disks for the Detection of β -Lactamase Enzymes* package insert

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

APPROVAL:

Date

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	04 Apr 19	Initial Release	L. Steven
2.0	30 Jun 21	Procedure reviewed and added to NTHSSA policy template	L. Steven

DRAFT

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.