Issuing Authority: Director of Health Services Next Review Date: Type: Laboratory Services Program SOP

Policy Number: Date Approved:

PROGRAM Standard Operating Procedure – Laboratory Services		
Title: MIC36400 -	Policy Number:	
Referral of Cat B Specimens to DL		
Program Name: Laboratory Services		
Applicable Domain: Lab, DI and Pharmacy Services		
Additional Domain(s):		
Effective Date:	Next Review Date:	
Issuing Authority:	Date Approved:	
Director of Health Services		
Accreditation Canada Applicable Standard:		

GUIDING PRINCIPLE:

Organisms that cannot be identified or require further susceptibility testing need to be referred to *Dyna*LIFE (DL). TDG regulations need to be followed including paperwork and packaging.

PURPOSE/RATIONALE:

To ensure microbiology Category B specimens are sent for referral testing to *Dyna*LIFE appropriately.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) and Medical Laboratory Assistants (MLAs) packaging Category B samples being referred to *DynaLIFE*.

SAMPLE INFORMATION:

The following category B specimens need to be sent to DL for referral testing:

- Any organism that cannot be identified
- Any susceptibility testing that cannot be performed at Stanton microbiology laboratory or that requires confirmation

SUPPLIES:

- Copan Transystem Culture swab transport system (with or without charcoal)
- DynaLIFE Referred-In Isolate Submission Form
- Category B box
- Category B package supplies
- DL RG2 Biosafety fax coversheet

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Policy Number: Date Approved: Page 1 of 5


Issuing Authority: Director of Health Services


Next Review Date:

Type: Laboratory Services Program SOP Policy Number:

Date Approved:

PROCEDURE INSTRUCTIONS:





Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.


Policy Number: Page 2 of 5

Issuing Authority: Director of Health Services

Next Review Date:

Type: Laboratory Services Program SOP

Policy Number: Date Approved:



Step	Action			
Comp	mplete organism information on the DL requisition:			
	Referred-In Isolate Submission Form (LTR47241) Edit Approved By:Hunter, Shane (11/20/2015)	$DynaLIFE_{Dx}$ Revision:3.00		
	DynaLIFEDx Suite 200, 10150 - 102 Street, Edmonton, AB T5J 5E5 Toll Free: 1-800-661-9876 Fax: (780) 453-9425 (Microbiology)			
	REFERRED IN ISOLATE SUBMI	SSION FORM [Complete all sections]		
	Patient Name: Patient's PHN: Patient's DOB: Patient Address (Required for Notifiable Organisms):	Submitting Facility Address: Stanton Terriforial Hospital 548 Bryne Road, P.O. Box 10 Yellowknife, NT X1A ZN1 Phone: (876) 767-9300 x 46691 Fax: (867) 669-4141		
	City: Postal Code:	Submitter ID 11122		
	Phone #:	Physician's Code: Report Location: Submitting Lab Number:		
	Clinical History:	Specimen Source:		
4	Current Antibiotics:	Specimen Collection Date:		
4	Medical Officer of Health and/or Infection Control notified :	Submission Date:		
	☐ Yes ☐ No (if required)			
	SUPPORTING DATA [C	omplete all relevant sections]		
	Test(s) requested ☐ Gram Smear Interpretation ☐ Ident	ification		
		eeze Only **Organism Identification:		
	Was this a pure culture? ☐ Yes ☐ No (list other isolates)			
	Direct Smear Results:	Growth Requirements: ☐ Acrobe ☐ Arnaerobe		
	Oxidase: Catalase: Other Tests:	Gram Stain:		
	Commercial Identification System & Susceptibility Results:			
	Identification (Suspected) :			
	Vegreville only Pre-incubated urine: Colony count Colony type			

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Issuing Authority: Director of Health Services

Next Review Date:

Type: Laboratory Services Program SOP Policy Number:

Date Approved:

Step	Action
Comp	lete highlighted areas on the DL RG2 coversheet:
	Stanton Territorial Hospital – Microbiology Laboratory 548 Byrne Road, PO Box 10 - Yellowknife, NT X1A 2N1 Phone: 867 767-9300 [ext.46691] Fax: 867 669-4141
	To: DynaLife Dx Microbiology Laboratory From: STH Microbiology Laboratory
	Fax: 1 (780) 453 9425 Pages:
5	Phone: 1 (800) 661 9876 Date:
	Re: Transfer of Human Pathogens (L-R2-09731-22-BK-00) CC: sth biosafety@gov.nt.ca
	STH Microbiology laboratory is transferring a package containing a Risk Group 2 Human pathogen to your facility The sample is to be directed to the

Step	Action	
Order a referred test in result entry screen:		
6	 Select Add Test Choose appropriate reference code: ?REFD- DynaLIFE Select Ok Add a period (.) in the results line Final report the test line 	

Step	Action			
Comp	Complete Category B send-out process:			
7	 Print off any Vitek results for organism being referred Submit organism on a labelled agar plate sealed with parafilm or on a labelled C&S swab Photocopy all papers and staple together. Hole punch photocopy stack and place in the <i>Dyna</i>LIFE and ProvLab Pending Referrals binder. Ensure the report is filed under the correct tab (<i>Dyna</i>LIFE, Cat B) Place a requisiton label on the DL requisition and scan into SoftMedia Pack specimen up according to TDG Category B regulations, send original paperwork with specimen Email completed coversheet to email locations listed Place Category B box in the Category B overpack box in the lab assitant area for shipment to <i>Dyna</i>LIFE 			

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Policy Number: Page 4 of 5

Title: MIC36400-Referral of Category B Specimens to DL

Issuing Authority: Director of Health Services

Next Review Date:

Type: Laboratory Services Program SOP
Policy Number:
Date Approved:

REFERENCES:

• DynaLIFE DX Referred-In Isolate Submission Form (LTR47241), revision 3.00

APPROVAL:		
Date		

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	28 Apr 17	Initial Release	L. Steven
2.0	16 Jun 19	Updated to reflect new hospital address	L. Steven
3.0	10 Mar 20	Procedure reviewed	L. Steven
4.0	30 Jun 22	Procedure reviewed and added to NTHSSA policy template	L. Steven

ref

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Policy Number: Date Approved: Page 5 of 5