

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC32100 – Yeast Culture	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director, Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

GUIDING PRINCIPLE:

Yeast cultures are performed to identify yeast in select specimens. Yeast isolates are identified and reported based on their full identification and clinical significance in the location of isolation.

PURPOSE/RATIONALE:

To determine the presence or absence of yeast in anal, penis, cervical and vaginal specimens.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for yeast culture.

SAMPLE INFORMATION:

Type	Swab <ul style="list-style-type: none"> • Amie’s with or without charcoal
Source	<ul style="list-style-type: none"> • Anus • Penis • Cervix • Vagina
Stability	If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none"> • Add specimen quality comment: “Delayed transport may adversely affect pathogen recovery”
Storage Requirements	Room temperature
Criteria for rejection	<ol style="list-style-type: none"> 1. Unlabeled/mislabeled swabs 2. Specimen container label does not match patient identification on requisition

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REAGENTS and/or MEDIA:

- Sabouraud Dextrose (SAB) agar

SUPPLIES:

- Disposable inoculation needles
- Microscope slides
- Sterile saline
- Wooden sticks
- Coverslips

EQUIPMENT

- Biosafety cabinet
- Vitek 2 and supplies

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to Test Manual for reagent quality control procedures

PROCEDURE INSTRUCTIONS:

Step	Action
Processing specimens for yeast culture	
1	In the biosafety cabinet: <ul style="list-style-type: none">• Inoculate SAB with the swab• Ensure all surfaces of swab make contact with the agar• Streak for isolated growth using a disposable inoculation needle• Write on the plate the date of the 48-hour read
2	Incubate the media: <ul style="list-style-type: none">• Place SAB on the workbench

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

INTERPRETATION OF RESULTS:

Step	Action	
1	<ul style="list-style-type: none"> Observe SAB plate at 48 hours Examine for white, creamy colonies resembling yeast 	
2	IF	THEN
	Colonies resembling yeast NOT seen	<ul style="list-style-type: none"> Workup complete Yeast not isolated
3	IF	THEN
	Wet prep NEGATIVE	<ul style="list-style-type: none"> Workup complete Yeast not isolated
	Wet prep POSITIVE	<ul style="list-style-type: none"> Perform Vitek 2 YST card

REPORTING INSTRUCTIONS:

IF	REPORT
Yeast not isolated	<ul style="list-style-type: none"> Report: "No Yeast Isolated"
Yeast Isolated, Candida albicans	<ul style="list-style-type: none"> Report: "Candida albicans" List quantitation as "Isolated"
Yeast Isolated, not Candida albicans	<ul style="list-style-type: none"> Report: "Yeast (NOT Candida albicans)" List quantitation as "Isolated"

REFERENCES:

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4thed.) Washington, D.C.: ASM Press
2. Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology*, 11th edition. Washington, D.C: ASM Press

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

APPROVAL:

Date

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	30 Jan 21	Initial Release	L. Steven

DRAFT