

| PROGRAM Standard Operating Procedure – Laboratory Services | |
|-------------------------------------------------------------------|-------------------|
| Title: MIC32200 – Nasal Culture | Policy Number: |
| Program Name: Laboratory Services | |
| Applicable Domain: Lab, DI and Pharmacy Services | |
| Additional Domain(s): | |
| Effective Date: | Next Review Date: |
| Issuing Authority: Director of Health Services | Date Approved: |
| Accreditation Canada Applicable Standard: N/A | |

GUIDING PRINCIPLE:

Nasal swabs are performed to identify nasal carriers of *Staphylococcus aureus*.

PURPOSE/RATIONALE:

To determine the presence or absence of *Staphylococcus aureus* in nasal specimens.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for nasal culture.

SAMPLE INFORMATION:

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Type | Swab <ul style="list-style-type: none"> Amie's with or without charcoal |
| Source | Nose |
| Stability | If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none"> Add specimen quality comment: "Delayed transport may adversely affect pathogen recovery" |
| Storage Requirements | Room temperature |
| Criteria for rejection | <ol style="list-style-type: none"> Unlabeled/mislabeled swabs Specimen container label does not match patient identification on requisition Duplicate specimens obtained with same collection method within 24 hours |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REAGENTS and/or MEDIA:

- Blood agar (BA) and Mueller Hinton agar (MHP)
- Identification reagents: catalase, Staph latex test and Cefoxitin antibiotic disks

SUPPLIES:

- Disposable inoculation needles
- Wooden sticks

EQUIPMENT:

- Biosafety cabinet
- 35° CO₂ incubator

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to Test Manual for reagent quality control procedures

PROCEDURE INSTRUCTIONS:

| Step | Action |
|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Processing specimens for nasal culture | |
| 1 | In the biosafety cabinet: <ul style="list-style-type: none">• Inoculate BA with the swab• Ensure all surfaces of swab make contact with the agar• Streak for isolated growth using a disposable inoculation needle |
| 2 | Incubate the media: <ul style="list-style-type: none">• Place BA in the CO₂ incubator |

INTERPRETATION OF RESULTS:

| Step | Action | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | <ul style="list-style-type: none"> Observe BA plate at 24 hours and 48 hours Examine for colonies resembling <i>Staphylococcus aureus</i> | |
| 2 | IF | THEN |
| | No <i>S.aureus</i> colonies seen at 24 hours | <ul style="list-style-type: none"> Record observations in the LIS Re-incubate plate in CO₂ incubator on the "Old wound culture" shelf |
| 3 | IF | THEN |
| | No <i>S.aureus</i> colonies seen at 48 hours | <ul style="list-style-type: none"> Record observations in the LIS Workup complete <i>S.aureus</i> not isolated |
| 3 | IF | THEN |
| | <i>S.aureus</i> is present on BA but not isolated | <ul style="list-style-type: none"> Subculture colonies to BA-S Perform Staph latex test |
| 4 | IF | THEN |
| | <i>S.aureus</i> is present on BA and isolated | <ul style="list-style-type: none"> Perform Staph latex test |
| | IF | THEN |
| | Staph latex test NEGATIVE | <ul style="list-style-type: none"> Record observations in the LIS Workup complete <i>S.aureus</i> not isolated |
| | Staph latex test POSITIVE | <ul style="list-style-type: none"> Record observations in the LIS <i>S.aureus</i> isolated Perform cefoxitin disk diffusion test |
| | IF | THEN |
| Cefoxitin screen SENSITIVE | <ul style="list-style-type: none"> Record observations in the LIS Methicillin sensitive <i>S.aureus</i> isolated | |
| Cefoxitin screen RESISTANT | <ul style="list-style-type: none"> Record observations in the LIS Methicillin resistant <i>S.aureus</i> isolated | |

REPORTING INSTRUCTIONS:

| IF | REPORT |
|-------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>Staphylococcus aureus</i> not isolated | <ul style="list-style-type: none"> Report: "No Staphylococcus aureus isolated" |
| Methicillin sensitive <i>Staphylococcus aureus</i> isolated | <ul style="list-style-type: none"> Add organism: "Staphylococcus aureus" List quantification as "Isolated" Report organism with isolate comment &MSSA |
| Methicillin resistant <i>Staphylococcus aureus</i> isolated | <ul style="list-style-type: none"> Add organism: "Staphylococcus aureus" List quantification as "Isolated" Report organism with isolate comment &cx01 In order entry, copy report to OCPHO (HPU1) Check the home address of the patient. If from Nunavut: Copy report to the applicable NU CPHO In order entry, copy report to appropriate IPAC ward if ER or In-patient In order entry add ESO code "MRSA" |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

NOTE:

STH IPAC ward is **SIPAC**. IRH IPAC ward is **IIPAC**. Territorial IPAC ward is **TIPAC**

REFERENCES:

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4thed.) Washington, D.C.: ASM Press
2. Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology, 11th edition*. Washington, D.C: ASM Press

APPROVAL:

Date

REVISION HISTORY:

| REVISION | DATE | Description of Change | REQUESTED BY |
|----------|-----------|--------------------------------------------------------|--------------|
| 1.0 | 03 Mar 19 | Initial Release | L. Steven |
| 2.0 | 22 Feb 21 | Procedure reviewed and added to NTHSSA policy template | L. Steven |
| | | | |
| | | | |
| | | | |
| | | | |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.