

<b>PROGRAM Standard Operating Procedure – Laboratory Services</b>	
Title: MIC32500 – Eye Culture-Superficial	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director of Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

**GUIDING PRINCIPLE:**

Inflammatory eye conditions can be separated into several clinical syndromes. A variety of microorganisms can play major roles, both in acute and in chronic conditions.

**PURPOSE/RATIONALE:**

To determine the presence or absence of bacterial pathogens in superficial eye specimens.

**SCOPE/APPLICABILITY:**

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for superficial eye culture.

**SAMPLE INFORMATION:**

<b>Type</b>	Swab <ul style="list-style-type: none"><li>• Amie’s with or without charcoal</li></ul>
<b>Source</b>	Conjunctiva: inflammation of the conjunctiva (the mucous membrane covering the sclera) <ul style="list-style-type: none"><li>➤ Swab of conjunctiva surface / pus</li></ul>
<b>Stability</b>	If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none"><li>• Add specimen quality comment: “Delayed transport may adversely affect pathogen recovery”</li></ul>
<b>Storage Requirements</b>	Room temperature

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

**Criteria for rejection**

1. Unlabeled/mislabelled swabs
2. Specimen container label does not match patient identification on requisition

**NOTE:**

- If gonorrhoeae culture is ordered on eye specimen, superficial eye culture along with gonorrhoeae culture will be performed. Refer to MIC33500-Gonorrhoeae Culture
- Refer to MIC34100-Body Fluid Culture for intraocular fluid
- Refer tissue or biopsy specimens for culture to *DynaLIFE*

**REAGENTS and/or MEDIA:**

- Blood agar (BA) and Chocolate agar (CHO)
- Identification reagents: catalase, oxidase, Staph latex test, Strep latex test, etc.

**SUPPLIES:**

- Disposable inoculation needles
- Microscope slides
- Wooden sticks

**EQUIPMENT:**

- Biosafety cabinet
- 35° ambient air and 35° CO<sub>2</sub> incubators
- Vitek 2 and supplies

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Ensure that appropriate hand hygiene practices be used.
- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

- Refer to Test Manual for reagent quality control procedures

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Processing specimens for superficial eye culture</b>	
<b>1</b>	In the biosafety cabinet: <ul style="list-style-type: none"> <li>• Inoculate BA and CHO with the swab</li> <li>• Ensure all surfaces of swab make contact with the agar</li> <li>• Streak for isolated growth using a disposable inoculation needle</li> <li>• Prepare smear by rolling the swab gently across the slide to avoid destruction of cellular elements and disruption of bacterial arrangements</li> </ul>
<b>2</b>	Incubate the media: <ul style="list-style-type: none"> <li>• Place BA and CHO in the CO<sub>2</sub> incubator</li> </ul>
<b>3</b>	Allow smear to dry and perform gram stain. Gram stain must be read before culture plates. Refer to MIC20115-Gram stain procedure.

<b>Probable Pathogens</b>	
<ul style="list-style-type: none"> <li>• <i>Haemophilus influenzae</i></li> <li>• <i>Staphylococcus aureus</i></li> <li>• <i>Streptococcus pneumoniae</i></li> <li>• <i>Streptococcus pyogenes</i></li> </ul>	<ul style="list-style-type: none"> <li>• <i>Moraxella</i> spp.</li> <li>• <i>Pseudomonas aeruginosa</i></li> <li>• <i>Neisseria gonorrhoeae</i></li> <li>• <i>Neisseria meningitidis</i></li> </ul>
<b>Potential Pathogens</b>	
<ul style="list-style-type: none"> <li>• Aerobic gram-negative-bacilli</li> <li>• <i>Enterococcus</i> spp.</li> </ul>	<ul style="list-style-type: none"> <li>• <i>Haemophilus parainfluenzae</i></li> </ul>
<b>Commensal Skin Flora</b>	
<ul style="list-style-type: none"> <li>• Coagulase-negative Staphylococcus</li> <li>• Micrococcus spp.</li> <li>• Corynebacterium spp.</li> </ul>	<ul style="list-style-type: none"> <li>• Bacillus spp. not listed above</li> <li>• Nonpathogenic Neisseria spp.</li> <li>• viridans Streptococcus grp.</li> <li>• Anaerobes</li> </ul>

**INTERPRETATION OF RESULTS:**

Step	Action
<b>1</b>	Ensure growth on culture media correlates with gram stain results. If discordant results are found between the gram stain and growth: <ul style="list-style-type: none"> <li>• Re-examine smear and culture plates</li> <li>• Check for anaerobic growth</li> <li>• Re-incubate media to resolve</li> <li>• Consider re-smearing or re-planting specimen</li> </ul>
<b>2</b>	<ul style="list-style-type: none"> <li>• Observe BA and CHO at 24 hours and 48 hours</li> </ul>
<b>3</b>	<ul style="list-style-type: none"> <li>• <u>If organism is a probable pathogen:</u> <ul style="list-style-type: none"> <li>➢ Perform and report full identification</li> <li>➢ Perform and report susceptibility testing as per ASTM</li> </ul> </li> </ul>
<b>4</b>	<ul style="list-style-type: none"> <li>• <u>If organism is a potential pathogen:</u> <ul style="list-style-type: none"> <li>➢ Perform and report full identification if there are ≤3 different potential pathogens</li> </ul> </li> </ul>

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

	<ul style="list-style-type: none"> <li>➤ Perform and report susceptibility testing on potential pathogens as per ASTM if any of the following is true:                             <ul style="list-style-type: none"> <li>○ 3-4+WBC were seen in the gram stain</li> <li>○ Growth is pure or predominant</li> <li>○ Patient is immunocompromised</li> </ul> </li> <li>➤ If &gt;3 potential pathogens are present, list and do not perform or report susceptibility testing</li> </ul>
<b>5</b>	Perform minimal identification and list commensal skin flora.

**REPORTING INSTRUCTIONS:**

IF	REPORT
No growth after 1 day	<b>PRELIM:</b> <ul style="list-style-type: none"> <li>• Report: <b>"No Growth After 1 Day"</b></li> <li>• Report: <b>"Further report to follow"</b></li> </ul>
No growth after 2 days	<b>FINAL:</b> <ul style="list-style-type: none"> <li>• Report: <b>"No Growth After 2 Days"</b></li> </ul>
Mix of commensal conjunctival flora	<ul style="list-style-type: none"> <li>• Report: <b>"Mixed commensal conjunctival flora"</b></li> <li>• List quantitation</li> </ul>
Mix of enteric Gram-negative bacilli	<ul style="list-style-type: none"> <li>• Report: <b>"Mixture of coliform organisms"</b></li> <li>• List quantitation</li> </ul>
Mix of anaerobic organisms	<ul style="list-style-type: none"> <li>• Report: <b>"Mixture of anaerobic organisms"</b></li> <li>• List quantitation</li> </ul>
Mix of non-pathogenic organisms	<ul style="list-style-type: none"> <li>• Report: <b>"Commensal flora"</b></li> <li>• List quantitation</li> </ul>
Growth of potential pathogens where minimal identification and listing is required	<ul style="list-style-type: none"> <li>• Report the minimal identification (i.e., Gram Negative Bacilli - Lactose Fermenter)</li> <li>• List quantitation</li> </ul>
Growth of potential pathogen where full identification is required	<ul style="list-style-type: none"> <li>• Report organisms identification</li> <li>• List quantitation</li> <li>• Report susceptibility as per interpretation of results</li> </ul>
Growth of pathogen	<ul style="list-style-type: none"> <li>• Report organisms full identification</li> <li>• List quantitation</li> <li>• Report susceptibility results as per ASTM</li> </ul>
<i>Neisseria gonorrhoeae</i> isolated and gonorrhoeae culture was not ordered	<ul style="list-style-type: none"> <li>• Add organism: <b>"Neisseria gonorrhoeae"</b></li> <li>• List quantification as: <b>"Isolated"</b></li> <li>• Add Beta-lactamase result if positive</li> <li>• Add isolate comment <b>&amp;REF6</b></li> <li>• Refer isolate to APL for susceptibility testing</li> <li>• Freeze isolate and log into stored isolates log</li> </ul>

**NOTE:**

- Refer to Reportable Diseases – Public Health Act as of September 2009 for reporting to OCPHO (HPU1)
- Refer to LQM70620-Laboratory Critical Results List-Microbiology for results that need to be phoned to ordering location
- Refer to MIC36100-Nosocomial Infection Notification Job Aid to determine if organism needs to be copied to Infection Prevention and Control

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

- Refer to MIC36200-Referral of Category A Specimens to APL for sending category A isolates to APL
- Refer to MIC36300-Referral of Category B Specimens to APL for sending isolates to APL
- Refer to MIC36400-Referral of Category B Specimens to DL for sending isolates to *DynaLIFE*.

#### **LIMITATIONS:**

1. False positive cultures can result from contamination of the specimen or plates with skin flora.
2. False negative results can occur if antimicrobial agents are given prior to collection of the specimen.
3. Even with the best techniques, culture often fails to yield the infecting organism.

#### **CROSS-REFERENCES:**

- LQM70620-Laboratory Critical Results List-Microbiology
- MIC20115-Gram stain procedure
- MIC33500-Gonorrhoeae Culture
- MIC34100-Body Fluid Culture
- MIC36100-Nosocomial Infection Notification Job Aid
- MIC36200-Referral of Category A Specimens to APL
- MIC36300-Referral of Category B Specimens to APL
- MIC36400-Referral of Category B Specimens to DL

#### **REFERENCES:**

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4<sup>th</sup>ed.) Washington, D.C.: ASM Press
2. Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology*, 11<sup>th</sup> edition. Washington, D.C: ASM Press

**APPROVAL:**

\_\_\_\_\_  
Date

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	28 May 18	Initial Release	L. Steven
2.0	26 Feb 21	Procedure reviewed and added to NTHSSA policy template	L. Steven

DRAFT

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.