

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC31100 – VRE Screen	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director, Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

GUIDING PRINCIPLE:

Specimens are submitted to identify carriers of vancomycin-resistant *Enterococcus faecium* and/or *Enterococcus faecalis* (VRE). Swabs may be submitted from any body site (other than nasal and axilla), but most commonly are collected from the rectum and are used to aid in the prevention and control of VRE in healthcare settings

PURPOSE/RATIONALE:

This standard operating procedure describes the screening for Vancomycin Resistant Enterococci (VRE) on admission and as part of Multi-Resistant Organism (MRO) screens.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for VRE screen.

SAMPLE INFORMATION:

Type	Swab <ul style="list-style-type: none"> • Amie’s with or without charcoal
Source	<ul style="list-style-type: none"> • Rectum • Stool • MRO screen: any site
Stability	If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none"> • Add specimen quality comment: “Delayed transport may adversely affect pathogen recovery”
Storage Requirements	Room temperature

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Criteria for rejection	<ol style="list-style-type: none">1. Unlabeled/mislabeled swabs2. Specimen container label does not match patient identification on requisition3. Duplicate specimens obtained with same collection method from same collection location within 24 hours
-------------------------------	--

REAGENTS and/or MEDIA:

- VREselect agar (VRE), Blood agar (BA) and Mueller Hinton agar (MH)
- Identification reagents: gram stain, catalase and vancomycin E-test

SUPPLIES:

- Disposable inoculation needles
- Wooden sticks

EQUIPMENT:

- Biosafety cabinet
- 35° ambient air incubator
- Microscope
- Vitek 2 and supplies

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures:

- Ensure that appropriate hand hygiene practices be used
- Lab gown must be worn when performing activities with potential pathogens
- Gloves must be worn when direct skin contact with infected materials is unavoidable
- Eye protection must be used when there is a known or potential risk of exposure of splashes
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC)
- The use of needles, syringes and other sharp objects should be strictly limited

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to MIC60040-Culture Media Quality Control procedure
- Refer to Test Manual for reagent quality control procedures

PROCEDURE INSTRUCTIONS:

Step	Action
Processing swabs for VRE screen	
1	Monday to Sunday, VRE swabs are processed by 15:00.

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

2	In the biosafety cabinet: <ul style="list-style-type: none"> • Inoculate VREselect agar with the swab • Ensure all surfaces of the swab make contact with the agar • Streak for isolated growth using a disposable inoculation needle
3	Label the VRE plate with: R (Date + 2 days).
4	Incubate the media: <ul style="list-style-type: none"> • Place VRE in the O₂ incubator in appropriate tray

INTERPRETATION OF RESULTS:

Step	Action	
1	<ul style="list-style-type: none"> • Observe VRE plate at 18 to 24 hours and 38 to 48 hours • Examine for blue or pink colonies 	
2	If no blue or pink colonies are seen at 18 to 24 hours: <ul style="list-style-type: none"> • Record observations in the LIS • Re-incubate plate in O₂ incubator on the "Old urine culture" shelf 	
3	If no blue or pink colonies are seen at 38 to 48 hours: <ul style="list-style-type: none"> • Record observations in the LIS • Workup complete, VRE not isolated 	
4	If blue colonies are seen: <ul style="list-style-type: none"> • Record observations in the LIS • Perform catalase test: 	
	IF	THEN
	Catalase POSITIVE	<ul style="list-style-type: none"> • Record observations in the LIS • Workup complete, VRE not isolated • Re-incubate plate for additional 38 to 48 hours if required
	Catalase NEGATIVE	<ul style="list-style-type: none"> • Record observations in the LIS • Subculture to BA plate
	From the BA sub plate: <ul style="list-style-type: none"> • Perform gram stain 	
	IF	THEN
	Gram stain NOT GPC	<ul style="list-style-type: none"> • Record observations in the LIS • Workup complete, VRE not isolated
	Gram stain GRAM POSITIVE COCCI	<ul style="list-style-type: none"> • Perform GPI • Set up vancomycin E-test
5	If pink colonies are seen: <ul style="list-style-type: none"> • Record observations in the LIS • Subculture colonies to BA sub plate • From the BA sub plate, perform gram stain 	
	IF	THEN
	Gram stain NOT GPC	<ul style="list-style-type: none"> • Record observations in the LIS • Workup complete, VRE not isolated
	Gram stain GRAM POSITIVE COCCI	<ul style="list-style-type: none"> • Perform GPI • Set up vancomycin E-test

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REPORTING INSTRUCTIONS:

IF	REPORT	
No Blue or Pink colonies	<ul style="list-style-type: none"> Report: "No Vancomycin Resistant Enterococci (VRE) isolated" 	
Blue and/or Pink colonies - Vitek ID: <i>E.gallinarum</i> and/or <i>E.casseliflavus</i>	<ul style="list-style-type: none"> Verify the organism ID Suppress GPI result in the isolates tab: <ul style="list-style-type: none"> ➤ Change the Isolate # to a letter ➤ Verify the result Enter and verify vancomycin E-test result Keep vancomycin E-test result suppressed Report: "No Vancomycin Resistant Enterococci (VRE) isolated" 	
Blue and/or Pink colonies - Vitek ID: <i>E.faecalis</i> and/or <i>E.faecium</i>	IF	REPORT
	Vancomycin E-test MIC = <4 µg/mL	<ul style="list-style-type: none"> Verify the organism ID Suppress GPI result in the isolates tab: <ul style="list-style-type: none"> ➤ Change the Isolate # to a letter ➤ Verify the result Enter and verify vancomycin E-test result Keep vancomycin E-test result suppressed Report: "No Vancomycin Resistant Enterococci (VRE) isolated"
Blue and/or Pink colonies - Vitek ID: <i>E.faecalis</i> and/or <i>E.faecium</i>	IF	REPORT
	Vancomycin E-test MIC = 4 µg/mL	<ul style="list-style-type: none"> Re-incubate vancomycin E-test <u>If after 48 hours MIC is still 4:</u> Verify the organism ID Suppress GPI result in the isolates tab: <ul style="list-style-type: none"> ➤ Change the isolate # to a letter ➤ Verify the result Enter and verify vancomycin E-test result Keep vancomycin E-test result suppressed Report: "No Vancomycin Resistant Enterococci (VRE) isolated"

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

		<p><u>If after 48 hours MIC \geq 8 μg/mL:</u></p> <ul style="list-style-type: none"> • Verify the organism ID • List quantitation as "Isolated" • Enter and verify vancomycin E-test result • Keep vancomycin E-test result suppressed • Report organism with isolate comment &VRE1 • Add test ?REFD and finalize with "." • In order entry, copy report to OCPHO (HPU1) • In order entry, copy report to appropriate IPAC ward • In order entry, add ESO code "VRE" • Freeze and record in isolate log • Forward isolate to <i>DynaLIFE</i> for vancomycin gene testing
<p>Blue and/or Pink colonies - Vitek ID: <i>E. faecalis</i> and/or <i>E. faecium</i></p>	<p style="text-align: center;">IF</p> <p>Vancomycin E-test MIC = 8-16 μg/mL</p>	<p style="text-align: center;">REPORT</p> <ul style="list-style-type: none"> • Repeat ID from vancomycin E-test plate • Verify the organism ID • List quantitation as "Isolated" • Enter and verify vancomycin E-test result • Keep vancomycin E-test result suppressed • Report organism with isolate comment &VRE2 • Add test ?REFD and finalize with "." • In order entry, copy report to OCPHO (HPU1) • In order entry, copy report to appropriate IPAC ward • In order entry, add ESO code "VRE" • Freeze and record in isolate log • Forward isolate to <i>DynaLIFE</i> for vancomycin gene testing

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

	IF	REPORT
Blue and/or Pink colonies - Vitek ID: <i>E. faecalis</i> and/or <i>E. faecium</i>	Vancomycin E-test MIC = ≥32 µg/mL	<ul style="list-style-type: none"> Verify the organism ID List quantitation as "Isolated" Enter and verify vancomycin E-test result Keep vancomycin E-test result suppressed The following isolate comment will be added: &VRE In order entry, copy report to OCPHO (HPU1) In order entry, copy report to appropriate IPAC ward In order entry, add ESO code "VRE" Freeze and record in isolate log

NOTE: STH IPAC ward is **SIPAC**. IRH IPAC ward is **IIPAC**

LIMITATIONS:

1. Organisms with atypical enzyme patterns may give anomalous results. The growth requirements of certain VRE can lead to their partial or total inhibition in culture.
2. Faecal specimens may cause some localized discolourization in the primary area of inoculation and should not be confused with a true chromogenic reaction wherein coloured colonies are visible. Interpret the colour of the isolate on well isolated colonies.
3. Strains of *E. faecalis* or *E. faecium* with intermediate resistance to vancomycin are infrequently encountered and may yield positive results.
4. Tightly clustered colonies of *Enterococcus gallinarum* and *Enterococcus casseliflavus* could appear as gray, slightly blue, or pink.
5. *Leuconostoc*, *Pediococcus* and *Lactobacillus* are intrinsically resistant to vancomycin, are inhibited or appear as colourless pinpoint colonies (blue or pink with tightly clustered colonies).
6. In some vancomycin-resistant *Enterococcus faecalis* strains such as ATCC 51299, which display low-level vancomycin resistance (MIC<32mg/L), the colonies may not develop the characteristic blue colour until 28h. If white colonies appear at 24h, a second reading must be performed at 28h.
7. Some Gram-negative bacteria may grow on VRESelect, but the size and colour of colonies enable differentiation from VREfm or VREfs.

CROSS-REFERENCES:

- MIC36400-Referral of Category B Specimens to DynaLIFE
- MIC60040-Culture Media Quality Control

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REFERENCES:

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4thed.) Washington, D.C.: ASM Press
2. Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology, 11th edition*. Washington, D.C: ASM Press
3. Bio-Rad. (2016/03). *VRESelect* package insert

APPROVAL:

Date

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	26 Apr 17	Initial Release	L. Steven
2.0	30 Nov 18	Updated to include new Vitek 2 instrument and two specimens per plate	L. Steven
3.0	30 Dec 21	Procedure reviewed and added to NTHSSA policy template	L. Steven
4.0	31 Aug 22	Updated to reflect new VRE agar VRESelect	L. Steven

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.