

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC32200 – Nasal Culture	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s):	
Effective Date:	Next Review Date:
Issuing Authority: Director of Health Services	Date Approved:
Accreditation Canada Applicable Standard: N/A	

GUIDING PRINCIPLE:

Nasal swabs are performed to identify nasal carriers of *Staphylococcus aureus*.

PURPOSE/RATIONALE:

This standard operating procedure describes the screening for *Staphylococcus aureus* in nasal specimens.

SCOPE/APPLICABILITY:

This procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for nasal culture.

SAMPLE INFORMATION:

Type	Swab <ul style="list-style-type: none">• Amie's with or without charcoal
Source	Nose
Stability	If the sample is received in the laboratory and processed greater than 48 hours from collection: <ul style="list-style-type: none">• Add specimen quality comment: "Delayed transport may adversely affect pathogen recovery"
Storage Requirements	Room temperature
Criteria for rejection	<ol style="list-style-type: none">1. Unlabeled/mislabeled swabs2. Specimen container label does not match patient identification on requisition3. Duplicate specimens obtained with same collection method within 24 hours

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REAGENTS and/or MEDIA:

- Blood agar (BA) and Mueller Hinton agar (MHP)
- Identification reagents: catalase, Staph latex test and cefoxitin antibiotic disks

SUPPLIES:

- Disposable inoculation needles
- Wooden sticks

EQUIPMENT:

- Biosafety cabinet
- 35° CO₂ incubator

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures:

- Ensure that appropriate hand hygiene practices be used
- Lab gown must be worn when performing activities with potential pathogens
- Gloves must be worn when direct skin contact with infected materials is unavoidable
- Eye protection must be used when there is a known or potential risk of exposure of splashes
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC)
- The use of needles, syringes and other sharp objects should be strictly limited

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to Test Manual for reagent quality control procedures

PROCEDURE INSTRUCTIONS:

Step	Action
Processing specimens for nasal culture	
1	In the biosafety cabinet: <ul style="list-style-type: none">• Inoculate BA with the swab• Ensure all surfaces of the swab make contact with the agar• Streak for isolated growth using a disposable inoculation needle
2	Incubate the media: <ul style="list-style-type: none">• Place BA in the CO₂ incubator

INTERPRETATION OF RESULTS:

Step	Action	
1	<ul style="list-style-type: none"> Observe BA plate at 24 hours and 48 hours Examine for colonies resembling <i>Staphylococcus aureus</i> 	
2	If no <i>S.aureus</i> colonies are seen at 24 hours: <ul style="list-style-type: none"> Record observations in the LIS Re-incubate plate in CO₂ incubator on the "Old wound culture" shelf 	
3	If no <i>S.aureus</i> colonies are seen at 48 hours: <ul style="list-style-type: none"> Record observations in the LIS Workup complete <i>Staphylococcus aureus</i> not isolated 	
4	If <i>S.aureus</i> colonies are seen: <ul style="list-style-type: none"> Record observations in the LIS Subculture colonies to BA plate if no isolated colonies are present If isolated colonies are present, perform Staph latex test 	
5	IF	THEN
	Staph latex test NEGATIVE	<ul style="list-style-type: none"> Record observations in the LIS Workup complete <i>S.aureus</i> not isolated
	Staph latex test POSITIVE	<ul style="list-style-type: none"> Record observations in the LIS <i>S.aureus</i> isolated Perform cefoxitin disk diffusion test
	IF	THEN
	Cefoxitin screen SENSITIVE	<ul style="list-style-type: none"> Record observations in the LIS Methicillin sensitive <i>S.aureus</i> isolated
	Cefoxitin screen RESISTANT	<ul style="list-style-type: none"> Record observations in the LIS Methicillin resistant <i>S.aureus</i> isolated

REPORTING INSTRUCTIONS:

IF	REPORT
<i>Staphylococcus aureus</i> not isolated	<ul style="list-style-type: none"> Report: "No Staphylococcus aureus isolated"
Methicillin sensitive <i>Staphylococcus aureus</i> isolated	<ul style="list-style-type: none"> Add organism: "Staphylococcus aureus" List quantification as "Isolated" Report organism with isolate comment &MSSA
Methicillin resistant <i>Staphylococcus aureus</i> isolated	<ul style="list-style-type: none"> Add organism: "Staphylococcus aureus" List quantification as "Isolated" Report organism with isolate comment &cx01 In order entry, copy report to OCPHO (HPU1) Check the home address of the patient. If from Nunavut: Copy report to the applicable NU CPHO In order entry, copy report to appropriate IPAC ward if ER or In-patient In order entry add ESO code "MRSA"

NOTE: STH IPAC ward is **SIPAC**. IRH IPAC ward is **IIPAC**

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REFERENCES:

1. Leber, A. (2016). *Clinical microbiology procedures handbook*. (4thed.) Washington, D.C.: ASM Press
2. Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology, 11th edition*. Washington, D.C: ASM Press

APPROVAL:

Date

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	03 Mar 19	Initial Release	L. Steven
2.0	22 Feb 21	Procedure reviewed and added to NTHSSA policy template	L. Steven
3.0	27 Feb 23	Procedure reviewed	L. Steven

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.