



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Microbiology Planting Job Aid

Distribution:
 Microbiology Specimen Processing Manual

Document Number:
 MIC10140

Date Issued:
 December 10, 2019

Uncontrolled When Printed

SPECIMEN TYPE	GRAM?	MEDIA	INCUBATION
Bacterial vaginosis	Yes	Only slide, no media	Only slide, no media
Blood culture: <ul style="list-style-type: none"> Positive bottles >24 hour bottles 	Yes	<u>Positive bottle:</u> BA, CHO, MAC, BRU <u>>24 hour bottle:</u> BA, CHO, BRU	BA and CHO = CO ₂ MAC = O ₂ BRU = AnO ₂
Body fluid culture	Yes	BA, CHO, MAC, BRU, THIO	BA and CHO = CO ₂ MAC and THIO = O ₂ BRU = ANO ₂
CSF culture	Yes	BA, CHO, MAC THIO if from shunt	BA and CHO = CO ₂ MAC (and THIO if indicated) = O ₂
Ear culture	Yes	BA, CHO, CNA, MAC	BA, CHO and CNA = CO ₂ MAC = O ₂
Eye culture-Superficial	Yes	BA, CHO	BA and CHO = CO ₂
Eye culture-Deep	Yes	BA, CHO, MAC, BRU, THIO	BA and CHO = CO ₂ MAC and THIO = O ₂ BRU = AnO ₂
Genital culture-Upper	Yes	BA, CHO, TM, MAC, BRU, THIO	BA, CHO and TM = CO ₂ MAC and THIO = O ₂ BRU = AnO ₂
Genital culture-Lower	Yes	BA, CHO, TM, MAC	BA, CHO and TM = CO ₂ MAC = O ₂
Gonorrhoeae culture	Male urethra	CHO, TM	CHO and TM = CO ₂
GBS screen	No	LIM, GBS	LIM = CO ₂
IUD culture	No	BRU, THIO	BRU = AnO ₂ THIO = O ₂
MRSA screen	No	MRS	MRS = O ₂
Nose culture	No	BA	BA = CO ₂

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Name:
Microbiology Planting Job Aid

Distribution:
Microbiology Specimen Processing Manual

Document Number:
MIC10140

Date Issued:
December 10, 2019

Uncontrolled When Printed

SPECIMEN TYPE	GRAM?	MEDIA	INCUBATION
Oral culture	No	SAB	SAB = room temperature
Respiratory culture	Yes	BA, CHO, MAC	BA and CHO = CO ₂ MAC = O ₂
Throat culture	No	BA	BA = AnO ₂
Tip culture	No	BA, MAC	BA = CO ₂ MAC = O ₂
Urine culture	No	URI	URI = O ₂
VRE screen	No	VRE	VRE = O ₂
Wound culture-Superficial	Yes	BA, CNA, MAC	BA and CNA = CO ₂ MAC = O ₂
Wound culture-Deep	Yes	BA, CHO, CNA, MAC, BRU, KV	BA, CHO and CNA = CO ₂ MAC = O ₂ BRU and KV = AnO ₂

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC10140-Microbiology Planting Job Aid

PRINT DATE: 1/11/2022 11:03 AM