

| PROGRAM Standard Operating Procedure – Laboratory Services | |
|---|-------------------|
| Title: MIC71000 – BACTEC FX Instrument | Policy Number: |
| Program Name: Laboratory Services | |
| Applicable Domain: Lab, DI and Pharmacy Services | |
| Additional Domain(s): NA | |
| Effective Date: | Next Review Date: |
| Issuing Authority: Director, Laboratory and Diagnostic Imaging Services | Date Approved: |
| Accreditation Canada Applicable Standard: NA | |

Uncontrolled When Printed

GUIDING PRINCIPLE:

The BACTEC FX Instrument is designed for the rapid detection of bacteria in clinical specimens. Blood cultures are collected from patients with suspected sepsis or bacteremia. Although primarily directed towards the processing of blood cultures, occasionally other specimen types (sterile fluids, abscess material, bone marrow, etc.) are received in blood culture vials. These vials may be processed in the same manner as blood cultures.

The BACTEC FX Instrument continuously monitors routine blood cultures for evidence of growth for 5 days. Negative results are auto-verified as follows:

- No growth after 48 hours (preliminary)
- No growth after 5 days (final)

PURPOSE/RATIONALE:

This standard operating procedure describes the BACTEC FX Instrument and its components.

SCOPE/APPLICABILITY:

This standard operating procedure applies to Medical Laboratory Technologists (MLTs) and Medical Laboratory Assistants (MLAs) processing specimens using the BACTEC FX Instrument.

SAMPLE INFORMATION:

| | |
|---------------|------------------------|
| Type | Blood culture vial |
| Source | Blood or sterile fluid |

REAGENTS and/or MEDIA:

| | |
|-------------------------------|---|
| Type | <ul style="list-style-type: none"> • BACTEC Plus Aerobic Culture Vials • BACTEC Lytic Anaerobic Culture Vials • BACTEC Peds Plus Culture Vials |
| Stability | Stable until date of expiration indicated on vial |
| Storage Requirements | Vial storage before blood collection: <ul style="list-style-type: none"> • Room temperature |
| Criteria for rejection | Do not use if: <ul style="list-style-type: none"> • The expiration date has passed • There are other signs of deterioration |

EQUIPMENT:

- BACTEC FX Instrument
- EpiCenter computer

ENVIRONMENTAL CONTROLS:

- Operating temperature: 18°C to 30°C
- Relative humidity: 25% to 80%, non-condensing

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures:

- Ensure that appropriate hand hygiene practices be used
- Lab gown must be worn when performing activities with potential pathogens
- Gloves must be worn when direct skin contact with infected materials is unavoidable
- Eye protection must be used when there is a known or potential risk of exposure of splashes
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC)
- The use of needles, syringes and other sharp objects should be strictly limited

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

PROCEDURE INSTRUCTIONS:


| Step | Action |
|------------------------------|--|
| Instrument components | |
| 1 | BACTEC FX Instrument: <ul style="list-style-type: none"> • The BACTEC FX Instrument is composed of 2 drawers: A and B • Drawers are divided into columns (numbered 1 to 10 from left to right) and rows of stations (lettered A to K excluding I) • There are a total of 100 stations in each drawer |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

2

Status display:

- The Status display is the main display shown when no other operation has been initiated or is in progress:




- Once a drawer is opened, you can initiate the major instrument activities from the Status display. Vial entry, remove positive vials, remove negative vials and identify anonymous vials can be initiated for any drawer in the instrument

3

System indicators:

- The system indicators are located on the front-center of the instrument:



- The system indicators inform you of various states in the instrument:

| Indicator Colour | State | Meaning |
|--|-------|--|
| Yellow=Light in unison for instrument | On | <ul style="list-style-type: none"> System Alert Indicator remains on until the condition is corrected/addressed |
| Green=One for each for left and right drawers | On | <ul style="list-style-type: none"> Negative vial Indicator remains lit until all negative vials are removed through the Remove Negative Vials activity |
| Red=One for each left and right drawers | On | <ul style="list-style-type: none"> Positive vial Indicator remains lit until all positive vials are removed through the Remove Positive Vials activity |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

4

Station indicators:

- Each station has a set of LED indicators that inform you of the stations or vials status
- The status indicator is located above each station:



- The color (red, green or yellow) and state (on, off or flashing) indicate the vial status:

| Colour | State | Meaning |
|------------|----------|-------------------------|
| Red | Flashing | Positive vial |
| Green | Flashing | Negative vial |
| Yellow | Flashing | Anonymous vial |
| Red Yellow | Flashing | Positive anonymous vial |
| Green | On | Available station |


5







Vial statuses:

- From the Status display, select the Drawer View button:




- The status of vials is displayed on the Drawer View Display:




| | <ul style="list-style-type: none"> The following Station statuses are shown: | | | |
|-----------|---|--|--|-----------------------------------|
| | Status | Icon | Meaning | How indicated |
| | Available |  | There is no vial in the station | Station indicator: GREEN |
| | Blocked |  | User has manually blocked the station | Station indicator: OFF |
| | Negative |  | Vial completed protocol with no evidence of positivity | Station indicator: FLASHING GREEN |
| | Ongoing |  | Vial is in the instrument and is in protocol | Station indicator: OFF |
| | Positive |  | Instrument has detected evidence of microbial growth | Station indicator: FLASHING RED |
| Anonymous |  | Vial was physically placed in instrument without its barcode sequence number being scanned | Station indicator: FLASHING YELLOW | |
| 6 | Audible tones and alarms: <ul style="list-style-type: none"> Numerous different sounds are generated by the BACTEC FX Instrument: | | | |
| | Type | Example | Sound | |
| | Activity complete | All negative vials were removed | High pitched tone repeated 3 times | |
| | Activity error | Did not scan accession barcode after scanning sequence barcode | Single high beep | |
| | Anonymous | Anonymous vial entered | Short buzz sound | |
| | Drawer ajar | The door is not closed | 2 tones, high then low frequency, repeating | |
| | Positive vial | A positive vial is detected | Pulsing fading sound, repeated | |
| | Vial entry | A vial was entered | High pitch chirp sound | |

| Step | Action |
|---|--|
| Loading vials into the BACTEC FX | |
| 1 | Open the drawer. Ensure that the Status screen is displayed, and the barcode scanner turns on. |
| 2 | Scan the vial sequence barcode (vial barcode) and the accession barcode (LIS label). |
| 3 | Place the vial into any available slot (solid green light) in the instrument. |
| 4 | Scan and place any other vials that need to be loaded. |
| 5 | Close drawer when finished. |


| | |
|----------|---|
| 6 | <p>Delayed entry of vials may lead to delayed results. Best practices dictate that vials should be placed into the BACTEC FX as soon as possible after collection:</p> <ul style="list-style-type: none"> All vials received >24 hours after collection need to be sub-cultured to ensure growth is not present Refer to MIC10100-Microbiology Specimen Processing |
| 7 | <p>Vials can be placed into available (GREEN indicator) stations without being scanned into the instrument. Vials that are not scanned are called "anonymous" vials.</p> |




| Step | Action |
|------------------------------------|--|
| Identifying anonymous vials | |
| 1 | <p>The Status display will indicate which drawer houses the anonymous vial:</p>  |
| 2 | <p>Open the drawer and locate station(s) with a flashing Yellow LED and remove the vial. NOTE: If the station flashes YELLOW and RED, then the vial is also POSITIVE</p> |
| 3 | <p>The ID Anonymous display appears, and the barcode scanner turns on.</p> |
| 4 | <p>Ensure the vial has an accession barcode (LIS label). If not, accession the sample and place the barcode on the vial.</p> |
| 5 | <p>Scan the vial sequence barcode (vial barcode) and the accession barcode (LIS label).</p> |
| 6 | <p>Place the vial in the FLASHING GREEN station (station from which the vial was pulled).</p> |
| 7 | <p>Repeat with any additional anonymous vials.</p> |

| Step | Action |
|--------------------------------|---|
| Removing positive vials | |
| 1 | <p>The Status display will indicate which drawer houses the positive vial:</p>  |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.




| | |
|----------|---|
| 2 | Open the drawer and locate the station(s) with a flashing Red LED and remove the vial(s). |
| 3 | The Positive Removal display appears. Scan the vial sequence barcode (vial barcode). NOTE: You must scan each positive vial that you remove in order for the instrument to re-light positive stations |
| 4 | When all positive vials are removed from the drawer, the Activity Complete tone sounds. |

| Step | Action |
|--------------------------------|---|
| Removing negative vials | |
| 1 | The Status display will indicate which drawer houses the negative vial:  |
| 2 | Open the drawer and locate the station(s) with a flashing Green LED and remove the vial(s). These vials do not have to be scanned. |
| 3 | When all negative vials are removed from the drawer, the Activity Complete tone sounds. |

| Step | Action |
|---|--|
| Extending incubation time of vials | |
| 1 | On the BACTEC Status display, select the Drawer View button:  |
| 2 | Select the desired station and select OK. |
| 3 | Select Modify and the Modify Protocol box will be displayed:  |
| 4 | Change the protocol length to 10 days and select OK. Select Save to save the changes. |
| 5 | In the LIS, under Results Entry, enter the order number to access the plate log/test comments <ul style="list-style-type: none"> In the plate log (Media Comments), look for the EXT media ID:  |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

- Double click in the resulted box (R) for the EXT Media ID so that a red check mark appears:
- | # | Media ID | Media Comment | R |
|---|----------|---|---|
| 1 | EXT | TO EXTEND PROTOCOL:RESULT THIS MEDIA-LOOK FOR RED CHECKMARK | ✓ |
| 2 | ANA | | |
| 3 | R&R | | |
- NOTE:** If a set of vials has been collected, both vials will need to be modified
- This alerts the LIS to stop the 5 day reporting
 - Save the order to save the changes made
- 6** If culture is negative:
- A 48 hour negative preliminary report will be automatically released
 - A 10 day no growth final report will automatically be released
- NOTE:** Manual reporting negative vials on extended protocol is not required

| Step | Action |
|-------------------------------|--|
| Resolving system alert | |
| 1 | A system alert is indicated by a yellow LED indicator on both drawers of the instrument. This alert usually indicates a power failure or communication interruption. |
| 2 | <ul style="list-style-type: none"> • Check the instrument for an error message • On the Status display, select the Systems Alert button to view the Alert List:  <ul style="list-style-type: none"> • Power interruptions will display the following alerts: <ul style="list-style-type: none"> ➢ Reboot Reason: Power fail ➢ The instrument has lost connectivity to the server ➢ EpiCenter Communications failure • Refer to the BACTEC FX Instrument User Manual for alert descriptions |
| 3 | <ul style="list-style-type: none"> • Log into EpiCenter computer: <ul style="list-style-type: none"> ➢ Log into Windows. When power goes out, Windows will re-boot and require re-login ➢ Icon will display:  Communication should RESYNC after logging into Windows (should take about one minute after logging in)  |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.


- A System Message window should pop-up displaying the errors. Click the x button to close or delete them
- The yellow system indicator lights on BACTEC FX door will stop glowing


| Step | Action |
|------|--------|
|------|--------|


| | |
|-----------------------------|--|
| Disassociating vials | |
|-----------------------------|--|

| | |
|----------|--|
| 1 | If a vial record contains an accession number, it is considered associated to that accession. The disassociate function enables you to break the link between a vial and an accession number. This can be useful when troubleshooting sample errors. |
|----------|--|

| | |
|----------|---|
| 2 | Remove the vial(s) from the BACTEC FX Instrument. |
|----------|---|

| | |
|----------|---|
| 3 | From the Status display, select the Culture button: <div style="text-align: center; margin-top: 10px;">  </div> |
|----------|---|




| | |
|----------|--|
| 4 | <p>The Culture-Patient display appears:</p> <div style="text-align: center; margin-top: 10px;">  </div> |
|----------|--|

| | |
|----------|--|
| 5 | <p>Select the Vial tab to access the Culture-Vial display:</p> <div style="text-align: center; margin-top: 10px;">  </div> |
|----------|--|

| | |
|----------|--|
| 6 | Scan the vial sequence barcode (vial barcode). |
|----------|--|

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

| | |
|-----------|--|
| 7 | Select the Disassociate button to disassociate the vial from the accession number. |
| 8 | Return to the Status display. |
| 9 | Open the drawer and proceed to load the vial into the instrument. |
| 10 | Scan the vial sequence barcode (vial barcode) and then scan the accession barcode (LIS label). |
| 11 | Place the vial in any available slot in the instrument. |

| Step | Action |
|--------------------------|---|
| Associating vials | |
| 1 | If a vial has been loaded into the instrument without an accession number, the vial needs to be associated with the accession number. |
| 2 | Remove the vial(s) from the BACTEC FX Instrument. |
| 3 | From the Status display, select the Culture button:  |
| 4 | The Culture-Patient display appears:  |
| 5 | Select the Specimen tab to access the Culture-Specimen display:  |

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

| | |
|----------|---|
| 6 | In the accession field, scan the accession barcode (LIS label). |
| 7 | Scan the vial sequence barcode (vial barcode) you want to attach. |
| 8 | Select the Save button to save the association. |

CROSS-REFERENCES:

- MIC10100-Microbiology Specimen Processing
- MIC20500-Gram stain resulting in LIS-Blood Cultures
- MIC60010-Microbiology Quality Control

REFERENCES:

1. Becton Dickinson and Company. (2016-12). *BD BACTEC FX Instrument User Manual*, 8005110(07)

APPROVAL:

Date

Director, Laboratory and Diagnostic Imaging Services

REVISION HISTORY:

| REVISION | DATE | Description of Change | REQUESTED BY |
|-----------------|-------------|---|---------------------|
| 1.0 | 06 Nov 17 | Initial Release | L. Steven |
| 2.0 | 26 Mar 19 | Updated to reflect addition of disassociating and associating vials | L. Steven |
| 3.0 | 16 Aug 21 | Procedure reviewed and added to NTHSSA policy template | L. Steven |
| 4.0 | 01 Oct 24 | Procedure reviewed | L. Steven |
| | | | |
| | | | |
| | | | |