

## Staphylococcus saprophyticus (LTR81529)

Edit Approved By: Van der Walt, Peet (09/27/2023)

Revision: 5.00

**Organism**                      **Staphylococcus saprophyticus**

---

**Clinical**                      This organism is associated with urinary tract infections, typically in sexually active younger females but also rarely in males, older females and children (often in association with indwelling urinary catheters or obstruction). It has been implicated in nongonococcal urethritis in males, prostatitis, and may be sexually transmitted. It has also been associated with bacteremia, septicemia, nosocomial pneumonia, and wound infections.

---

**Usual susceptibility pattern**                      *S. saprophyticus* is usually susceptible to amoxicillin/clavulanate, first generation cephalosporins, quinolones, nitrofurantoin, and tetracyclines. It is generally susceptible to TMP-SMX although some resistance has been described. This organisms is resistant to fosfomycin.

---

**Susceptibility method**                      VITEK2. Additional tests include disc diffusion and Etest method.

Disc diffusion		Mueller-Hinton agar incubated in ambient air at 35°C for 16-18 hours
Cefoxitin Screen disc		Mueller-Hinton agar incubated in ambient air at 35°C for 24 hours. Use 0.5 McFarland suspension in saline.
Etest	Oxacillin/ Cloxacillin	Mueller-Hinton agar with 2% NaCl incubated in ambient air at 35°C for 48 hours. Use 1.0 McFarland suspension in saline.
	Vancomycin	Mueller-Hinton agar incubated in ambient air at 35°C for 24 hours. Use 0.5 McFarland suspension in saline.
	Other	Mueller-Hinton agar incubated in ambient air at 35°C for 16-20 hours. Use 0.5 McFarland suspension in saline.

# Staphylococcus saprophyticus, Continued

## Susceptibility reporting

	Blood	Sterile Body Site	Urine (non-pregnant) See Note	Urine (pregnant) See Note	Comments
Ampicillin			✓	✓	Disc diffusion using <b>2 µg disc</b> Report as R if oxacillin R <b>See Special Considerations</b>
Cefazolin	✓	✓	✓	✓	Report same as ox/clox
Ciprofloxacin			✓		Do not report in patients < 18 y
Doxycycline		✓	✓		If tetra S - report doxy S If tetra I/R - do doxy disc Do not report if patient <8 y
Oxacillin/ Cloxacillin	✓	✓	✓	✓	<b>Refer to Staphylococcus Oxacillin Reporting Flowchart (Doc ID: MIC - 37934)</b>
Nitrofurantoin			✓	✓	
TMP-SMX		✓	✓	✓	Do not report if patient <2 months
Vancomycin	2	2			2 <sup>nd</sup> line if ox/clox R If vancomycin ≥4 µg/mL see <b>Special Considerations</b>

### Note

<b>Urine specimens</b>	Susceptibility testing not routinely performed, as urinary tract infections (UTI) due to this organism respond to concentrations achieved in urine of antibiotics routinely used to treat UTIs.	
	<b>Susceptibility testing may be performed at physician request or if patient pregnant as therapeutic options are limited.</b>	
	<b>IF ...</b>	<b>THEN ...</b>
<ul style="list-style-type: none"> <li>pregnant</li> <li>physician request</li> </ul>	Perform susceptibility testing according to reporting chart	
None of the above	Add comment: "Susceptibility testing is not routinely performed on Staphylococcus saprophyticus. If antibiotic therapy is indicated, they are normally susceptible to amoxicillin-clavulanate, cephalexin, trimethoprim-sulfamethoxazole, nitrofurantoin and fluoroquinolones." <b>&amp;A366</b>	

# Staphylococcus saprophyticus, Continued

## Special considerations

<p><u>Ampicillin:</u></p>	<p>Use ampicillin <b>2 µg discs</b>. Interpret using EUCAST disc diffusion breakpoints.</p> <p><b>EUCAST breakpoints:</b></p> <table border="1" data-bbox="399 386 854 506"> <thead> <tr> <th>Zone</th> <th>Interpretation</th> </tr> </thead> <tbody> <tr> <td>≤17</td> <td>R</td> </tr> <tr> <td>≥18</td> <td>S</td> </tr> </tbody> </table> <p>Add comment: "Interpretation is based upon EUCAST breakpoints." <b>#amp3</b></p>	Zone	Interpretation	≤17	R	≥18	S		
Zone	Interpretation								
≤17	R								
≥18	S								
<p><u>Vancomycin:</u></p>	<p>Isolates with VITEK2 MIC <b>≥4 µg/mL</b>, confirm MIC by Etest and consult microbiologist.</p> <table border="1" data-bbox="399 653 1373 1686"> <thead> <tr> <th>IF vancomycin is...</th> <th>THEN....</th> </tr> </thead> <tbody> <tr> <td data-bbox="399 653 716 1050"> <p><b>4 µg/mL</b> (confirmed by Etest)</p> </td> <td data-bbox="716 653 1373 1050"> <ul style="list-style-type: none"> <li>The clinical failure rate of vancomycin may be significant.</li> <li>Consult Technical Supervisor</li> <li>Add comment: "This isolate tests at the upper limit of susceptibility to vancomycin. Careful follow up to assess clinical response is required, or an alternate agent should be considered. Expert consultation is suggested." <b>#va04</b></li> </ul> </td> </tr> <tr> <td data-bbox="399 1050 716 1251"> <p><b>8-16 µg/mL</b> (confirmed by Etest)</p> </td> <td data-bbox="716 1050 1373 1251"> <ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as I</li> <li>Add comment: "This isolate exhibits resistance to vancomycin." <b>#va11</b></li> <li>Notify Infection Control &amp; MOH</li> </ul> </td> </tr> <tr> <td data-bbox="399 1251 716 1686"> <p><b>≥ 32 µg/mL</b> (confirmed by Etest)</p> </td> <td data-bbox="716 1251 1373 1686"> <ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as R</li> <li>Add comments: "Preliminary tests indicate this organism may be resistant to vancomycin" <b>#va12</b> "Referred to Public Health Laboratory, Alberta Precision Laboratories, for Van gene testing." <b>#va13</b></li> <li>Notify Infection Control &amp; MOH</li> <li>Send to reference laboratory for Van gene testing.</li> </ul> </td> </tr> </tbody> </table>	IF vancomycin is...	THEN....	<p><b>4 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>The clinical failure rate of vancomycin may be significant.</li> <li>Consult Technical Supervisor</li> <li>Add comment: "This isolate tests at the upper limit of susceptibility to vancomycin. Careful follow up to assess clinical response is required, or an alternate agent should be considered. Expert consultation is suggested." <b>#va04</b></li> </ul>	<p><b>8-16 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as I</li> <li>Add comment: "This isolate exhibits resistance to vancomycin." <b>#va11</b></li> <li>Notify Infection Control &amp; MOH</li> </ul>	<p><b>≥ 32 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as R</li> <li>Add comments: "Preliminary tests indicate this organism may be resistant to vancomycin" <b>#va12</b> "Referred to Public Health Laboratory, Alberta Precision Laboratories, for Van gene testing." <b>#va13</b></li> <li>Notify Infection Control &amp; MOH</li> <li>Send to reference laboratory for Van gene testing.</li> </ul>
IF vancomycin is...	THEN....								
<p><b>4 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>The clinical failure rate of vancomycin may be significant.</li> <li>Consult Technical Supervisor</li> <li>Add comment: "This isolate tests at the upper limit of susceptibility to vancomycin. Careful follow up to assess clinical response is required, or an alternate agent should be considered. Expert consultation is suggested." <b>#va04</b></li> </ul>								
<p><b>8-16 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as I</li> <li>Add comment: "This isolate exhibits resistance to vancomycin." <b>#va11</b></li> <li>Notify Infection Control &amp; MOH</li> </ul>								
<p><b>≥ 32 µg/mL</b> (confirmed by Etest)</p>	<ul style="list-style-type: none"> <li>Consult Technical Supervisor</li> <li>Report vancomycin as R</li> <li>Add comments: "Preliminary tests indicate this organism may be resistant to vancomycin" <b>#va12</b> "Referred to Public Health Laboratory, Alberta Precision Laboratories, for Van gene testing." <b>#va13</b></li> <li>Notify Infection Control &amp; MOH</li> <li>Send to reference laboratory for Van gene testing.</li> </ul>								

## Staphylococcus saprophyticus, Continued

**Interpretation** For Etest, report actual MIC result. For interpretation (S, I, or R) report according to the nearest higher doubling dilution (**Appendix 1**).

Use **CLSI** interpretive document for **Staphylococcus spp.** using Coagulase Negative Staphylococcal breakpoints.

For oxacillin and ceftiofloxacin – Use EUCAST breakpoints. Refer to *Staphylococcus Oxacillin Reporting Flowchart* (Doc ID: MIC - 37934)

For ampicillin – Refer to Special Considerations