

**Citrobacter freundii complex (LTR62253)**

Edit Approved By: Solomon, Natalia (05/03/2022)

Revision: 5.00

**Organism****Citrobacter freundii complex:**

- *C. braakii*
  - *C. freundii*
  - *C. werkmanii*
  - *C. youngae*
  - *C. murlinae*
- 

**Clinical**

These organisms may be found in water, soil, food and human/animal gastrointestinal tracts. They are associated with hospital-acquired infections, especially of the respiratory and urinary tracts. Sepsis with *Citrobacter* spp. is often polymicrobial and associated with high morbidity and mortality.

Organisms in this group are important opportunistic pathogens and tend to be more resistant to antibiotics.

---

**Usual  
susceptibility  
pattern**

These organisms produce a chromosomally mediated inducible cephalosporinase (AmpC) and are resistant to penicillins, and first/second generation cephalosporins. Although they often exhibit in vitro susceptibility to third generation cephalosporins, use of these agents clinically may result in selection of resistant strains. The beta-lactamase produced by these organisms is not inhibited by beta-lactamase inhibitors and as such, beta-lactam/beta-lactamase inhibitor combinations should not be reported. Although extended spectrum beta-lactamase (ESBL) may be found in these organisms, conventional ESBL testing is not reliable, due to interference by the chromosomal cephalosporinase. Testing of cefepime +/- clavulanic acid may detect an ESBL enzyme. Most strains are susceptible to aminoglycosides, TMP-SMX, quinolones, and nitrofurantoin.

---

**Susceptibility  
method**

VITEK2. Additional tests (Disc diffusion or Etest method) are performed using Mueller-Hinton agar incubated in ambient air at 35°C for 16-20 hours.

**Note:** For Etest use 0.5 McFarland suspension in saline.  
For mucoid strains use 1.0 McFarland.

## Citrobacter freundii complex, Continued

### Susceptibility reporting

	CSF/ Brain	Blood/Sterile Body Site/ Endovascular Catheter	Urine	Other	Comments
Amikacin		3	3	3	3 <sup>rd</sup> line if gent and tobra I/R Disc diffusion
Ampicillin	R	R	R	R	
Cefazolin		R	R	R	
Cefixime			R		
Ceftriaxone	R	R			
Ciprofloxacin		✓	✓	✓	Do not report in patients < 18 y
Doxycycline			2		2 <sup>nd</sup> line if cipro I/R For patients ≤17 y report 1 <sup>st</sup> line Disc diffusion If patient <8 y <b>See special considerations</b>
Ertapenem		✓	2	2	2 <sup>nd</sup> line if cipro or TMP-SMX I/R If S do not report in patients < 3 months
Gentamicin	*	✓**	✓	✓**	* Report only in neonates (< 1 month) <b>**See Special Considerations</b>
Imipenem *		✓	2	2	2 <sup>nd</sup> line if cipro or TMP-SMX I/R
Meropenem	✓	✓	2*	2*	2 <sup>nd</sup> line if cipro or TMP-SMX I/R * Report 1 <sup>st</sup> line in neonates (< 1 month)
Nitrofurantoin			✓		<b>Add comment:</b> For uncomplicated lower UTI only <b>#1</b>
TMP-SMX	*	✓	✓	✓	* Report only at physician request
Tobramycin		2*	2	2*	2 <sup>nd</sup> line if gent I/R <b>*See Special Considerations</b>

\* Do NOT report Imipenem from the VITEK

## Citrobacter freundii complex, Continued

### Special considerations

<p><u>Doxycycline:</u></p>	<p>If reporting doxycycline on patients &lt;8 years add the following comment:</p> <p>“Doxycycline can now be prescribed for children &lt;8y for short-course (&lt;21 d) therapy; OTHER tetracyclines are still contraindicated for this age group.” <b>(27664)</b></p>
<p><u>Gentamicin/tobramycin:</u></p>	<p>Organisms testing at upper limit of susceptibility (4µg/mL) may not achieve optimal pharmacokinetics/pharmacodynamics.</p> <p><b>For non-urine isolates:</b></p> <p>If MIC 4.0 µg/mL add comment:</p> <p>“This isolate tests at the upper limit of susceptibility for gentamicin. Clinical failure may occur despite in vitro susceptibility.” <b>#A312</b></p> <p style="text-align: center;"><b>or</b></p> <p>“This isolate tests at the upper limit of susceptibility for tobramycin. Clinical failure may occur despite in vitro susceptibility”. <b>#A313</b></p> <p style="text-align: center;"><b>or</b></p> <p>“This isolate tests at the upper limit of susceptibility for both gentamicin and tobramycin. Clinical failure may occur despite in vitro susceptibility”. <b>#A314</b></p>

### Interpretation

For Etest, report actual MIC result. For interpretation (S, I, or R) report according to the nearest higher doubling dilution (**Appendix 1**).

Use **CLSI** interpretive document for **Enterobacterales**.

**For Beta-lactam drugs** – Refer to Beta-lactam Resistance Detection Charts.

**For gentamicin and tobramycin** – Refer to Special Considerations