

Proteus mirabilis (LTR62809)

Edit Approved By: Solomon, Natalia (06/08/2022) Revision: 8.00

Organism

Proteus mirabilis

Clinical

This organism is found in a variety of environmental sources and may be part of the normal flora of the gastrointestinal tract. It is a cause of both community and nosocomial urinary tract infections, usually catheter associated infections or in chronically catheterized patients. It accounts for 90% of *Proteus spp.* infections. Infection with this organism indicates upper urinary tract involvement.

P. mirabilis has also been associated with post-operative wound infections, skin/soft tissue infections, and bacteremia (usually secondary to a urinary tract focus).

Usual susceptibility pattern

This organism does not produce any chromosomal beta-lactamase and is usually susceptible to ampicillin and cephalosporin. However, plasmid mediated beta-lactamases [penicillinase, cephalosporinase, extended spectrum beta-lactamases (ESBL) and carbapenemase (metalloenzyme, KPC or Class D OXA enzyme)] are increasingly being described. ESBL production may be more difficult to detect. Beta-lactam resistance may also be mediated by permeability mutations. *P. mirabilis* is usually susceptible to aminoglycosides, quinolones, TMP-SMX and carbapenems. (Exception: May exhibit decreased susceptibility to imipenem (decreased affinity to PBP2) that does not affect other carbapenems.) *P. mirabilis* is resistant to nitrofurantoin, tetracycline and colistin.

Note: Although generally susceptible to beta-lactam antibiotics, non-beta-lactam antibiotics appear to be more efficacious in treatment of urinary tract infections with this organism.

Susceptibility method

VITEK2. Additional tests (Disc diffusion or Etest method) are performed using Mueller-Hinton agar incubated in ambient air at 35°C for 16-20 hours.

Note: For Etest use 0.5 McFarland suspension in saline. For mucoid strains use 1.0 McFarland.

Proteus mirabilis, Continued

Susceptibility reporting

	CSF/ Brain	Blood/ Endo- vascular Catheter	Sterile Body Site	Urine	Other	Comments
Amikacin		3	3	3	3	3 rd line if gent and tobra I/R Disc diffusion
Amoxicillin/ Clavulanateoral				✓	✓	If Amp S and Amox-Clav I/R – Do Amp and AMC disc diffusion, consult tech 2 with results See Special Considerations
Amoxicillin/ ClavulanateIV		2*	2*		2*	2nd line if ampicillin I/R, cefazolin I/R and ceftriaxone S *Report 1st line if ampicillin I/R and ceftriaxone S and anaerobes, Enterococcus or S. aureus (MSSA) co-isolated Report same as AMC oral
Ampicillin	*	✓	✓	✓	✓	* Report only in neonates (<1 month)
Cefazolin		✓	✓	√ *	✓	*If MIC ≤16 Do not report (offer); ≥32 report as R. Refer to Beta-Lactam Resistance Detection Charts
Cefixime				✓		
Ceftriaxone	√ *	√ *	2*	2*	2*	Always report if I/R 2 nd line if cefazolin I/R *If patient < 1 mo - report cefotaxime instead of ceftriaxone using the same interpretation.
Cephalexin				✓		Add comment: For uncomplicated lower UTI only #CIX1
Ciprofloxacin		√	√	2*	✓	Do not report in patients < 18 y 2 nd line if cefixime and TMP-SMX I/R. Always report if I/R *See Special Considerations
Ertapenem		3	3	3	3	3 rd line if ceftriaxone or ceftazidime I/R If S do not report in patients < 3 months See Special Consideration s
Gentamicin	*	√ **	√ **	✓	√ **	* Report only in neonates (< 1 month) **See Special Considerations
Meropenem	2	3	3	3	3	2 nd or 3 rd line if ceftriaxone or ceftazidime I/R See Special Considerations
Nitrofurantoin				R		Add comment: For uncomplicated lower UTI only #f1
Piperacillin/ Tazobactam		3* ╬	3*		3*	3 rd line if AMC IV R and ceftriaxone S * Report if <i>P.aeruginosa</i> co-isolated and ceftriaxone S † For bloods report if I/R and ceftriaxone S See Special Considerations
TMP-SMX	*	✓	✓	✓	✓	*Report only at physician request
Tobramycin		2*	2*	2	2*	2 nd line if gent I/R *See Special Considerations

Proteus mirabilis, Continued

Special considerations

Amoxicillin/	A MIC of 8/4 µg/mL is at upper limit of susceptibility. This may be adequate to					
<u>Clavulanate</u>	achieve reasonable pharmacodynamics in urine but may not be optimal for non-					
oral:	urinary sites.					
	For all non-urinary sites if MIC 8/4 µg/mL and interpretation is S add comment: "This isolate tests at the upper limit of susceptibility for amoxicillin/clavulanate. Clinical failure may occur despite in vitro susceptibility". #A315					
	I/R VITEK 2 results are not recommended due to a card limitation. Perform an					
	alternate method prior to reporting of results.					
Ciprofloxacin:	For urine cultures add the following comment when not reporting ciprofloxacin					
<u>Cipionoxaciii.</u>	(patients ≥ 18 y):					
	"Ciprofloxacin is not routinely reported, given the potential for significant adverse events and increasing antimicrobial resistance." &3206					
Ertapenem /	If reporting ertapenem and meropenem add comment:					
Meropenem:	"Imipenem has intrinsically low activity against this organism." #A375					
Gentamicin/ tobramycin:	Organisms testing at upper limit of susceptibility (4µg/mL) may not achieve optimal pharmacokinetics/pharmacodynamics.					
	For non-urine isolates: If MIC 4.0 µg/mL add comment: "This isolate tests at the upper limit of susceptibility for gentamicin. Clinical failure may occur despite in vitro susceptibility".#A312					
	"This isolate tests at the upper limit of susceptibility for tobramycin. Clinical failure may occur despite in vitro susceptibility". #A313					
	Or "This is plate tooks at the upper limit of susceptibility for both gentamicin and					
	"This isolate tests at the upper limit of susceptibility for both gentamicin and tobramycin. Clinical failure may occur despite in vitro susceptibility". #A314					
Piperacillin/	This antibiotic is frequently used as empiric therapy for polymicrobial infections (i.e.					
tazobactam:	co-infections with <i>S. aureus, Enterococcus, Pseudomonas aeruginosa</i> and/or anaerobes), febrile neutropenia or sepsis syndromes.					
	Note: Do not report as S if ceftriaxone or ceftazidime I/R (≥2 μg/mL) as piperacillin/ tazobactam is not recommended if either an extended spectrum beta-lactamase (ESBL) and/or a cephalosporinase is present.					

Proteus mirabilis, Continued

Interpretation

For Etest, report actual MIC result. For interpretation (S, I, and R) report according to nearest higher doubling dilution (Appendix 1).

Use **CLSI** interpretive document for **Enterobacterales**.

For Beta-lactam drugs – Refer to Beta-lactam Resistance Detection Charts. For Amoxicillin/Clavulanate, gentamicin and tobramycin – Refer to Special Considerations