



Detection of Beta Lactam Resistance Based on Ertapenem, Imipenem and Meropenem Susceptibility

**Chart 1**





Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g/mL}$

**Chart 2A**


*Escherichia spp.*  
*Edwardsiella spp.*  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae complex*  
*Klebsiella variicola*  
*Proteus mirabilis*  
*Raoultella spp.*  
*Salmonella spp.*  
*Shigella spp.*

MAST AmpC and ESBL discs (D68C)

Interpret:

- Amox/clav AMC oral - R
- Amox/clav IV AMC IV - R - Do not report
- Ampicillin AM - R
- Cephalexin CEX - R (Urine Only) **Note 1**
- Cefazolin CZ - R
- Cefixime CFM - R (Urine Only)
- Ceftriaxone CRO - R
- Pip/taz TZP - Do not report

Note 1: Report for *E. coli*, *K. pneumo* and *P. mirabilis* only



**NOTES:**


For *E. coli* and *Shigella spp.*, it is not possible, with current phenotypic methods to differentiate between plasmid mediated AmpC cephalosporinase and hyperproduction of intrinsic chromosomal AmpC enzyme.

*Edwardsiella spp.*, *Klebsiella spp.*, *Proteus mirabilis* and *Salmonella spp.* do not produce chromosomal beta-lactamases. Hence, all beta-lactamase detection implies acquisition of a transmissible beta-lactamase.

Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g/mL}$

**Chart 2B**

*Citrobacter amalonaticus*  
*Citrobacter farmeri*  
*Citrobacter koseri*  
*Citrobacter sedlakii*  
*Klebsiella oxytoca*




**NOTE:** Hyperproduction of the class A enzymes (K1 in *K. oxytoca*, Cko in *Citrobacter koseri* and CdiA *Citrobacter amalonaticus*) typically results in resistance to aztreonam, ceftriaxone (cefotaxime less affected than ceftriaxone) and cefepime. Although ceftazidime usually tests susceptible its use is not recommended as it may select resistant mutants.

Ceftriaxone **and/or** Ceftazidime  $\geq 2 \mu\text{g/mL}$

**Chart 2C**

*Proteus hauseri*  
*Proteus penneri*  
*Proteus vulgaris*



These organisms produce an inducible class A cephalosporinase (cefuroximase) that results in resistance to penicillins, and 1<sup>st</sup>/2<sup>nd</sup> generation cephalosporins. They remain susceptible to beta-lactamase inhibitor combination drugs. Hyperproduction/derepression of this enzyme typically results in cefotaxime/ceftriaxone resistance (ceftazidime remains susceptible).

Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g/mL}$

## Chart 2D

Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriazone CRO - R  
Pip/taz TZP - Do not report

*Cedecea spp.*  
*Citrobacter freundii complex*  
*Enterobacter spp.*  
*Hafnia spp.*  
*Klebsiella aerogenes*  
*Morganella morganii*  
*Providencia spp.*  
*Serratia spp.* (other than *S. fonticola*)

Comment  
&1746

These organisms produce an inducible chromosomal Amp C cephalosporinase that typically results in resistance to penicillins and 1<sup>st</sup>/2<sup>nd</sup> generation cephalosporins, including cefoxitin (exceptions: *Serratia* spp, cefuroxime may test more resistant than cefoxitin; *Hafnia* and *Providencia* species, both may still test susceptible to cefoxitin and cefuroxime). Induction and/or derepression of the Amp C enzyme results in broad resistance including 3<sup>rd</sup> generation cephalosporins, especially cefotaxime and ceftriaxone. Cefepime, a 4<sup>th</sup> generation cephalosporin, is poorly hydrolyzed by the Amp C enzyme (typical MICs  $\leq 1\mu\text{g/mL}$ ).

**NOTE:** Resistance to cefepime implies either acquisition of an extended spectrum beta lactamase (ESBL) or an extended spectrum Amp C cephalosporinases (ESAC) resulting from further derepression of the chromosomal enzyme. Acquisition of an ESBL is common but difficult to detect by standard ESBL confirmatory tests due to interference from the chromosomal Amp C enzyme that is not inhibited by clavulanate. Cefepime ESBL discs may be useful in detecting ESBL production in the presence of an Amp C cephalosporinase.

It is not possible to detect an ESBL in the presence of an ESAC as the higher level of Amp C cephalosporinase interferes with the cefepime ESBL test. Nosocomial outbreaks of cefepime resistant organisms warrant molecular confirmation of resistance mechanism(s) (ESBL or ESAC).

Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g/mL}$

## Chart 2E

*Erwinia spp.*  
*Kluyvera spp.*  
*Rahnella spp.*

Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO - R

Comment

&1744

These organisms produce a low level chromosomal Class A extended spectrum beta-lactamase (ESBL) conferring resistance to penicillins, and 1<sup>st</sup>/ 2<sup>nd</sup> generation cephalosporins, but not cefoxitin. Like other ESBLs, resistance to 3<sup>rd</sup> generation cephalosporins may not be apparent in vitro, but should be implied. ESBL testing is not useful as all should test positive due to the chromosomal ESBL.

Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g/mL}$

## Chart 2F

*Cronobacter spp.*  
*Lelliottia amnigena*  
*Pluralibacter gergoviae*  
*Pantoea spp.*

Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO - R

Comment  
&1746

These organisms differ from *E. cloacae* and *K. aerogenes* in that the chromosomal Amp C cephalosporinase is produced at very low levels. Ampicillin, 1<sup>st</sup> generation cephalosporins and ceftiofur often test susceptible but still should be reported as resistant. For these organisms, the Amp C cephalosporinase may or may not be inducible. Third generation cephalosporins and/or aztreonam resistance implies acquisition of ESBL or hyperproduction of Amp C cephalosporinase.



Ceftriaxone **and/or** Ceftazidime  $\geq 2\mu\text{g}/\text{mL}$

## Chart 2G

*Serratia fonticola*

Multiple  
Beta - lactamases



Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO - R




**Comment**  
**&BL11**

*Serratia fonticola* produces both an inducible chromosomal cephalosporinase (Amp C) and an inducible, extended spectrum Class A enzyme conferring resistance to penicillins, amoxicillin/clavulanate and 1<sup>st</sup> generation cephalosporins.

If 3<sup>rd</sup> generation cephalosporins are resistant, all penicillins, 1<sup>st</sup>/2<sup>nd</sup>/3<sup>rd</sup> generation cephalosporins and beta-lactamase inhibitor combinations should be considered resistant. Acquisition of ESBL enzymes may occur, but would be difficult to detect.

Ceftriaxone and Ceftazidime < 2 µg/mL and Cefoxitin ≥ 32 µg/mL

**Chart 3**



## Chart 3A

Ceftriaxone and/or Ceftazidime < 2 µg/mL  
and  
Cefoxitin ≥ 32 µg/mL

### Interpret:

Ampicillin AM - R  
Cephalexin CEX - R (Urine Only) <sup>Note 1</sup>  
Cefazolin CZ - R  
Cefixime CFM - As Tested (Urine Only)  
Ceftriaxone CRO - S  
Amox/clav AMC oral - R  
Amox/clav IV AMC IV - R  
Pip/taz TZP - As Tested

*Citrobacter amalonaticus*  
*Citrobacter farmeri*  
*Citrobacter koseri*  
*Citrobacter sedlakii*  
*Edwardsiella spp.*  
*Escherichia spp.*  
*Klebsiella oxytoca*  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae*  
complex  
*Klebsiella variicola*  
*Proteus mirabilis*  
*Raoultella spp.*  
*Salmonella spp.*  
*Shigella spp.*

Note 1: Report for *E. coli*, *K. pneumo* and *P. mirabilis* only

*E. coli* and *Shigella* spp. produce a noninducible Amp C cephalosporinase in minute amount, not resulting in any significant beta-lactam resistance (although the MICs of cephalothin/cephalexin and occasionally amoxicillin/clavulanate may be slightly elevated). Cefoxitin resistance is due to either a permeability mutation or expression of Amp C cephalosporinase (hyperproduction of chromosomal Amp C or acquisition of plasmid Amp C cephalosporinase enzyme).

Hyperproduction of the chromosomal enzymes of *Klebsiella* spp., *Citrobacter koseri*, *Citrobacter amalonaticus* does not affect cefoxitin. Cefoxitin resistance indicates a permeability mutation or acquisition of AmpC cephalosporinase.

**Note:** Cefazolin and amoxicillin/clavulanate should be reported as R if either permeability mutation or cephalosporinase suspected.

Ceftriaxone and/or Ceftazidime < 2 µg/mL  
and  
Cefoxitin ≥ 32 µg/mL

## Chart 3B

*Proteus hauseri*  
*Proteus penneri*  
*Proteus vulgaris*

### Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO – Do Not Report  
Amox/clav AMC oral - R  
Amox/clav IV AMC IV - R  
Pip/taz TZP - As Tested

These organisms produce an inducible chromosomal cephalosporinase (Class A cefuroximase). Hyperproduction of this enzyme should not affect cefoxitin. Cefoxitin resistance indicates a permeability mutation or acquisition of AmpC cephalosporinase.

Note: Cefazolin and amoxicillin/clavulanate should be reported as R if either permeability mutation or cephalosporinase suspected.

Ceftriaxone and/or Ceftazidime < 2 µg/mL  
**and**  
Cefoxitin ≥ 32 µg/mL

## Chart 3C

*Erwinia spp.*  
*Kluyvera spp.*  
*Rahnella spp.*

Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO - R

Comment  
&1744

These organisms produce a low level chromosomal Class A extended spectrum beta-lactamase (ESBL) conferring resistance to penicillins, and 1<sup>st</sup>/ 2<sup>nd</sup> generation cephalosporins, but not cefoxitin. Even hyperproduction of this enzyme does not affect cefoxitin. Cefoxitin resistance indicates a permeability mutation (porin mutation/upregulated efflux). Typically, acquisition of a plasmid-mediated Amp C cephalosporinase enzyme results in elevated MIC (≥ 2 µg/mL) to one or more 3<sup>rd</sup> generation cephalosporins. Like other ESBLs, resistance to 3<sup>rd</sup> generation cephalosporins may not be apparent in vitro, but should be implied. ESBL testing is not useful as all should test positive due to the chromosomal ESBL enzyme.

Ceftriaxone and/or Ceftazidime < 2µg/mL  
and  
Cefoxitin ≥ 32 µg/mL

## Chart 3D

### Interpret:

Ampicillin AM - R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO - R for CSF/Brain, Blood/Sterile body site/  
Endovascular catheter sites  
- Do not report for other sites  
Pip/taz TZP - R Exception: *Morganella morganii* - As Tested\*


*Cedecea spp.*  
*Citrobacter freundii complex*  
*Cronobacter spp.*  
*Enterobacter spp.*  
*Hafnia spp.*  
*Klebsiella aerogenes*  
*Lelliottia amnigena*  
*Morganella morganii*  
*Pantoeae spp.*  
*Pluralibacter gergoviae*  
*Providencia spp.*  
*Serratia spp.*

All of these organisms produce a chromosomal Amp C cephalosporinase. In many, it is inducible, and in some it is constitutively expressed at low or high levels. The use of 3<sup>rd</sup> generation cephalosporins, even if susceptible (< 2 µg/mL), should be avoided to prevent selection of derepressed mutants. Selection of derepressed mutants is more common in certain clinical settings, such as pneumonia, bacteremia and abscesses. For these infections, 3<sup>rd</sup> generation cephalosporins should not be reported. For urine isolates, where therapeutic options are limited, cefotaxime or ceftriaxone may be reported if the MIC is < 2 µg/mL, as the risk of selecting derepressed mutants is low, given the high drug concentration of these agents in urine. Ceftazidime is a strong selector of derepressed mutants and should not be reported on any of these isolates (regardless of MIC). Permeability mutations in association with Amp C cephalosporinase may result in carbapenem resistance.

\*As the AmpC cephalosporinase is poorly inhibited by beta-lactamase inhibitors, amoxicillin/clavulanate and piperacillin/tazobactam should be reported as resistant. **Exception:** The cephalosporinase of *Morganella morganii* is inhibited by tazobactam (but not clavulanate). Therefore piperacillin/tazobactam can be reported as tested.


Ceftriaxone and Ceftazidime < 2 µg/mL  
**and**  
 Cefoxitin ≤ 16 µg/mL  
**and**  
 Cefazolin ≥ 8 µg/mL

**Chart 4**



# Chart 4A

*Escherichia coli*  
*Escherichia fergusonii*  
*Escherichia vulneris*  
*Proteus mirabilis*



**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**

For *E. coli* and *P. mirabilis*, a cefazolin MIC  $\geq 8\mu\text{g/mL}$  implies a high level production of penicillinase, rather than a cephalosporinase or a low level penicillinase. Low level penicillinase would affect only ampicillin and cephalothin/cephalexin, but not cefazolin, while a cephalosporinase should result in elevated cefoxitin MIC.

**Comment #clx1** For uncomplicated lower UTI only. (auto adds to cephalixin) **Add to cephalixin**


**Comment #A335** Cefazolin susceptibility result is available upon request. (auto adds to organism for *E. coli*, *K. pneumoniae*, *P. mirabilis*) **Add to cephalixin**



Ceftriaxone and Ceftazidime < 2 µg/mL  
**and**  
 Cefoxitin ≤ 16 µg/mL  
**and**  
 Cefazolin ≥ 8 µg/mL

**Chart 4B**

*Proteus hauseri*  
*Proteus penneri*  
*Proteus vulgaris*



**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**

These organisms all have an inducible chromosomal cephalosporinase (Class A cefuroximase) that typically results in resistance to ampicillin and 1<sup>st</sup>/2<sup>nd</sup> generation cephalosporins, but not cefoxitin. Amoxicillin/clavulanate and piperacillin/tazobactam typically remain susceptible.

Ceftriaxone and Ceftazidime < 2 µg/mL  
and  
Cefoxitin ≤ 16 µg/mL  
and  
Cefazolin ≥ 8 µg/mL

## Chart 4C

*Erwinia spp.*  
*Kluyvera spp.*  
*Rahnella spp.*

### Interpret

Ampicillin AM – R  
Cefazolin CZ - R  
Cefixime CFM - R (Urine Only)  
Ceftriaxone CRO – R

Comment  
&1744

These organisms have a chromosomal ESBL that typically results in cefazolin resistance. Beta lactamase inhibitor combination drugs may test susceptible, but should be used with caution as the level of beta-lactamase inhibitor may not be sufficient to inhibit the ESBL enzyme produced. It is safest not to report these agents at all.

Ceftriaxone and Ceftazidime < 2 µg/mL

and

Cefoxitin ≤ 16 µg/mL

and

Cefazolin ≥ 8 µg/mL

## Chart 4D

*Cedeceae* spp.  
*Citrobacter freundii* complex  
*Cronobacter* spp.  
*Enterobacter* spp.  
*Hafnia* spp.  
*Klebsiella aerogenes*  
*Lelliottia amnigena*  
*Morganella* spp.  
*Pantoea* spp.  
*Pluralibacter gergoviae*  
*Providencia* spp.  
*Serratia* spp.

### Interpret

Ampicillin AM - R

Cefazolin CZ - R

Cefixime CFM - R (Urine Only)

Ceftriaxone CRO - R for CSF/Brain, Blood/Sterile body site/Endovascular catheter sites  
- Do Not Report for other sites


Pip/taz TZP - R Exception: *Morganella morganii* - As Tested\*

All of these organisms produce a chromosomal Amp C cephalosporinase. In most, it is inducible, and in some it is constitutively expressed at low or high levels. The use of 3<sup>rd</sup> generation cephalosporins, even if susceptible (< 2 µg/mL), should be avoided to prevent selection of derepressed mutants. Selection of derepressed mutants is common in certain clinical settings, such as pneumonia, bacteremia and abscesses. For these infections, 3<sup>rd</sup> generation cephalosporins should not be reported. For urine isolates, where therapeutic options are limited, cefotaxime or ceftriaxone may be reported if the MIC is < 2 µg/mL, as the risk of selecting derepressed mutants is low, given the high drug concentration of these agents in urine. Ceftazidime is a strong selector of derepressed mutants and should never be reported on any of these isolates (regardless of MIC). Permeability mutations in association with Amp C cephalosporinase may result in carbapenem resistance.

\*As this cephalosporinase is poorly inhibited by beta-lactamase inhibitors, amoxicillin/clavulanate and piperacillin/tazobactam should be reported as resistant. **Exception:** The cephalosporinase of *Morganella morganii* is inhibited by tazobactam (but not clavulanate). Therefore, piperacillin/tazobactam can be reported as tested.

# Chart 4E

*Citrobacter amalonaticus*  
*Citrobacter farmeri*  
*Citrobacter koseri*  
*Citrobacter sedlakii*  
*Edwardsiella spp.*  
*Salmonella spp.*  
*Shigella spp.*




**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**

*C. koseri* and *C. amalonaticus* produce a chromosomal class A penicillinase which confers resistance to penicillins (amoxicillin, piperacillin) and, to a lesser degree, to cephalothin/cephalexin. Cefazolin MIC ≥ 8µg/mL implies high level penicillinase production.

# Chart 4F

*Escherichia hermanii*  
*Klebsiella oxytoca*  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae*  
 complex  
*Klebsiella variicola*  
*Raoultella* spp.



**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**


*K. pneumoniae* and *K. oxytoca* produce a chromosomal class A penicillinase which confers resistance to penicillins (amoxicillin, piperacillin) and, to a lesser degree, to cephalothin/cephalexin. Cefazolin MIC ≥ 8µg/mL implies high level penicillinase production.

**Comment #clx1** For uncomplicated lower UTI only **Add to cephalixin**

**Comment #A335** Cefazolin susceptibility result is available upon request **Add to cephalixin**


Ceftriaxone and Ceftazidime < 2 µg/mL  
**and**  
 Cefoxitin ≤ 16 µg/mL  
**and**  
 Cefazolin ≤ 4 µg/mL

**Chart 5**



# Chart 5A

*Escherichia coli*  
*Escherichia fergusonii*  
*Escherichia vulneris*  
*Proteus mirabilis*



\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.

*E. coli* produce a non-inducible basal Amp C cephalosporinase that typically does not result in significant beta-lactam resistance (although in vitro the MICs of cephalothin/cephalexin and occasionally amoxicillin/clavulanate may be slightly elevated).


Comment #clx1 For uncomplicated lower UTI only. Add to cephalixin

Comment #A335 Cefazolin susceptibility result is available upon request. Add to cephalixin

Comment 21104: Cefazolin result to follow. If susceptible, this would be the preferred (or most narrow spectrum) parenteral cephalosporin.

# Chart 5B

*Citrobacter amalonoticus*  
*Citrobacter farmeri*  
*Citrobacter koseri*  
*Citrobacter sedlakii*  
*Edwardsiella spp.*



**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**

The chromosomal penicillinase of these organisms has variable activity against oral/narrow spectrum cephalosporins (cephalexin). It is safest to report them all as I/R but not S.

**Comment #A335** Cefazolin susceptibility result is available upon request.



Ceftriaxone and Ceftazidime < 2 µg/mL  
**and**  
Cefoxitin ≤ 16 µg/mL  
**and**  
Cefazolin ≤ 4 µg/mL

## Chart 5C

*Salmonella spp*  
*Shigella spp*

### Interpret:

Ampicillin AMP - As Tested  
Cefixime CFM - As Tested  
Ceftriaxone CRO - S

# Chart 5D

*Escherichia hermanii*  
*Klebsiella oxytoca*  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae*  
*complex*  
*Klebsiella variicola*  
*Raoultella spp.*

Ceftriaxone and Ceftazidime < 2 µg/mL  
**and**  
 Cefoxitin ≤ 16µg/mL  
**and**  
 Cefazolin ≤ 4 µg/mL

Urine

Yes

No

*E. hermanii*  
*K. oxytoca*  
*K. pneumoniae complex*  
*K. variicola*  
*Raoultella spp.*

*K. pneumoniae*

**Cefazolin Disc**  
 If Blood culture add comment: **21104**

S

I/R

**Interpret:**

**Interpret:**

**Interpret:**  
 Ampicillin AM – R  
 Cephalixin CEX - Do not report  
 Cefazolin CZ - Do not report **(21105)**  
 Cefixime CFM - As Tested  
 Amox/clav AMC - As Tested\*  
 Pip/taz TZP - As Tested\*

**Interpret:**  
 Ampicillin AM – R  
 Cephalixin CEX - S **(27402)**  
 Cefazolin CZ - Do not report **(21105)**  
 Cefixime CFM - As Tested  
 Amox/clav AMC - As Tested\*  
 Pip/taz TZP - As Tested\*

Ampicillin AM – R  
 Cefazolin CZ – S  
 Ceftriaxone CRO - S  
 Amox/clav AMC oral – As Tested\*  
 Amox/clav IV AMC IV – same as AMC oral  
 Pip/taz TZP – As Tested\*


Ampicillin AM – R  
 Cefazolin CZ – I/R  
 Ceftriaxone CRO - S  
 Amox/clav AMC oral – As Tested\*  
 Amox/clav IV AMC IV – same as AMC oral  
 Pip/taz TZP – As Tested\*

**\*Confirm results if piperacillin/tazobactam I/R and amoxicillin/clavulanate oral S. If TZP confirms as I/R report AMC oral as R.**

**Comment #clx1** For uncomplicated lower UTI only. **Add to cephalixin**  
**Comment #A335** Cefazolin susceptibility result is available upon request. **Add to cephalixin**  
**Comment** Cefazolin result to follow. If susceptible, this would be the preferred (or most narrow spectrum) parenteral cephalosporin.

Ertapenem  $\geq 1$   $\mu\text{g/mL}$  (I/R) **or** Meropenem  $\geq 2$   $\mu\text{g/mL}$  (I/R) **or**  
Imipenem  $\geq 2$   $\mu\text{g/mL}$  (I/R)

## Chart 6



**Comment 21143:** Further testing is being done to determine the mechanism of carbapenem resistance (CRO).

**IPC:** Infection Prevention and Control

**PR:** Preliminary report

**\***: Phone call to IPC not required if only I/R to ertapenem.


**Note 1:** Chromosomal Class A carbapenemases have been found in *Enterobacter* spp. (NMC, IMI), *Serratia marcescens* (SME) and *Serratia fonticola* (SFC). These genes have also been found on transmissible elements. It is necessary to perform molecular testing to definitely determine the type of carbapenemase.

**Note 2:** Certain OXA enzymes (eg. OXA-48) hydrolyze carbapenems but not cephalosporins.

**Note 3:** Vitek N390 card has a card limitation for *Serratia marcescens* and imipenem. Test imipenem by disc diffusion before reporting imipenem result.

# Chart 6A

*Morganella spp.*  
*Proteus spp.*  
*Providencia spp.*



**Comment 21143:** Further testing is being done to determine the mechanism of carbapenem resistance (CRO).

**Comment 28325:** Imipenem has intrinsically low activity against this organism.

**IPC:** Infection Prevention and Control

**PR:** Preliminary report


**Note 1:** *Proteus*, *Providencia* and *Morganella* spp. may have low-level imipenem resistance (2-8  $\mu\text{g}/\text{mL}$ ) due to altered penicillin binding protein (PBP2).

**Note 2:** Consult Tech II if discrepancy between disc diffusion and automated system for ertapenem or meropenem.

**Interpretation of NG-Test CARBA 5 for Enterobacterales other than Proteus species**

**Chart 6B**

**NG-Test CARBA 5 results**



- Probable Carbapenemase
- Notify Microbiologist to phone IPC & physician
- Comments: 23415 + 21286 + 21393
- Freeze isolate as indefinite save (IS)
- Send isolate to Reference Laboratory for carbapenemase gene testing. See Note 3.
- Wait for carbapenem disc diffusion results before reporting further carbapenem results
- When carbapenem disc diffusion results available, report all carbapenems as tested (Notes 1 & 2). (For *Morganella* and *Providencia* spp. do not report imipenem, add comment 28325 )


- Notify Microbiologist
- Wait for carbapenem disc diffusion results before reporting further carbapenems
- When carbapenem disc diffusion results available, report all carbapenems as tested (Notes 1 & 2). (For *Morganella* and *Providencia* spp. do not report imipenem, add comment 28325 )
- Comments: 21143 + 21286 + 21393
- Freeze isolate as indefinite save (IS)
- Send isolate to Reference Laboratory for carbapenemase gene testing. See Note 3.

**Note 1:** Some carbapenemase producing isolates are categorized as susceptible to one or more carbapenems and therefore carbapenems should be reported as tested.  
**Note 2:** Consult Tech II if discrepancy between disc diffusion and automated system for ertapenem, imipenem and/or meropenem.  
**Note 3:** Check for previous positive gene testing. If gene testing was done within 1 year on the same organism/from same body site, report the organism and refer to the gene testing results. If different body site consult Microbiologist.

**Comment 23415:** Phenotypic testing suggests probable CARBAPENEMASE.  
**Comment 21143:** Further testing is being done to determine the mechanism of carbapenem resistance (CRO).  
**Comment 21286:** Sent for carbapenemase gene testing to National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.  
**Comment 21393:** Treatment of this organism is complex and may require high dose combination antibiotic therapy. Infectious Diseases consultation (by phone consultation if on site specialty services don't exist) is advisable.  
**Comment 28325:** Imipenem has intrinsically low activity against this organism.  
**IPC:** Infection Prevention and Control

**Interpretation of NG-Test CARBA 5 NEGATIVE for Enterobacterales other than Proteus species**

**Chart 6C**



**Comments:**

- If ESBL: 21744 + 23339
- If AmpC: 21746 + 23339
- If ESBL + AmpC: 21744 + 23339
- If MAST Negative or Undetermined

**Mechanism:**

Set up cefepime ESBL discs.  
If cefepime ESBL positive: 21744 + 23339  
If cefepime ESBL negative: Consult

**Note:** Interpret beta-lactams according to Chart 2

**Comments:**

- If ESBL: 21744 + 23340 + \*
- If AmpC: 21746 + 23340 + \*
- If ESBL + AmpC: 21744 + 23340 + \*
- If MAST Negative or Undetermined

Mechanism: Set up cefepime ESBL discs.

If cefepime ESBL positive: 21744 + 23340 + \*

If cefepime ESBL negative: Consult

**Note:** Interpret beta-lactams according to Chart 2

- Notify Microbiologist
- Report all carbapenems as tested (Notes 1 & 2)  
(For *Morganella* and *Providencia* spp. do not report imipenem, add comment 28325)
- Comments: 21143 + 21286 + 21393
- Freeze isolate as indefinite save (IS)
- Phone IPC if not previously phoned
- Send isolate to Reference Laboratory for carbapenemase gene testing  
**Note:** Check for previous positive gene testing. If gene testing was done within 1 year on the same organism/from same body site, report the organism and refer to the gene testing results. If different body site consult Microbiologist.


**Chart 6E**

**Comment 21143:** Further testing is being done to determine the mechanism of carbapenem resistance (CRO).  
**Comment 21286:** Sent for carbapenemase gene testing to National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.  
**Comment 21393:** Treatment of this organism is complex and may require high dose combination antibiotic therapy. Infectious Diseases consultation (by phone consultation if on site specialty services don't exist) is advisable.  
**Comment 23339:** This organism exhibits resistance to ertapenem (CRO) likely mediated by impermeability. Susceptible carbapenems should be used with caution.  
**Comment 23340:** This organism exhibits resistance to carbapenems (CRO) likely mediated by impermeability.  
**Comment 28325:** Imipenem has intrinsically low activity against this organism.  
 \* add: Susceptible carbapenems should be used with caution.  
 IPC: Infection Prevention and Control

*Proteus spp.*


**Interpretation of NG-Test CARBA 5 for Proteus species**

NG-Test CARBA 5 results  
Note 3



**Comment 23415:** Phenotypic testing suggests probable CARBAPENEMASE.  
**Comment 21143:** Further testing is being done to determine the mechanism of carbapenem resistance (CRO).  
**Comment 21286:** Sent for carbapenemase gene testing to National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.  
**Comment 21393:** Treatment of this organism is complex and may require high dose combination antibiotic therapy. Infectious Diseases consultation (by phone consultation if on site specialty services don't exist) is advisable.  
**Comment 28325:** Imipenem has intrinsically low activity against this organism.  
**IPC:** Infection Prevention and Control

**Note 1:** Some carbapenemase producing isolates are categorized as susceptible to one or more carbapenems and therefore carbapenems should be reported as tested.  
**Note 2:** Consult Tech II if discrepancy between disc diffusion and automated system for ertapenem and/or meropenem.  
**Note 3:** False negative NG-Test CARBA 5 results may occur with Proteus species.  
**Note 4:** Check for previous positive gene testing. If gene testing was done within 1 year on the same organism/from same body site, report the organism and refer to the gene testing results. If different body site consult Microbiologist.



**Comment 23340:** This organism exhibits resistance to carbapenems (CRO) likely mediated by impermeability.  
**IPC:** Infection Prevention and control

**Note 1:** Most likely a permeability mutation. Permeability mutations are most common in *Klebsiella* and *Enterobacter* spp. Although they typically affect only ertapenem, they can also affect other carbapenems especially if concurrent ESBL or AmpC present. In addition to ertapenem, permeability mutations may preferentially affect meropenem (*Klebsiella* spp.) or imipenem (*Enterobacter* spp.).



**Table 1**  
**Carbapenemase Gene Testing Result Codes**

CODE	COMMENT
<b>21252+21314</b>	Positive for KPC gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21251+21314</b>	Positive for IMP gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21266+21314</b>	Positive for VIM gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21255+21314</b>	Positive for NDM gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21249+21314</b>	Positive for GES gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21261+21314</b>	Positive for OXA-48 gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21265+21314</b>	Positive for SME gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21256+21314</b>	Positive for NMC gene. Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2.
<b>21201+21314</b>	Negative for carbapenemase genes (GES, IMP, KPC, NDM, NMC, OXA and VIM). Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2. <b>Note: add "-48" to OXA.</b>
<b>21209+21314</b> (For <i>Serratia</i> spp.)	Negative for carbapenemase genes (GES, IMP, KPC, NDM, NMC, OXA, SME and VIM). Testing performed by the National Microbiology Laboratory (NML), 1015 Arlington Street, Winnipeg, MB, Canada, R3E 3R2. <b>Note: add "-48" to OXA.</b>

### BETA-LACTAM RESISTANCE INTERPRETATION CHART

SIS Code	Comment	Ampicillin	Cephalexin	Cefazolin	Cefotaxime	Ceftriaxone	Ceftazidime	Cefixime	Cefepime	Aztreonam	Piperacillin/ Tazobactam	Amoxicillin/ Clavulanate	Ertapenem	Imipenem	Meropenem	Probable mechanism of resistance	Phone Infection Prevention and Control?
&1744	This organism is considered resistant to penicillins and cephalosporins. Amoxicillin-clavulanate and piperacillin-tazobactam are generally NOT recommended.	R	R	R	R	R	R	R	R	R	Do Not Report	R	S	S	S	ESBL	No
&1746	This organism is considered resistant to penicillins and cephalosporins (not including cefepime). Amoxicillin-clavulanate and piperacillin-tazobactam are generally NOT recommended.	R	R	R	R	R	R	R	Use With Caution	R	Do Not Report	R	S	S	S	AmpC cephalosporinase	No
&1744	This organism is considered resistant to penicillins and cephalosporins. Amoxicillin-clavulanate and piperacillin-tazobactam are generally NOT recommended.	R	R	R	R	R	R	R	R	R	Do Not Report	R	S	S	S	ESBL + AmpC cephalosporinase	No
&BL07	Beta-lactam resistance of this organism is mediated by an undetermined mechanism of resistance. All penicillins, cephalosporins and aztreonam should be considered resistant.	R	R	R	R	R	R	R	Do Not Report	R	Do Not Report	R	S	S	S	Indeterminate	No

**BETA-LACTAM RESISTANCE INTERPRETATION CHART**

SIT Code	Comment	Ampicillin	Cephalexin	Cefazolin	Cefotaxime	Ceftriaxone	Ceftazidime	Cefixime	Cefepime	Aztreonam	Piperacillin/ Tazobactam	Clavulanate	Amoxicillin/ Clavulanate	Ertapenem	Imipenem	Meropenem	Probable mechanism of resistance	Phone Infection Prevention and Control?
&BL19	Beta-lactam resistance of this organism is likely mediated by hyperproduction of its chromosomal beta-lactamase conferring resistance to all penicillins and cephalosporins, aztreonam and beta-lactamase inhibitor combination drugs.	R	R	R	R	R	Do Not Report	R	R	R	Do Not Report	R	S	S	S	Hyper-production of chromosomal enzymes	No	
&BL11	Beta-lactam resistance of this organism is likely mediated by multiple beta-lactamases conferring resistance to all penicillins and cephalosporins, aztreonam and beta-lactamase inhibitor combination drugs.	R	R	R	R	R	R	R	R	R	Do Not Report	R	S	S	S	Multiple beta-lactamases	No	

CODE	COMMENT
&BL11	Beta-lactam resistance of this organism is likely mediated by multiple beta-lactamases conferring resistance to all penicillins and cephalosporins, aztreonam and beta-lactamase inhibitor combination drugs.
&BL07	Beta-lactam resistance of this organism is mediated by an undetermined mechanism of resistance. All penicillins, cephalosporins and aztreonam should be considered resistant.
#cz1	Cefazolin result to follow. If susceptible, this would be the preferred (or most narrow spectrum) parenteral cephalosporin.
#A335	Cefazolin susceptibility result is available upon request.
&2143	Further testing is being done to determine the mechanism of carbapenem resistance (CRO).
&2130	Susceptibility testing for this organism was performed by a non-reference method and/or required modifications to the standard test conditions.
&1744	This organism is considered resistant to penicillins and cephalosporins. Amoxicillin-clavulanate and piperacillin-tazobactam are generally NOT recommended.
&1746	This organism is considered resistant to penicillins and cephalosporins (not including cefepime). Amoxicillin-clavulanate and piperacillin-tazobactam are generally NOT recommended.
&BL19	Beta-lactam resistance of this organism is likely mediated by hyperproduction of its chromosomal beta-lactamase conferring resistance to all penicillins and cephalosporins, aztreonam and beta-lactamase inhibitor combination drugs.

<b>DRUG NAME</b>	<b>CODE</b>
Amoxicillin/Clavulanate	AMC
Ampicillin	AM
Aztreonam	ATM
Cefazolin	CZ
Cefepime	FEP
Cefixime	CFM
Cefotaxime	CTX
Cefoxitin	FOX
Ceftazidime	CAZ
Ceftriaxone	CRO
Cephalexin	CEX
Ertapenem	ETP
Imipenem	IPM
Meropenem	MEM
Piperacillin-tazobactam	TZP