

Haemophilus influenzae (LTR79351)

Edit Approved By: Dragan, Tatiana (07/31/2023)

Revision: 4.00

Organism **Haemophilus influenzae**

Clinical *H. influenzae* may colonize the nasopharynx, conjunctiva, and occasionally the genital tract of humans.

This organism may also colonize the lower respiratory tract of individuals with chronic obstructive pulmonary disease (COPD). *H. influenzae* may cause meningitis, epiglottitis, cellulitis (including periorbital cellulitis), septic arthritis, osteomyelitis, otitis media, sinusitis, acute exacerbations of chronic bronchitis, pneumonia, empyema, conjunctivitis, bacteremia, neonatal and maternal sepsis, and rarely urinary tract infections.

**Usual
susceptibility
pattern**

These organisms have variable susceptibility to ampicillin, and TMP-SMX. Ampicillin resistance is most often mediated by β -lactamase production and less frequently by mutation in the penicillin binding protein (PBP3). Alteration of PBP3 usually affects only ampicillin, cefaclor and minimally cefuroxime. One specific mutation may also affect third generation cephalosporins. Combination of both β -lactamase and altered PBP results in very high level resistance to ampicillin. Although resistance has been described, most strains remain susceptible to third generation cephalosporins, rifampin, quinolones, aminoglycosides, tetracycline and chloramphenicol. **These organisms are resistant to first generation cephalosporins.** Of the oral cephalosporins, cefuroxime and cefixime have the best activity, while the activity of cefprozil and cefaclor is less reliable. Although azithromycin has the best activity, the activity of macrolides is not reliable. Resistance to macrolides is mediated mostly by efflux mechanisms. Rare imipenem or meropenem resistant isolates have been observed mainly among β -lactamase negative ampicillin resistant isolates.

**Susceptibility
method**

Modified Kirby-Bauer or Etest method using Haemophilus Test Medium (HTM) incubated in 5% CO₂ at 35°C for 16-18 hours (20-24 hours for Etest).

Note: For Etest use 0.5 McFarland suspension in broth. For mucoid strains use 1.0 McFarland

* 1. Eye samples: Follow ASTM

* 2. Resp/Ear samples:

- do BL
- if BL negative, add drug comment #BLN
- if BL positive, add drug comment #BLP

Haemophilus influenzae, Continued

Susceptibility reporting

	CSF/ Brain	Blood	Sterile Body Site	Eye * (See Note)	Other*	Comments
β-lactamase	*	*	*	*	*	Test but do not report
Amox-Clav					2	2 nd line if β-lactamase positive or anaerobes/Moraxella/ Staph aureus present on smear or culture If amox-clav I/R see Special Considerations
Ampicillin	✓*	✓*	✓*	✓	✓	*Etest method If β-lactamase positive – report amp R If ampicillin I/R and β-lactamase negative see Special Considerations
Cefotaxime	*	*	*	*	3*	*Report if patient <1 month instead of ceftriaxone Etest method 3 rd line if cefuroxime I/R See Special Considerations
Ceftriaxone	✓	✓	✓	✓	3	Do not report if patient <1 month Etest method 3 rd line if cefuroxime I/R See Special Considerations
Cefuroxime		2	2	2	2	2 nd line if beta lactamase positive If cefuroxime I/R see Special Considerations
Ciprofloxacin		✓	✓	✓		Do not report in patients < 18 y Etest method See Special Considerations
Meropenem	2	2	2		4	Etest method 2 nd or 4 th line if cefotaxime/ceftriaxone not susceptible See Special Considerations
Rifampin	*	*	*			*Report if rifampin I/R If rifampin S see Special Considerations
Tetracycline					✓	Do not report in patients < 8 y
TMP-SMX			✓	✓	✓	

Note

Deep eye specimens:	Perform susceptibility test if:		
	• vitreous fluid	• canaliculitis	• corneal ulcer / scrapings
	• chamber aspirate	• endophthalmitis	• contact lens related infections
	• intraocular fluid	• donor sclera	• ophthalmology clinic/ward
	• keratitis	• chorioretinitis	• history of failure of therapy
	• injury/surgery	• cornea	• preseptal/orbital cellulitis
Superficial eye specimens:	Susceptibility testing of superficial eye specimens not routinely performed. Add comment: “Susceptibility testing of topical antibiotics is not standardized and is not routinely performed on superficial eye specimens”.&A89		

Haemophilus influenzae, Continued

Special considerations

<u>Amoxicillin-clavulanate / Ampicillin / Cefuroxime:</u>	<p>Ampicillin resistance is usually due β-lactamase production. In some strains that are ampicillin resistance and β-lactamase negative (BLNAR) resistance is due to alteration of penicillin binding proteins.</p>	
	IF...	THEN...
	<p>Ampicillin I/R and β-lactamase negative</p>	<ul style="list-style-type: none"> Do not report amoxicillin/clavulanate or cefuroxime Repeat β-lactamase testing Confirm ampicillin by Etest If confirmed as ampicillin I/R and β-lactamase negative, this indicates an altered penicillin binding protein mechanism of resistance (BLNAR). Consult microbiologist Ampicillin, cefuroxime and amoxicillin/clavulanate should be reported as R.
	<p>Amoxicillin-clavulanate I/R and β-lactamase positive (ampicillin R)</p>	<ul style="list-style-type: none"> Do not report amoxicillin/clavulanate or cefuroxime Repeat β-lactamase testing Confirm amoxicillin/clavulanate by Etest Consult microbiologist If confirmed as β-lactamase positive and amoxicillin/clavulanate I/R, this indicates multiple mechanisms of resistance including both β-lactamase production and altered penicillin-binding proteins (BLPACR). Ampicillin, cefuroxime and amoxicillin/clavulanate should be reported as R.
	<p>Amoxicillin-clavulanate I/R and Ampicillin S</p>	<ul style="list-style-type: none"> Unusual susceptibility pattern Do not report ampicillin or amoxicillin/clavulanate Confirm ampicillin and amoxicillin/clavulanate by Etest Consult microbiologist
	<p>Cefuroxime I/R</p>	<ul style="list-style-type: none"> Confirm cefuroxime by Etest Consult microbiologist
<u>Cefotaxime/ Ceftriaxone/ Meropenem:</u>	<p>Haemophilus sp. should be susceptible to these antibiotics. Consult Supervisor if not susceptible. If I/R, the identification of the organism and its susceptibility should be confirmed by repeat testing. Consider sending confirmed isolates to reference laboratory.</p>	

Haemophilus influenzae, Continued

Special considerations (continued)

<u>Ciprofloxacin:</u>	Current CLSI breakpoints may underestimate resistance. Use EUCAST breakpoints for interpretation: <table border="1" data-bbox="410 327 833 438"><thead><tr><th>MIC</th><th>Interpretation</th></tr></thead><tbody><tr><td>≤ 0.06 µg/mL</td><td>S</td></tr><tr><td>≥ 0.12 µg/mL</td><td>R</td></tr></tbody></table> Add comment: "Interpretation is based upon EUCAST breakpoints." #2117	MIC	Interpretation	≤ 0.06 µg/mL	S	≥ 0.12 µg/mL	R
MIC	Interpretation						
≤ 0.06 µg/mL	S						
≥ 0.12 µg/mL	R						
<u>Rifampin:</u>	If rifampin susceptible report only upon microbiologist's request. Add comment: "Rifampin is appropriate ONLY for prophylaxis of case contacts. These interpretations do NOT apply to therapy of patients with invasive disease." (23751)						

Interpretation For Etest, report actual MIC result. For interpretation (S, I, or R) report according to the nearest higher doubling dilution (**Appendix 1**).

Use **CLSI** interpretive document for *Haemophilus influenzae* and *Haemophilus parainfluenzae*.

For ciprofloxacin - Refer to **Special Considerations**