



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Name:
Microbiology Planting Job Aid

Distribution:
Microbiology Specimen Processing Manual

Document Number:
MIC10140

Date Issued:
December 10, 2019

Uncontrolled When Printed

SPECIMEN TYPE	GRAM?	MEDIA	INCUBATION
Bacterial vaginosis	Yes	Only slide, no media	Only slide, no media
Blood culture: <ul style="list-style-type: none">Positive bottles>24 hour bottles	Yes	<u>Positive bottle:</u> BA, CHO, MAC, BRU <u>>24 hour bottle:</u> BA, CHO, BRU	BA and CHO = CO ₂ MAC = O ₂ BRU = AnO ₂
Body fluid culture	Yes	BA, CHO, MAC, BRU, THIO	BA and CHO = CO ₂ MAC and THIO = O ₂ BRU = ANO ₂
CSF culture	Yes	BA, CHO, MAC THIO if from shunt	BA and CHO = CO ₂ MAC (and THIO if indicated) = O ₂
Ear culture	Yes	BA, CHO, CNA, MAC	BA, CHO and CNA = CO ₂ MAC = O ₂
Eye culture-Superficial	Yes	BA, CHO	BA and CHO = CO ₂
Eye culture-Deep	Yes	BA, CHO, MAC, BRU	BA and CHO = CO ₂ MAC = O ₂ BRU = AnO ₂
Genital culture-Upper	Yes	BA, CHO, CNA, TM, MAC, BRU, KV	BA, CHO, CNA and TM = CO ₂ MAC = O ₂ BRU and KV = AnO ₂
Genital culture-Lower	Yes	BA, CHO, CNA, TM, MAC	BA, CHO, CNA and TM = CO ₂ MAC = O ₂
Gonorrhoeae culture	Male urethra	CHO, TM	CHO and TM = CO ₂
GBS screen	No	LIM, GBS	LIM = CO ₂
IUD culture	No	BRU, THIO	BRU = AnO ₂ THIO = O ₂
MRO	No	MRS, VRE	MRS = O ₂ VRE = O ₂
MRSA screen	No	MRS	MRS = O ₂

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Name:
Microbiology Planting Job Aid

Distribution:
Microbiology Specimen Processing Manual

Document Number:
MIC10140

Date Issued:
December 10, 2019

Uncontrolled When Printed

SPECIMEN TYPE	GRAM?	MEDIA	INCUBATION
Nose culture	No	BA	BA = O ₂
Oral culture	No	YST	YST = O ₂
Respiratory culture	Yes	BA, CHO, CNA, MAC	BA, CHO, CNA = CO ₂ MAC = O ₂
Throat culture	No	BA	BA = AnO ₂
Tip culture	No	BA, MAC	BA = CO ₂ MAC = O ₂
Urine culture	No	URI	URI = O ₂
VRE screen	No	VRE	VRE = O ₂
Wound culture-Superficial	Yes	BA, CNA, MAC	BA and CNA = CO ₂ MAC = O ₂
Wound culture-Deep	Yes	BA, CHO, CNA, MAC, BRU, KV	BA, CHO and CNA = CO ₂ MAC = O ₂ BRU and KV = AnO ₂
Yeast culture	No	YST	YST = O ₂

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC10140-Microbiology Planting Job Aid

PRINT DATE: 1/11/2022 11:03 AM