

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC32000 – Oral Culture	Policy Number: 15-137-V2
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s): NA	
Effective Date: 06/05/2025	Next Review Date: 06/05/2027
Issuing Authority: Director, Laboratory and Diagnostic Imaging Services	Date Approved: 06/05/2025
Accreditation Canada Applicable Standard: NA	

**Uncontrolled When Printed**

**GUIDING PRINCIPLE:**

Mouth swabs are usually obtained to identify oral yeast infections (thrush). Thrush is a common infection of the mouth and tongue caused by *Candida* spp. such as *Candida albicans*, *Candida glabrata*, *Candida tropicalis*, and *Candida krusei*.

**PURPOSE/RATIONALE:**

This standard operating procedure describes the screening for *Candida* spp. in oral (mouth) specimens.

**SCOPE/APPLICABILITY:**

This standard operating procedure applies to Medical Laboratory Technologists (MLTs) processing specimens for oral culture.

**SAMPLE INFORMATION:**

<b>Type</b>	Swab • Amie's with or without charcoal
<b>Source</b>	Mouth or tongue
<b>Stability</b>	If the sample is received in the laboratory and processed greater than 48 hours from collection: • Add specimen quality comment: "Delayed transport may adversely affect pathogen recovery"
<b>Storage Requirements</b>	Room temperature
<b>Criteria for rejection</b>	1. Unlabeled/mislabeled swabs 2. Specimen container label does not match patient identification on requisition

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

### REAGENTS and/or MEDIA:

- CandiSelect agar (YST)

### SUPPLIES:

- Disposable inoculation needles
- Microscope slides
- Sterile saline
- Wooden sticks
- Coverslips

### EQUIPMENT:

- Biosafety cabinet
- Microscope
- VITEK 2 and supplies

### SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures:

- Ensure that appropriate hand hygiene practices be used
- Lab gown must be worn when performing activities with potential pathogens
- Gloves must be worn when direct skin contact with infected materials is unavoidable
- Eye protection must be used when there is a known or potential risk of exposure of splashes
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC)
- The use of needles, syringes and other sharp objects should be strictly limited

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

### QUALITY CONTROL:

- Refer to MIC60040-Culture Media Quality Control procedure
- Refer to Test Manual for reagent quality control procedures

### PROCEDURE INSTRUCTIONS:

Step	Action
<b>Processing specimens for oral culture</b>	
<b>1</b>	In the biosafety cabinet: <ul style="list-style-type: none"><li>• Inoculate YST agar with the swab</li><li>• Ensure all surfaces of the swab make contact with the agar</li><li>• Streak for isolated growth using a disposable inoculation needle</li></ul>
<b>2</b>	Label the YST plate with: R (Date + 2 day).

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

<b>3</b>	Incubate the media: <ul style="list-style-type: none"> <li>Place YST plate in the O<sub>2</sub> incubator on the "YST SCREEN" section on the old cultures shelf</li> </ul>
----------	--

#### INTERPRETATION OF RESULTS:

Step	Action			
<b>1</b>	<ul style="list-style-type: none"> <li>Observe YST plate at 48 hours</li> <li>Examine for the presence of yeast growth, characterized by colored colonies on the agar</li> </ul>			
<b>2</b>	If no colonies resembling yeast are seen at 48 hours: <ul style="list-style-type: none"> <li>Record observations in the LIS</li> <li>Workup complete</li> <li>No yeast isolated</li> </ul>			
<b>3</b>	If colonies resembling yeast are seen: <ul style="list-style-type: none"> <li>Record observations in the LIS</li> <li>Perform a wet prep of the colonies</li> </ul>			
<b>4</b>	<b>IF</b>			
	<b>THEN</b>			
	<table border="1"> <tr> <td style="text-align: center;">Wet prep NEGATIVE</td><td> <ul style="list-style-type: none"> <li>Workup complete</li> <li>Yeast not isolated</li> </ul> </td></tr> <tr> <td style="text-align: center;">Wet prep POSITIVE</td><td> <ul style="list-style-type: none"> <li>Perform VITEK 2 YST card</li> </ul> </td></tr> </table>	Wet prep NEGATIVE	<ul style="list-style-type: none"> <li>Workup complete</li> <li>Yeast not isolated</li> </ul>	Wet prep POSITIVE
Wet prep NEGATIVE	<ul style="list-style-type: none"> <li>Workup complete</li> <li>Yeast not isolated</li> </ul>			
Wet prep POSITIVE	<ul style="list-style-type: none"> <li>Perform VITEK 2 YST card</li> </ul>			

#### REPORTING INSTRUCTIONS:

IF	REPORT
Yeast not isolated	<ul style="list-style-type: none"> <li>Report: <b>"No Yeast Isolated"</b></li> </ul>
Yeast Isolated, <i>Candida</i> spp.	<ul style="list-style-type: none"> <li>Report organism identification</li> <li>List quantitation as <b>"Isolated"</b></li> </ul>
Yeast Isolated, not <i>Candida</i> spp.	<ul style="list-style-type: none"> <li>Add isolate: <b>"Yeast (NOT Candida species)"</b></li> <li>Use LIS OrgID: <b>ystnot</b></li> <li>List quantitation as <b>"Isolated"</b></li> </ul>

#### LIMITATIONS:

- The VITEK 2 YST card provides an acceptable identification for *Candida* spp. and results can be reported as tested.
- If yeast other than *Candida* spp. are isolated, the identification should not be reported.

#### CROSS-REFERENCES:

- MIC60040-Culture Media Quality Control

#### REFERENCES:


- Leber, A. (2016). *Clinical microbiology procedures handbook*. (4<sup>th</sup>ed.) Washington, D.C.: ASM Press
- Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. (2015). *Manual of Clinical Microbiology*, 11<sup>th</sup> edition. Washington, D.C: ASM Press

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

## APPROVAL:

May 06, 2025

Date



Acting Director, Laboratory and Diagnostic Imaging Services

## REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	06 Nov 17	Initial Release	L. Steven
2.0	04 Dec 19	Updated to include new VITEK 2 instrument	L. Steven
3.0	31 Dec 21	Procedure reviewed and added to NTHSSA policy template	L. Steven
4.0	31 Aug 23	Procedure reviewed	L. Steven
5.0	31 Mar 25	Procedure reviewed and updated to reflect new yeast agar CandiSelect	L. Steven

**Disclaimer Message:** This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.