

PROGRAM Standard Operating Procedure – Laboratory Services	
Title: MIC71000 – BACTEC FX Instrument	Policy Number:
Program Name: Laboratory Services	
Applicable Domain: Lab, DI and Pharmacy Services	
Additional Domain(s): NA	
Effective Date:	Next Review Date:
Issuing Authority: Director, Laboratory and Diagnostic Imaging Services	Date Approved:
Accreditation Canada Applicable Standard: NA	

Uncontrolled When Printed

GUIDING PRINCIPLE:

The BACTEC FX Instrument is designed for the rapid detection of bacteria in clinical specimens. Blood cultures are collected from patients with suspected sepsis or bacteremia. Although primarily directed towards the processing of blood cultures, occasionally other specimen types (sterile fluids, abscess material, bone marrow, etc.) are received in blood culture bottles. These bottles may be processed in the same manner as blood cultures.

The BACTEC FX Instrument continuously monitors routine blood cultures for evidence of growth for 5 days. Negative results are auto-verified as follows:

- No growth after 48 hours (preliminary)
- No growth after 5 days (final)

PURPOSE/RATIONALE:

This standard operating procedure describes the BACTEC FX Instrument and its components.

SCOPE/APPLICABILITY:

This standard operating procedure applies to Medical Laboratory Technologists (MLTs) and Medical Laboratory Assistants (MLAs) processing specimens using the BACTEC FX Instrument.

SAMPLE INFORMATION:

Type	Blood culture bottle
Source	Blood or sterile fluid

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

REAGENTS and/or MEDIA:

Type	<ul style="list-style-type: none">• BACTEC Plus Aerobic Culture Bottles• BACTEC Lytic Anaerobic Culture Bottles• BACTEC Peds Plus Culture Bottles
Stability	Stable until date of expiration indicated on bottle
Storage Requirements	Bottle storage before blood collection: <ul style="list-style-type: none">• Room temperature
Criteria for rejection	Do not use if: <ul style="list-style-type: none">• The expiration date has passed• There are other signs of deterioration

EQUIPMENT:

- BACTEC FX Instrument
- EpiCenter computer

ENVIRONMENTAL CONTROLS:

- Operating temperature: 18°C to 30°C
- Relative humidity: 25% to 80%, non-condensing

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures:

- Ensure that appropriate hand hygiene practices be used
- Lab gown must be worn when performing activities with potential pathogens
- Gloves must be worn when direct skin contact with infected materials is unavoidable
- Eye protection must be used when there is a known or potential risk of exposure of splashes
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC)
- The use of needles, syringes and other sharp objects should be strictly limited

All patient specimens are assumed to be potentially infectious. Routine Practices must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

PROCEDURE INSTRUCTIONS:


Step	Action
Instrument components	
1	BACTEC FX Instrument: <ul style="list-style-type: none">• The BACTEC FX Instrument is composed of 2 drawers: A and B• Drawers are divided into columns (numbered 1 to 10 from left to right) and rows of stations (lettered A to K excluding I)• There are a total of 100 stations in each drawer

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

2

Status display:

- The Status display is the main display shown when no other operation has been initiated or is in progress:



- Once a drawer is opened, you can initiate the major instrument activities from the Status display. Vial entry, remove positive vials, remove negative vials and identify anonymous vials can be initiated for any drawer in the instrument

3

System indicators:

- The system indicators are located on the front-center of the instrument:



- The system indicators inform you of various states in the instrument:

Indicator Colour	State	Meaning
Yellow=Light in unison for instrument	On	<ul style="list-style-type: none"> System Alert Indicator remains on until the condition is corrected/addressed
Green=One for each for left and right drawers	On	<ul style="list-style-type: none"> Negative bottle Indicator remains lit until all negative bottles are removed through the Remove Negative Vials activity
Red=One for each left and right drawers	On	<ul style="list-style-type: none"> Positive bottle Indicator remains lit until all positive bottles are removed through the Remove Positive Vials activity

4

Station indicators:

- Each station has a set of LED indicators that inform you of the stations or bottles status
- The status indicator is located above each station:



- The color (red, green or yellow) and state (on, off or flashing) indicate the bottle status:

Colour	State	Meaning
Red	Flashing	Positive bottle
Green	Flashing	Negative bottle
Yellow	Flashing	Anonymous bottle
Red Yellow	Flashing	Positive anonymous bottle
Green	On	Available station


5







Bottle statuses:

- From the Status display, select the Drawer View button:



- The status of bottles is displayed on the Drawer View Display:




	<ul style="list-style-type: none"> The following Station statuses are shown: 		
	Status	Icon	Meaning
	Available		There is no bottle in the station
	Blocked		User has manually blocked the station
	Negative		Bottle completed protocol with no evidence of positivity
	Ongoing		Bottle is in the instrument and is in protocol
	Positive		Instrument has detected evidence of microbial growth
6	Anonymous		Bottle was physically placed in instrument without barcode sequence number being scanned
	Audible tones and alarms: <ul style="list-style-type: none"> Numerous different sounds are generated by the BACTEC FX Instrument: 		
	Type	Example	Sound
	Activity complete	All negative bottles were removed	High pitched tone repeated 3 times
	Activity error	Did not scan accession barcode after scanning sequence barcode	Single high beep
	Anonymous	Anonymous bottle entered	Short buzz sound
	Drawer ajar	The door is not closed	2 tones, high then low frequency, repeating
	Positive bottle	A positive bottle is detected	Pulsing fading sound, repeated
	Bottle entry	A bottle was entered	High pitch chirp sound


Step	Action
Loading bottles into the BACTEC FX	
1	Open the drawer. Ensure that the Status screen is displayed, and the barcode scanner turns on.
2	Scan the bottle sequence barcode (bottle barcode) and the accession barcode (LIS label).
3	Place the bottle into any available slot (solid green light) in the instrument.
4	Scan and place any other bottles that need to be loaded.
5	Close drawer when finished.


Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.



6	Delayed entry of bottles may lead to delayed results. Best practices dictate that bottles should be placed into the BACTEC FX as soon as possible after collection: <ul style="list-style-type: none"> All bottles received >24 hours after collection need to be sub-cultured to ensure growth is not present Refer to MIC10100-Microbiology Specimen Processing
7	Bottles can be placed into available (GREEN indicator) stations without being scanned into the instrument. Bottle that are not scanned are called "anonymous" bottles.

Step	Action
Identifying anonymous bottles	
1	The Status display will indicate which drawer houses the anonymous bottle: 
2	Open the drawer and locate station(s) with a flashing Yellow LED and remove the bottle. NOTE: If the station flashes YELLOW and RED, then the bottle is also POSITIVE
3	The ID Anonymous display appears, and the barcode scanner turns on.
4	Ensure the bottle has an accession barcode (LIS label). If not, accession the sample and place the barcode on the bottle.
5	Scan the bottle sequence barcode (bottle barcode) and the accession barcode (LIS label).
6	Place the bottle in the FLASHING GREEN station (station from which the bottle was pulled).
7	Repeat with any additional anonymous bottles.



Step	Action
Removing positive bottles	
1	The Status display will indicate which drawer houses the positive bottle:


	
2	Open the drawer and locate the station(s) with a flashing Red LED and remove the bottle(s).
3	The Positive Removal display appears. Scan the bottle sequence barcode (bottle barcode). NOTE: You must scan each positive bottle that you remove in order for the instrument to re-light positive stations
4	When all positive bottles are removed from the drawer, the Activity Complete tone sounds.

Step	Action
Removing negative bottles	
1	The Status display will indicate which drawer houses the negative bottle: 
2	Open the drawer and locate the station(s) with a flashing Green LED and remove the bottle (s). These bottles do not have to be scanned.
3	When all negative bottles are removed from the drawer, the Activity Complete tone sounds.



Step	Action
Extending incubation time of bottles	
1	On the BACTEC Status display, select the Drawer View button: 
2	Select the desired station and select OK.
3	Select Modify and the Modify Protocol box will be displayed: 



Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

4	Change the protocol length to 10 days and select OK. Select Save to save the changes.
5	<p>In the LIS, under Results Entry, enter the order number to access the plate log/test comments</p> <ul style="list-style-type: none"> In the plate log (Media Comments), look for the EXT media ID:  <ul style="list-style-type: none"> Double click in the resulted box (R) for the EXT Media ID so that a red check mark appears:  <p>NOTE: If a set of bottles has been collected, both bottles will need to be modified</p> <ul style="list-style-type: none"> This alerts the LIS to stop the 5 day reporting Save the order to save the changes made
6	<p>If culture is negative:</p> <ul style="list-style-type: none"> A 48 hour negative preliminary report will be automatically released A 10 day no growth final report will automatically be released <p>NOTE: Manual reporting negative bottles on extended protocol is not required</p>


Step	Action
Resolving system alert	
1	A system alert is indicated by a yellow LED indicator on both drawers of the instrument. This alert usually indicates a power failure or communication interruption.
2	<ul style="list-style-type: none"> Check the instrument for an error message On the Status display, select the Systems Alert button to view the Alert List:  <ul style="list-style-type: none"> Power interruptions will display the following alerts: <ul style="list-style-type: none"> ➤ Reboot Reason: Power fail ➤ The instrument has lost connectivity to the server ➤ EpiCenter Communications failure Refer to the BACTEC FX Instrument User Manual for alert descriptions


Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.



3	<ul style="list-style-type: none"> Log into EpiCenter computer: <ul style="list-style-type: none"> Log into Windows. When power goes out, Windows will re-boot and require re-login Icon will display:  Communication should RESYNC  after logging into Windows (should take about one minute after logging in) A System Message window should pop-up displaying the errors. Click the x button to close or delete them The yellow system indicator lights on BACTEC FX door will stop glowing
----------	--

Step	Action
Disassociating bottles	
1	If a bottle record contains an accession number, it is considered associated to that accession. The disassociate function enables you to break the link between a bottle and an accession number. This can be useful when troubleshooting sample errors.
2	Remove the bottle (s) from the BACTEC FX Instrument.
3	From the Status display, select the Culture button: 
4	The Culture-Patient display appears: 
5	Select the Vial tab to access the Culture-Vial display:

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

	
6	Scan the bottle sequence barcode (bottle barcode).
7	Select the Disassociate button to disassociate the bottle from the accession number.
8	Return to the Status display.
9	Open the drawer and proceed to load the bottle into the instrument.
10	Scan the bottle I sequence barcode (bottle barcode) and then scan the accession barcode (LIS label).
11	Place the bottle in any available slot in the instrument.

Step	Action
Associating bottles	
1	If a bottle has been loaded into the instrument without an accession number, the bottle needs to be associated with the accession number.
2	Remove the bottle (s) from the BACTEC FX Instrument.
3	From the Status display, select the Culture button: 
4	The Culture-Patient display appears:

	
5	<p>Select the Specimen tab to access the Culture-Specimen display:</p> 
6	In the accession field, scan the accession barcode (LIS label).
7	Scan the bottle sequence barcode (bottle barcode) you want to attach.
8	Select the Save button to save the association.

CROSS-REFERENCES:

- MIC10100-Microbiology Specimen Processing
- MIC20500-Gram stain resulting in LIS-Blood Cultures
- MIC60010-Microbiology Quality Control

REFERENCES:

1. Becton Dickinson and Company. (2016-12). *BD BACTEC FX Instrument User Manual*, 8005110(07)

APPROVAL:

 Date

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.

Director, Laboratory and Diagnostic Imaging Services

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	06 Nov 17	Initial Release	L. Steven
2.0	26 Mar 19	Updated to reflect addition of disassociating and associating bottles	L. Steven
3.0	16 Aug 21	Procedure reviewed and added to NTHSSA policy template	L. Steven
4.0	01 Oct 24	Procedure reviewed	L. Steven

Disclaimer Message: This is a **CONTROLLED** document for internal use only. Any documents appearing in paper form are not controlled and should be checked against the electronic file version prior to use.