

PROCEDURE: ANTIBIOTIC BATTERIES

The laboratory reviews CLSI M100 annually to update procedure for relevant organism/antimicrobial changes to ensure appropriateness of reporting. A Director actively participates with the Antimicrobial Stewardship Committee to address the needs of the patient population.

Non-standard susceptibility testing should not be performed. Always refer to CLSI M100 for standards.

All providers requesting recommendations on antimicrobial therapy should be directed to the Antimicrobial Stewardship Pharmacist at: 350-2205.

- *Haemophilus* species – perform cefinase only and report as the appropriate isolate comment: **&HBLN** (*Beta Lactamase Negative. Haemophilus isolates negative for beta-lactamase are likely to be susceptible to Amoxicillin, Macrolides and Cephalosporin antibiotics.*) or **&HBLP** (*Beta-Lactamase positive. Haemophilus Isolates producing beta lactamase are resistant to Amoxicillin.*)
- *Moraxella catarrhalis* - include isolate comment: **&BCAT** (*All isolates should be considered Beta-Lactamase positive.*)
- *Bacteroides fragilis* group – include isolate comment: **&BCAT** (*All isolates should be considered Beta-Lactamase positive.*)
- Antibiotics for MDR *Acinetobacter* screens are not routinely reported. The provider is looking for the presence or absence of the resistant organism. Refer to *Acinetobacter Screen Procedure* for guidance
- Additional antimicrobial requests – check CLSI standards for appropriateness (Table 1 & Table 2) & bring up on ROUNDS
ARUP is the reference laboratory for send-out testing.
Document all requests in worksheet.
(Reference laboratories require organisms to be submitted growing on an agar slant)
- Send isolate to RIDOH if the following reported result is: Ertapenem resistant *Enterobacterales*; Meropenem/Imipenem resistant *Enterobacterales*, non-mucoid *P. aeruginosa*, or *Acinetobacter baumannii* complex; Vancomycin confirmed as non-susceptible (MIC is 4) *Staphylococcus aureus*; pan-resistant organisms

KIRBY BAUER: GRAM NEGATIVE RODS (Non-CSF SOURCES)

<i>Routine Method</i>	<i>Back-up Method</i>	
<i>PSEUDOMONAS AERUGINOSA</i>	<i>ENTEROBACTERIACEAE</i>	<i>SALMONELLA / SHIGELLA</i>
Amikacin (AN)	Amikacin (AN)	Ampicillin (AMP)
Aztreonam (ATM)	Ampicillin (AMP)	Ceftriaxone (CRO) ^d
Cefepime (FEP)	Ampicillin/Sulbactam (SAM)	Ciprofloxacin (CIP)
Ceftazidime (CAZ)	Cefazolin (CZ)	Trimeth/Sulfa (SXT)
Ciprofloxacin (CIP)	Cefuroxime (CXM)	
Gentamicin (GM)	Ceftriaxone (CRO)	
Levofloxacin (LV)	Cefepime (FEP)	
Meropenem (MEM)	Ciprofloxacin (CIP)	
Piperacillin/Tazobactam (TZP)	Gentamicin (GM)	
Tobramycin (NN)	Meropenem (MEM)	
Ampicillin/Sulbactam (SAM) ^e	Trimeth/Sulfa (SXT)	
Trimeth/Sulfa (SXT) ^e	Nitrofurantoin (FD) ^a	
	Piperacillin/Tazobactam (TZP) ^b	

a - Urines cultures only

b - Additional antibiotic reported for inpatients

c - ≤25mm – set up ESBL confirmatory disks

d – Only extra-intestinal isolates

e – *Acinetobacter baumannii* only

KIRBY-BAUER: GRAM POSITIVE COCCI (Non-CSF SOURCES)

Routine Method		Back-up Method	
BETA-HEMOLYTIC STREP ^{g,h}	VIRIDANS STREP ^{b,g}	ENTEROCOCCUS	STAPHYLOCOCCUS ^{e,h,i}
URINES ^c	URINES	URINES	URINES
Ampicillin (AM)	Ceftriaxone (CRO)	Nitrofurantoin (FD)	Gentamicin (GM)
Levofloxacin (LV)	Penicillin (E TEST)	Penicillin (P)	Levofloxacin (LV)
Ceftriaxone (CRO)		Tetracycline (TE)	Trimeth/Sulfa (SXT)
	NON-URINE SOURCES	Vancomycin (VA)	Cefoxitin (FOX for Oxacillin) ^e
NON-URINE SOURCES	Clindamycin (CC) ^b	Levofloxacin (LV)	Vancomycin (E TEST)
Ampicillin (AM)	Erythromycin (E) ^b		
Clindamycin (CC) ^b	Ceftriaxone (CRO)	NON-URINE SOURCES	NON-URINE SOURCES
Erythromycin (E) ^b	Penicillin (E TEST)	Penicillin (P)	Clindamycin (CC) ^b
Levofloxacin (LV)	Vancomycin (VA)	Vancomycin (VA)	Erythromycin (E) ^b
Ceftriaxone (CRO)		Gent-500 (QUAD PLATE)	Gentamicin (GM)
Vancomycin (VA)		Strep-2000 (QUAD PLATE)	Levofloxacin (LV) ^d
		Linezolid (E TEST) ^a	Trimeth/Sulfa (SXT)
SCREENS ^c			Cefoxitin (FOX for Oxacillin)
Clindamycin (CC)			Vancomycin (E TEST)
Dtest reported if positive ^f			

a - Release if organism is resistant to all antibiotics reported from a sterile site

b - DTEST performed only upon special request

c - When requested by provider, otherwise report isolate comment: &GBS

d - Quinolones are not reported for *Staphylococcus aureus* in blood cultures

e - *Staphylococcus lugdunensis* and *Staphylococcus pseudointermedius* have special reporting criteria. Refer to CLSI Standards (Table 2C)

f - Ceftriaxone and Vancomycin can be released if DTEST is positive

g - *Streptococcus anginosus* group and Group D non-enterococcus are treated as viridans streptococcus

h - Daptomycin Etest is appropriate to add for non-respiratory sources

i - Ceftaroline is appropriate to skin/soft tissue/blood sources

KIRBY BAUER: CSF ONLY

Routine Method		Back-up Method		
VIRIDANS STREP.	PSEUDOMONAS AERUGINOSA	STAPHYLOCOCCUS^b	GRAM NEGATIVE RODS	STREP. PNEUMONIAE
Ceftriaxone (CRO) ^a	Amikacin (AN)	Gentamicin (GM)	Amikacin (AN)	Ceftriaxone (CRO) ^a
Meropenem (MEM)	Aztreonam (ATM)	Vancomycin (E TEST)	Ampicillin (AMP)	Meropenem (MEM)
Penicillin (E TEST)	Cefepime (FEP)	Cefoxitin (FOX for Oxacillin) ^b	Amp/Sulbactam (SAM)	Penicillin (E TEST)
Vancomycin (VA)	Ceftazidime (CAZ)		Cefepime (FEP)	Vancomycin (VA)
	Gentamicin (GM)		Ceftriaxone (CRO)	
	Meropenem (MEM)		Gentamicin (GM)	
	Pip/Tazobactam (TZP)		Meropenem (MEM)	
	Tobramycin (NN)		Pip/Tazobactam (TZP)	
			Tobramycin (NN)	

a - Streptococcus pneumoniae that are resistant to ceftriaxone by disk diffusion need to have an MIC performed. If no MIC is available, the isolate must be sent to the RIDOH

b - Staphylococcus lugdunensis and Staphylococcus pseudointermedius have special reporting criteria. Refer to CLSI Standards (Table 2C).

VITEK PANELS¹

<i>Routine Method</i>		
GN-84	GP-67	GP-67
GRAM NEGATIVE RODS	STAPHYLOCOCCUS ^{h, j}	ENTEROCOCCUS
Amoxicillin/CA	Ciprofloxacin ^c	NON-URINE SOURCES
Ampicillin	Clindamycin ^m	Penicillin
Aztreonam	Erythromycin ^m	Vancomycin
Cefazolin	Gentamicin	Gentamicin High Level
Cefepime	Inducible Clindamycin Resistance	Streptomycin High Level
Ceftriaxone	Levofloxacin ^c	Linezolid ^f
Ciprofloxacin	Moxifloxacin ^{c, d}	Nitrofurantoin ^a
Ertapenem ^k	(Cefoxitin) Oxacillin ^{g, h, i}	
Gentamicin	Tetracycline	URINES
Levofloxacin	Trimethoprim/Sulfamethoxazole	Nitrofurantoin
Meropenem ^k	Vancomycin ^e	Penicillin
Piperacillin/Tazobactam ^{b, n}		Tetracycline
Tetracycline		Vancomycin
Trimethoprim/Sulfamethoxazole		
Nitrofurantoin ^a		

a – Urine cultures only

b – Add TZP disk if Vitek TZP is not performed (Inpatients)

c – Quinolones are not reported for *Staphylococcus aureus* in blood cultures or for MRSA isolates in wound cultures

d – Moxifloxacin only reported for MSSA

e – Vancomycin results of $\geq 2\mu\text{g/ml}$ must be confirmed by GPN3F; $\geq 4\mu\text{g/ml}$ sent to CDC through RIDOH.

f – Release if organism is resistant to all antibiotics reported from a sterile site

g) 1– If Cefoxitin screen is negative and switches oxacillin result of 0.5 $\mu\text{g/ml}$ resistant setup a PBP2a test. If PBP2A is positive report out isolate as MRSA, if result is negative result just as SA. Bring up both cases on rounds.

2. If Cefoxitin screen is negative and oxacillin is 1 $\mu\text{g/ml}$ or 2 $\mu\text{g/ml}$, perform a PBP2a test. If positive, report as MRSA and bring up on rounds. If result in negative, report out as just a SA.

h – When Staph lugdunensis and Staph pseudointermedius results for oxacillin and cefoxitin are discrepant repeat the vitek and perform a KB. Also perform a oxacillin screen for comparison. If the 2nd vitek result matches the KB, report the vitek result. If the vitek and the KB disagree, bring up on rounds.

i – *Staphylococcus lugdunensis* and *Staphylococcus pseudointermedius* have special reporting criteria. Refer to CLSI Standards (Table 2C).

j – Acceptable to release Rifampin for ID provider

k – Ertapenem and Meropenem results that are I or R need to be confirmed with the GNX2F Trek Sensititre panel.

l – Reported antimicrobials vary based on organism identification and product limitations

m – Not routinely reported on urine cultures

n – Resistant Pip/Tazo results from *Enterobacterales* isolated from blood cultures should be confirmed with the GNX2F Sensititre panel

TREK SENSITITRE PANEL

Routine Method		Supplemental Method	
GNX2F^{a,g}	GNX2F^{a,g}	STP6F^b	GPN3F^c
NLF GRAM NEGATIVE RODS	BURKHOLDERIA CEPECIA	STREPTOCOCCUS PNEUMONIAE^f	STAPHYLOCOCCUS AUREUS
Amikacin	Ceftazidime	Ceftriaxone (meningitis)	VANCOMYCIN ^d
Aztreonam	Levofloxacin	Ceftriaxone (non-meningitis)	
Cefepime	Meropenem	Chloramphenicol	
Cefotaxime	Minocycline	Clindamycin ^e	
Ceftazidime	Trimeth/Sulfa	Erythromycin	
Ciprofloxacin		Levofloxacin	
Doxycycline	STENOTROPHOMONAS MALTOPHILIA	Meropenem	
Gentamicin	Ceftazidime (R)	Penicillin(oral)	
Levofloxacin	Levofloxacin ^d	Penicillin(meningitis)	
Meropenem	Minocycline	Penicillin(non-meningitis)	
Minocycline	Trimeth/Sulfa	Vancomycin	
Pip/Tazobactam			
Tobramycin			
Trimeth/Sulfa			

a – Never release tigecycline from this panel

b – Available upon request: amoxicillin/clav, azithromycin, cefepime, cefotaxime, cefuroxime, daptomycin, ertapenem, linezolid, moxifloxacin & trimeth/sulfa

c – For confirmation of possible VISA/VRSA isolates.

d – Suppressed if non-susceptible

e – Not reported if erythromycin is non-susceptible & clindamycin is susceptible

f – Only report: ceftriaxone, meropenem, penicillin & vancomycin for *Strep. pneumoniae* isolates from CSF sources.

g – Reported antimicrobials vary based on organism identification and product limitations

REVISIONS:

4/30/2020 Added Levofloxacin to PA and suppression of Quinolones for MRSA isolates in wound cultures

03/22/2023 Updated RIDOH submission guidelines and confirmatory testing requirements for *Enterobacterales* with Pip/Tazo resistance