

PROCEDURE: ORGANISM ID & AST**I. Aerobic****A. Gram-Positive Cocci - Refer to [Appendix AP33](#) and [Appendix AP34](#)**

1. *Abiotrophia* species/ *Granulicatella* species
 - a. These are the nutritionally variant *Streptococci*. They usually satellite around a *Staphylococcus* streak. Gram stain should be performed on satellite colonies to rule out *Haemophilus* species
 - b. If differentiation of genera is not achievable, report as *Abiotrophia*/ *Granulicatella*
2. *Aerococcus* species
 - a. *A. viridians* is alpha-hemolytic, LAP negative & PYR positive
 - b. *A. urinae* is alpha-hemolytic, LAP positive & PYR negative
 - c. Speciation performed by Vitek MS/ Vitek GP ID card
 - d. Susceptibilities are not routinely performed on non-urine isolates.
 - i. For isolates of *A. urinae* from urine specimens add comment **&AERU**:
"Isolates are typically susceptible to Beta-lactam antibiotics."
3. Alpha-Hemolytic/Non-hemolytic *Streptococcus* species (viridans group):
 - a. Consists of *S. mitis* group, *S. anginosus* group (microaerophilic), *S. mutans* group, *S. salivarius* group, *S. bovis* group:
 - i. Microaerophilic *Streptococci*
 - a) Consists of *Streptococcus anginosus*, *Streptococcus constellatus*, and *Streptococcus intermedius*
 - b) Identified group-wise by their atmospheric requirements
 1. Sub-culture suspected organism onto three BAP plates. Incubate one plate in ambient, another in CO₂, and another in anaerobic conditions. If growth is observed on plates from only CO₂ and anaerobic condition, then this rules in microaerophilic *Streptococci*.
 2. Some strains may grow a little in ambient atmosphere
 - b. LAP positive & PYR negative
 - c. Speciation performed by Vitek MS/ Vitek GP ID card
 - d. Susceptibility is performed by Kirby-Bauer disk diffusion method on Blood Mueller-Hinton
 - i. Note: a penicillin E-test must be tested, since the penicillin disk method is unreliable
 - ii. If unable to grow for susceptibilities add the isolate comment **&UNAB**:
 - a) "Unable to perform susceptibility testing due to the atmospheric growth requirements of this organism."
 4. *Streptococcus pneumoniae*
 - a. Alpha-hemolytic, catalase negative
 - b. Preliminary identification by optochin disk (PTAB) with a zone inhibition of ≥ 14 mm
 - c. Confirmation of identification by Vitek MS/Vitek GP ID
 - d. Susceptibilities performed by MIC
 5. Beta-Hemolytic *Streptococcus* species
 - a. Speciation performed by strep grouping kit, including all sera types.
 - b. Susceptibilities performed by Kirby-Bauer disk diffusion method on Blood Mueller Hinton agar, and incubated in CO₂ at 35°

- c. Susceptibilities are routinely performed on the following isolates:
 - i. All Beta-Hemolytic *Streptococcus* isolated from sterile body sites
 - ii. Beta-Hemolytic *Streptococcus* Group A and B isolated from wounds
 - iii. Group A *Streptococcus* isolated from urine specimens
 - d. Susceptibilities are not routinely performed on Group B *Streptococcus* isolated from urine or genital specimens
 - i. Add isolate comment **&GBS**:
 - “Susceptibility testing not routinely performed. Group B Streptococci are predictably susceptible to ampicillin and penicillin. Call laboratory for further testing if patient is allergic to penicillin.”
 - ii. Susceptibilities performed upon physician request only.
 - e. Microaerophilic, minute Beta-Hemolytic *Streptococci* are identified group-wise by their atmospheric requirements
 - i. Identified group-wise by their atmospheric requirements
 - a) Sub-culture suspected organism onto three BAP plates. Incubate one plate in ambient, another in CO₂, and another in anaerobic conditions. If growth is observed on plates from only CO₂ and anaerobic condition, then this rules in microaerophilic *Streptococci*.
 - b) Some strains may grow a little in ambient atmosphere
 - ii. Speciation is performed using Vitek MS/Vitek GP card
 - iii. Susceptibility testing for microaerophilic, minute Beta-Hemolytic *Streptococci* is performed by Kirby-Bauer disk diffusion method on Blood Mueller-Hinton, including a penicillin E-test
 - iv. Microaerophilic *Streptococci* can group serologically as A, C, F, or G
 - f. All Group F Beta-Hemolytic *Streptococcus* are microaerophilic
6. *Enterococcus* species
- a. Organism is PYR positive & LAP positive
 - b. Speciate all *Enterococcus* isolates reported as pathogens
 - c. Release Linezolid susceptibility results on vancomycin resistant isolates from blood cultures and sterile body sites – Do not release for urines.
 - d. Add High Level Aminoglycoside Resistance statement to non-urine reports **&HLR**:
 - “Synergy (from use of an aminoglycoside plus either Penicillin or Vancomycin) cannot be predicted unless both antibiotics used in combination are susceptible.”
 - e. If sensitivities are performed by Kirby-Bauer method, a synergy quad plate must be performed for non-urines.
 - f. Add nitrofurantoin susceptibility disk to Kirby-Bauer performed on urines
 - g. *E. raffinosus* susceptibilities must be performed by Kirby-Bauer
 - h. *E. casseliflavus* and *E. gallinarum* are intrinsically resistant to Vancomycin and are not considered VREs
7. *Staphylococcus* species/ *Micrococcus* species
- a. *Micrococcus* species
 - i. Bacitracin (Atab) disk sensitive- demonstrating a zone of inhibition of ≥10mm
 - ii. Speciation can be performed by Vitek MS/Vitek GP ID, if necessary
 - b. *Staphylococcus* species
 - i. Sensitivities performed by Vitek GP AST. Alternatively, Kirby-Bauer may be performed.
 - ii. Clindamycin inducible enzyme results may be reported from Vitek directly, otherwise, if a Kirby-Bauer is reported, save plates for 7 days in case a D-test is requested

- iii. If sensitivities are performed by Kirby-Bauer method, set-up an E-Test for vancomycin. Performing vancomycin by disk diffusion is unreliable for *Staphylococcus* spp.
- iv. Coagulase-Negative *Staphylococcus* species
 - a) Latex agglutination can be performed on suspicious *Staphylococcus* species to r/o *S. aureus*
 - b) Speciation can be performed by Vitek MS/Vitek GP ID, if required:
 - 1. *S. lugdunensis*
 - a. PYR positive
 - b. If isolate is reported as a pathogen with susceptibilities, speciation must be done due to differences in oxacillin interpretation between other Coagulase-negative *Staphylococcus* species
 - 2. *S. intermedius*
 - a. May demonstrate delayed latex positivity. If there is a suspicion of *S. intermedius* identify by Vitek MS/Vitek GP ID card
 - 3. *S. saprophyticus*
 - a. Speciation required when isolated from urine
 - b. Susceptibilities not routinely performed
 - i. Add isolate comment **&SAP:**
 "Routine susceptibility testing is not performed. Infections respond to concentrations achieved in urine of antimicrobial agents commonly used to treat acute uncomplicated urinary tract infections."
- v. *Staphylococcus aureus*
 - a) Circular, smooth, light-golden colonies, 1-4mm diameter
 - b) Often Beta-hemolytic (some strains may appear non-hemolytic)
 - c) Speciated by latex agglutination or Vitek MS/Vitek GP ID
 - d) MRSA/MSSA may be determined by performing PBP2a, if requested or there is a discordance with sensitivities. Refer to [Procedure: Alere PBP2a SA Culture Colony Test](#) for further guidance.
 - e) If vancomycin MIC is ≥ 2 by Vitek, an MIC must be done for confirmation – Refer to [Appendix AP41](#)

B. Gram-positive bacilli - Refer to [Appendix AP35](#)

- 1. *Arcanobacterium haemolyticum*
 - a. Beta-hemolytic, irregular rod
 - b. Catalase negative
 - c. Esculin positive
 - d. Identification performed by Vitek MS/Vitek ANC
- 2. *Bacillus* species
 - a. Catalase positive
 - b. *Bacillus anthracis* and *Bacillus cereus* biovar *anthracis* are potential bioterrorism agents and must be ruled out on all isolates
 - i. *Bacillus anthracis* is non-motile, non-hemolytic and has a ground glass colony morphology.
 - ii. *Bacillus cereus* biovar *anthracis* has variable motility, is non-hemolytic and has a ground glass colony morphology
 - c. Potential bioterrorism isolates must be reported to RIDOH Epidemiology and an isolate sent to RIDOH laboratory for confirmation
 - d. If *B. anthracis* and *B. cereus* biovar *anthracis* are ruled out, report generically as *Bacillus* species non-anthraxis

- e. Non-bioterrorism *Bacillus* isolates may be identified by Vitek MS upon request
3. *Corynebacterium* species
 - a. Club-shaped rods
 - b. Catalase positive
 - c. Specification performed by Vitek MS/ Vitek ANC, if appropriate:
 - i. *Ex. Corynebacterium* species isolated in pure culture from multiple blood cultures, catheter tips, or multiple sterile body sites
 - d. *C. ureolyticum* is ruled out in urine specimens if appropriate quantitation present
 - e. Susceptibilities not routinely performed
 4. *Listeria* species
 - a. Coccobacillary rods, singly or in chains
 - b. Catalase positive
 - c. Demonstrates tumbling motility in wet prep, especially at 30°C
 - d. Gray colonies exhibiting weak Beta-hemolysis
 - e. Send to RIDOH for serotyping
 - f. Susceptibilities not routinely performed
 5. *Nocardia* / *Rhodococcus* / Aerotolerant *Actinomyces* / Rapid-Growing AFB
 - a. Send to the Mycobacteriology/ Mycology laboratory for further workup
 - b. *Nocardia* species are PAF positive & AFB negative
 - c. *Mycobacterium* species are PAF positive & AFB positive
 - d. Other species are PAF & AFB variable
 - e. *Rhodococcus* is generally very mucoid and orange
- C. Gram-negative bacilli - Refer to [Appendix AP36](#)
1. *Aeromonas*/ *Plesiomonas*
 - a. Identification from non-stool specimens performed by Vitek MS/Vitek GN ID
 - b. Oxidase positive
 - c. Susceptibilities performed by MIC
 2. *Acinetobacter* species
 - a. Identification can be performed by Vitek MS/Vitek GN card
 - b. For all *Acinetobacter* species perform a Kirby Bauer
 - c. Add amox-clav comment to all *Acinetobacter* species on all benches
&ACIN:
 "Amoxicillin-clavulanate is not active against *Acinetobacter*."
 3. *Achromobacter* species
 - a. If identified by MALDI as *Achromobacter dentrificans*/*xylosoxidans*, send organism to Vitek for GN card for better speciation
 - b. Susceptibility testing performed by MIC
 4. *Burkholderia cepacia*
 - a. Identification performed by Vitek MS/Vitek GN card
 - b. Susceptibility testing performed by MIC
 - c. BCSA plates are inoculated for screening of cystic fibrosis (CF) patient respiratory specimens
 - i. Plates must be incubated ambiently at 35°C for 72 hours
 - ii. *B. cepacia* will grow well and have a yellow to pink-yellow zone on BCSA
 - iii. For other species of *Burkholderia*, refer to miscellaneous GNRs

- d. If *B. cepacia* or member of *B. cepacia* complex is isolated from a CF patient, isolate must be reported as probable *B. cepacia* (*B. cepacia* complex) then sent to reference for identification confirmation
5. *Campylobacter* species
 - a. Curved, thin rods, usually "S" or sea-gull shaped rods
 - b. Oxidase & catalase positive
 - c. Grows best at 42°C in microaerophilic atmosphere
 - d. Susceptibilities not routinely performed
6. *Eikenella corrodans*
 - a. Oxidase positive, usually catalase negative
 - b. Unable to grow on MacConkey agar
 - c. Creates depressions or "pits" in the agar. There is a discoloration of a greenish pigment of the agar.
 - d. Is part of the normal flora of the oral cavity, and therefore can be found in bites, clenched-fist wounds and subacute bacterial endocarditis
 - e. Member of HACEK group
 - f. Identification performed by Vitek MS/Vitek NH card
 - g. Susceptibilities not routinely performed
7. *Enterobacteriaceae* species (*Escherichia*, *Klebsiella*, *Enterobacter*, *Citrobacter*, *Cronobacter*, *Serratia*, *Proteus*, and other *Enterobacteriaceae*)
 - a. Identification performed by Vitek MS/Vitek GN ID
 - b. Susceptibilities performed by Vitek GN AST, or MIC if necessary
 - i. ESBL (Extended-Spectrum-Beta-Lactamases) – [Refer to Appendix AP39](#)
 - ii. CRE (Carbapenem Resistant *Enterobacteriaceae*) – [Refer to Appendix AP40](#)
 - c. *Escherichia coli*
 - i. Spot indole can be performed for rapid identification for flat, dry, lactose fermenting colonies on the urine specimens only
 - ii. If Ceftriaxone ≥ 4/R and spot indole was used for identification, confirmation of the identification must be performed by Vitek MS/Vitek GN card
 - d. *Klebsiella* species
 - i. Identification can be performed by Vitek MS/Vitek GN card
 - ii. Susceptibilities are performed by Vitek GN AST card
 - iii. *K. pneumoniae* and *K. oxytoca* can be differentiated by spot indole
 - a) *K. pneumoniae* is indole positive and *K. oxytoca* is indole negative
 - e. *Pantoea* species
 - i. Identification can be performed by Vitek MS/Vitek GN card
 - ii. Susceptibilities are performed by Vitek GN AST card
 - f. *Serratia marcescens*
 - i. A Pip/tazo disk must be performed on all in-house patients.
 - g. *Salmonella* species
 - i. Speciation performed by Vitek MS/ Vitek GN card
 - ii. Susceptibilities may be done by Kirby-Bauer
 - a) If Levofloxacin is needed, a MIC must be performed for lack of interpretations by Kirby-Bauer method
 - iii. Send isolates to RIDOH for typing.

- h. *Shigella* species
 - i. Identification performed by Vitek GN ID card
 - ii. Vitek MS should not be performed because of its inability to differentiate from clear colonies that could be *E.coli*
 - iii. Susceptibilities performed by Kirby-Bauer
 - iv. Send to RIDOH for confirmation
- 8. *Vibrio* species
 - a. Identification from non-stools specimens performed by Vitek MS/Vitek GN ID
 - b. TCBS media and *Vibrio* species colony morphology:
 - i. *V. cholerae* - large yellow colonies
 - ii. *V. parahaemolyticus* - colonies with blue to green centers
 - iii. *V. alginolyticus* - Large yellow colonies
 - c. Send to RIDOH for confirmation if possible *V. cholera*, *V. parahemolyticus*, or *V. vulnificus*
 - d. Sensitivities are not routinely performed from stool isolates
 - e. Susceptibilities are performed by MIC when recovered from non-stool specimens
- 9. *Yersinia* species
 - a. Identification is performed by Vitek MS/ Vitek GN ID card
 - b. On CIN media *Y. enterocolitica* is characterized as "bull's eye", with a deep-red center, surrounded by transparent border
 - c. Sensitivities are done by Kirby-Bauer method
 - d. Send to RIDOH for confirmation
- 10. *Haemophilus* species
 - a. Identification performed by Vitek MS/ Vitek NH ID
 - b. Member of HACEK group
 - c. Identified by Gram-stain, colony morphology and satellites around *S. aureus*
 - d. A Haem-Quad plate can be performed to rule in certain *Haemophilus* spp. Refer to *Haemophilus ID Agar Procedure*.
 - e. If isolate is *Haemophilus influenzae*, perform B-lactamase and result accordingly using the following isolate comments:
 - i. B-Lactamase Positive- (&HBLP) "*Haemophilus* isolates producing beta lactamase are resistant to Amoxicillin."
 - ii. B-Lactamase Negative- (&HBLN) "*Haemophilus* isolates negative for beta-lactamase are likely to be susceptible to Amoxicillin, Macrolides and Cephalosporin antibiotics."
- 11. *Moraxella* species/ *Neisseria* species
 - a. *N. meningitidis* or *Moraxella catarrhalis* can colonize the nasopharynx and is generally considered as part of normal flora. Work-up only if it is the predominant organism and/ or the Gram-stain demonstrates intracellular Gram-negative diplococci.
 - b. Bring up any *N. meningitides* isolates on ROUNDS
 - c. Susceptibilities are not performed
 - d. Molecular confirmation is required on isolates from patients ≤ 17 years old for *Neisseria gonorrhoeae* and any legal cases
 - e. *Moraxella catarrhalis* is always considered beta-lactamase positive.
- 12. *Pasteurella* species
 - a. Oxidase positive, indole positive
 - b. Unable to grow on MacConkey agar
 - c. Identification performed by Vitek MS/ Vitek GN ID card
 - d. Add isolate comment (&NOSU), "*Susceptibilities not routinely performed.*"

13. *Pseudomonas aeruginosa*
 - a. Identification may be determined if organism is oxidase positive, demonstrates a green pigment, and exhibits as a flat, fuzzy colony
 - b. If mucoid, supplement report with isolate comment &MUC "*Mucoid isolate*"
 - c. Mucoid *P. aeruginosa* can be identified if oxidase positive, green pigment is variable, oxidative metabolism of glucose, and grows at 42°C
 - d. Susceptibilities are performed by Kirby-Bauer with Muller-Hinton at 35°C, in ambient atmosphere
 - e. A MIC can be performed on resistant isolates or for additional antibiotic requests
14. *Stenotrophomonas maltophilia*
 - a. Identification performed by Vitek MS/ Vitek GN ID card
 - b. Susceptibility testing is performed by MIC method
 - i. *S. maltophilia* should be susceptible to SXT. When reading the MIC for SXT, choose the well where the growth is 80-90% inhibited compared to the positive control well
15. Miscellaneous Gram-negative Bacilli/ Non-Fermenting Gram-negative Bacilli - [Refer to Appendix AP37](#)
 - a. Identification performed by Vitek MS/ Vitek GN ID card
 - b. Susceptibilities are performed by the MIC
 - c. *Burkholderia species* (not *B. cepacia*)
 - i. *B. cenocepacia* and *B. multivorans* grow well with pink to red zone on BCSA

II. Anaerobic - [Refer to Appendix AP38](#)

A. Cocci

1. Speciation is done by Vitek MS/ Vitek ANC ID card
 - a. Vitek MS result is acceptable if it is consistent with colony morphology, fluorescence (if applicable), Gram-stain, or any rapid biochemicals.
 - b. If a Vitek ANC ID card is used, perform Gram-stain, fluorescence, and any rapid biochemicals, aerotolerance plate and special potency disk diffusion by using: kanamycin, vancomycin, and colistin for affirmation

B. Bacilli

1. Speciation by Vitek MS/ Vitek ANC ID card
 - a. Vitek MS result is acceptable if it is consistent with colony morphology, fluorescence (if applicable), Gram-stain, or any rapid biochemicals.
 - b. If a Vitek ANC card is used, perform Gram-stain, fluorescence, and any rapid biochemicals, aerotolerance plate and special potency disk diffusion by using: kanamycin, vancomycin, colistin, and bile (for Gram-negative) disks for affirmation
 - c. *Cutibacterium* (formerly *Propionibacterium*) *acnes*
 - i. Gram positive rods
 - ii. Catalase positive
 - iii. Spot indole variable (most are positive)
 - d. *Actinomyces* species
 - i. Gram positive rods
 - ii. Catalase negative
 - a) Note: rare *Actinomyces* species are catalase positive
 - e. *Clostridium* species
 - i. Gram positive bacillus with spores
 - ii. Catalase negative

III. Fungi

A. Yeast

1. Generally considered part of mixed respiratory flora, and should not be reported, unless it is *Cryptococcus neoformans*
2. Urine isolates are speciated if isolated in significant amounts. See urine protocol
 - a. Urine isolates showing “feet” can be reported as Probable *C. albicans*
3. PNA FISH smear must be performed on a patient’s first positive blood culture when Gram-stain shows yeast
4. Isolates from catheter tips are speciated, regardless of quantity. Correlate with blood culture results
5. Sensitivities are performed on select *Candida* spp. by MIC
6. *Candida auris*
 - a. Can be identified by Vitek YST card
 - b. Currently is difficult to identify due to its lack of distinct phenotypic characteristics from other *Candida* species
 - c. Budding yeast, which almost never forms short pseudohyphae and does not form germ tube
 - d. May be misidentified as:
 - i. *C. haemulonii* (Vitek MS)
 - ii. *C. haemulonii* or *C. duobushaemulonii*, (Vitek YST card)
 - iii. Repeat testing using opposite method for confirmation
 - e. Report out *C. auris* if identified
 - f. If *C. auris* is not identified, then it must be sent to RI Department of Health for further identification
 - g. Multidrug resistant
 - h. Isolate must be reported to Infection Control and Department of Health if *C. auris* is either identified or not ruled out.

B. Filamentous fungi

1. Refer to AFB/ Mycology laboratory for workup
 - a. If only one colony grows, do not report. Refer culture to AFB/ Mycology laboratory for its significance.

IV. Attachments

- A. [Appendix AP33 – Catalase Negative, GPC Flowchart](#)
- B. [Appendix AP34 – Catalase Positive, GPC Flowchart](#)
- C. [Appendix AP35 – Aerobic Gram-Positive Rods Flowchart](#)
- D. [Appendix AP36 – Aerobic GNR, Good Growth on BAP](#)
- E. [Appendix AP37 – Aerobic GNR, Poor Growth on BAP](#)
- F. [Appendix AP38 – Anaerobe Identification Flowchart](#)

V. Revisions

- A. 01/29/20 Alterations to testing for organisms to accommodate changes to accessible testing materials or changes to procedures. Addition of *Candida auris*.
- B. 02/03/2025 Updated testing protocols for *Enterobacterales*, *Salmonella/Shigella sp.*, and *Acinetobacter sp.* resulting from Vitek AST-N812 implementation.