

### Carbapenem Resistant Organisms

- I. If an organism is an *Enterobacterales* and resistant to Meropenem (MIC $\geq$ 4) or Ertapenem (MIC $\geq$ 2) by Vitek, confirm result by performing Kirby-Bauer disk diffusion.
- II. Isolates exhibiting Carbapenem resistance are forwarded to RIDOH for genetic testing to determine presence or absence of carbapenem-resistance genes when the following situations occur:
  1. Isolate recovered is:
    - a. *Enterobacterales* **AND** isolate is resistant to at least one of the following:
      - i. Imipenem, Meropenem, Doripenem, Ertapenem by standard susceptibility testing methods (i.e., minimum inhibitory concentrations of  $\geq$ 4  $\mu$ g/mL for Imipenem, Meropenem, Doripenem, or  $\geq$ 2  $\mu$ g/mL for Ertapenem)
    - b. Non-mucoid *Pseudomonas aeruginosa* or *Acinetobacter baumannii* **AND** isolate is resistant to at least one of the following:
      - i. Imipenem, Meropenem, or Doripenem by standard susceptibility testing methods (i.e., minimum inhibitory concentrations of  $\geq$ 8  $\mu$ g/mL)
- III. Reporting:
  1. If the organism is carbapenem resistant add the isolate comment: **&CRE**  
  
**“Carbapenem Resistance Detected. Isolate forwarded to RIDOH for confirmation.”**
  2. Subculture isolate to send to RIDOH for further testing.
  3. Stock isolate in appropriate freezer box, refer to [Isolate Stocking Procedure](#).
  4. The provider/unit must be notified if isolate is a Carbapenem Resistant *Enterobacterales*, refer to [Critical Results Notification](#). (Infection Control will be notified electronically via Theradoc, no phone report is needed)
  5. Culture may be finalized.
- IV. References
  1. January 2015 CDC surveillance definition for CRE  
<http://www.cdc.gov/hai/organisms/cre/definition.html>
  2. RIDOH specimen submissions – Carbapenem Resistant Bacteria  
<http://www.health.ri.gov/programs/laboratory/biological/about/specimensubmission/>