

PROCEDURE: WOUND PROTOCOL**I. PRINCIPLE**

- A. A wide variety of microorganisms that reside on the skin and mucous membranes of the body, as well as those found in the environment can cause skin and soft tissue infections. Superficial wound and abscess specimens usually grow primary pathogens that can cause infections. However, interpretation of these cultures taken from open skin or abscesses may be compromised due to the fact that these lesions are often colonized with many indigenous microbiota. The presence of inflammatory cells and the resultant pus is a hallmark of local infection. Evidence of this process can be documented by the presence of PMNs in the gram stain. The presence of epithelial cells indicates contamination of the specimen with skin or mucous membrane microbiota and may compromise the significance of the culture results. Therefore, the quality of the wound specimen should be determined by the gram stain results.

II. AVAILABILITY

- A. 7 days a week; all shifts

III. SPECIMEN

- A. Refer to [Laboratory Guide – Brown University Health AMC Laboratories](#) for complete acceptance/rejection criteria.
- B. Refer to [Planting Procedure](#) for complete specimen processing information and media selection guidelines.

IV. PROCEDURE

- A. Examine plates at 24-hours and 48-hours incubation.
- B. No or poor growth after initial 24-hour incubation should be re-incubated for evaluation at 48-hours.

V. INTERPRETATION

- A. Any growth of the following organisms should be worked up fully:
1. *Staphylococcus aureus*
 2. *Pseudomonas aeruginosa*
 3. Beta-hemolytic *Streptococcus* Group A or B
- B. Additional organism less frequently encountered that should be worked up fully:
1. *Neisseria gonorrhoeae*
 2. *Listeria monocytogenes*
 3. *Erysipelothrix rhusiopathiae*
 4. *Nocardia sp.*
 5. Rapidly growing *Mycobacteria sp.*
- C. Agents of Bioterrorism should always be ruled out.
- D. Potential pathogens are site dependent and may include *H. influenzae*, Microaerophilic Strep., Non-fermenting GNRs, Enteric pathogens, *B. cereus*, and *Actinomyces*.

- E. Staphylococcus:
1. All *Staphylococcus* isolates morphologically consistent with *S. aureus* must have a Staphaurex or Maldi performed and documented.
 - a. Refer to [Procedure: Organism ID/AST](#) for notes on colony morphology.
 2. Coagulase-negative *Staphylococcus* are fully worked-up with if isolated from multiple cultures.
 3. Pure cultures of Coagulase-negative *Staphylococcus* are reported generically with a consult statement:
 - a. If *S. lugdenensis*, fully work-up with ID/AST.
 4. Only isolates of Coagulase-negative *Staphylococcus species* reported as a potential pathogen need a Maldi or Staphaurex performed and documented in the worksheet.
- F. Yeast:
1. Yeast should be considered commensal flora unless it is isolated in pure culture.
 2. If pure and “feet” present, report as probable *C. albicans*.
 3. If pure and no “feet” present, set up Maldi and report ID. Refer unusual/uncommon identifications to the Mycology lab.
 4. Any isolate identified as *Cryptococcus* species should be reported.
 5. All molds that are not obvious contaminants (growing off the streaked area) should be forwarded to Mycology for evaluation. **Do not report** - The Mycology lab will decide if clinically significant and whether to report.
- G. Commensal Flora: [Refer to Appendix AP19: Wound Culture - Commensal flora](#)
1. Organisms considered commensal flora are dependent on source/site.
 2. Organisms isolated from sources other than those included in the list of commensal flora should be brought up on Rounds.
 3. Organisms reported as commensal flora are to be listed generically in the worksheet.
- H. Pathogen Workup: [Refer to Appendix AP18: Wound Culture Work-up Flowchart](#)
1. Determine if a single morphology or multiple morphologies are isolated.
 2. Single morphology isolated:
 - a. Potential Pathogen: Perform ID/AST
 - b. Commensal Flora:
 - i. Perform ID on isolates from invasively collected specimens and list with consult statement.
 - ii. Isolates from non-invasively collected specimens should be brought up on Rounds to determine significance and extent of workup.
 3. Multiple morphologies isolated:
 - a. Determine the number of potential pathogens:
 - i. **≤2 potential pathogens**
 - a) Perform ID/AST on potential pathogens greater than commensal flora.
 - b) Non-predominant potential pathogens are either listed generically or grouped as commensal flora based on quantitation, gram stain and culture correlation.
 - ii. **>2 potential pathogens**
 - a) Potential pathogens are either reported as mixed aerobic flora or grouped as commensal flora based on quantitation, gram stain and culture correlation.

VI. REPORTING RESULTS

- A. All gram stains, biochemical testing, and/or Maldi identification performed to rule-out suspected pathogens should be documented in the worksheet.
1. Examples of required testing:

Note: This is not a complete list of testing that may be needed to rule out potential pathogens

 - a. MALDI or Strep grouping for Beta-hemolytic *Streptococcus sp.*
 - b. MALDI or Oxidase for NLF Gram-negative rods
 - c. MALDI or Staphaurex for atypical morphologies of *S. aureus*
- B. Quantitate all significant isolates and report with appropriate susceptibility results.
- C. Organisms reported as commensal flora should be listed generically in the worksheet (gram stain and biochemical testing is not required unless performed to rule-out suspected pathogen).
- D. Use the appropriate organism code based on site when reporting commensal flora:

mixcut	Mixed Cutaneous Flora
mixres	Mixed Respiratory Flora
mixora	Mixed Oral Flora
mixent	Mixed Enteric Flora
mixgu	Mixed Genito-urinary Flora
mixvag	Mixed Vaginal Flora
mixaer	Mixed Aerobic Flora

- E. If commensal flora is reported, a designated isolate comment must be included.

1. *S. aureus*, *P. aeruginosa*, or Beta-hemolytic *Streptococcus* Group A or B

ISOLATED:

- a. Add isolate comment: **&NFW – No further workup**

Example:

* Wound Cult and GS	FINAL 01/31/19 10:39
1+ Staphylococcus aureus	
<i>Oxacillin-susceptible staphylococci are susceptible to nafcillin, oxacillin, cephalosporins, carbapenems (meropenem) and B-lactam/B-lactamase inhibitor combinations.</i>	
2+ Mixed Aerobic Flora	
<i>No further workup</i>	

2. *S. aureus*, *P. aeruginosa*, or Beta-hemolytic *Streptococcus* Group A or B

NOT ISOLATED:

- a. Add isolate comment: **&NSPS – No Staphylococcus aureus, Pseudomonas aeruginosa or Beta-hemolytic Streptococcus Group A and B isolated.**

Example:

* Wound Cult and GS	FINAL 01/31/19 10:37
2+ Mixed Aerobic Flora	
<i>No Staphylococcus aureus, Pseudomonas aeruginosa or Beta-hemolytic Streptococcus Group A and B isolated.</i>	

Example:

* Wound Cult and GS	FINAL 01/31/19 10:35
2+ Serratia marcescens	
2+ Mixed Cutaneous Flora	
<i>No Staphylococcus aureus, Pseudomonas aeruginosa or Beta-hemolytic Streptococcus Group A and B isolated.</i>	

VII. PROCEDURE NOTES

- A. Technologist discretion and professional judgement should be used to ensure all necessary testing (ie. Gram stain, Staphaurex, oxidase, Strep Grouping, Maldi) has been performed and documented to either rule out or identify any potential pathogens present.
- B. Hold plates without full ID/AST performed in the 7-day hold racks, not including mixed commensal flora.
- C. [Appendix AP18 - Wound Culture Work-up Flowchart](#) is to be utilized to guide decisions regarding extent of workup for potential pathogens in conjunction with specimen source, site, gram stain, and patient history. It is not all-inclusive and technologist discretion is expected in evaluating potential clinical significance of isolated organisms.
- D. Any workup not consistent with stated guidelines should be brought up on Rounds for further discussion.

VIII. LIMITATIONS

- A. The significance of culture results is dependent on appropriate specimen collection and processing.
- B. Vaginal wounds requested for GC should be tested by Aptima for optimum detection.

IX. REFERENCES

- A. A.L. Leber. 2016. Wound/Abscess and Soft Tissue Cultures p.3.13.1.1-3.13.120. Clinical Microbiology Procedures Handbook, 4th Edition, Vol 1 ASM Press, Washington, D.C.
- B. Mount Sinai Hospital, Department of Microbiology Wounds/Tissues/Aspirates Culture Manual, Version 17.

X. REVISIONS

- A. 8/30/2023 – Procedure updated to clarify required testing and documentation of commensal flora vs potential pathogens present
- B. 3/18/2026 – Updated protocol to reflect changes to CoNS and mixed culture workup