

## TRAINING UPDATE

**Lab Location:** SGMC & WOMC  
**Department:** Core Lab

**Date Distributed:** 1/2/2020  
**Due Date:** 1/31/2020  
**Implementation:** 12/17/2019

### DESCRIPTION OF REVISION

<b>Name of procedure:</b>
<b>Specimen Processing for Microbiology SGAH.M04 v10</b>
<b>Description of change(s):</b>
<p>Header: Changed WAH to WOMC</p> <p>5.3.4 Deleted dialysis and pharmacy cultures; Added sending lab water to Chantilly</p> <p>5.3.3 Only need thio for CSF and synovial fluid</p> <p><b>This revised SOP was implemented December 17, 2019</b></p>

Document your compliance with this training update by taking the quiz in the MTS system.

Non-Technical SOP

<b>Title</b>	<b>Specimen Processing for Microbiology</b>	
<b>Prepared by</b>	Ronald Master	Date: 4/14/2009
<b>Owner</b>	Ronald Master	Date: 4/14/2009

<b>Laboratory Approval</b>		
<b>Print Name and Title</b>	<b>Signature</b>	<b>Date</b>
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Effective Date:		

**TABLE OF CONTENTS**

1. PURPOSE.....	2
2. SCOPE.....	1
3. RESPONSIBILITY.....	1
4. DEFINITIONS.....	1
5. PROCEDURE.....	2
6. RELATED DOCUMENTS .....	4
7. REFERENCES .....	4
8. REVISION HISTORY.....	4
9. ADDENDA AND APPENDICES.....	5

**1. PURPOSE**

To describe the process for microbiology specimen setup, plating and management.

**2. SCOPE**

The scope of this SOP is to ensure the pre analytic processes for microbiology specimens are outlined. These procedures are imperative in determining what pathogenic organisms are present in specimens obtained from patients.

**3. RESPONSIBILITY**

It is the responsibility of all personnel assigned to Microbiology to read, understand and to perform all procedures as described in this SOP.

**4. DEFINITIONS**

**Plating** – inoculation of plated/tubes media with clinical specimen for microbiology culture.

**Inoculation** – to implant microorganisms or infectious material onto a culture medium.

**Streaking** – The use of a loop or other plating tool to inoculate a specimen in order to differentiate microorganisms by color or texture from its surroundings on a culture medium.

## 5. PROCEDURE

### 5.1 Routine Procedure for Plating Cultures:

**All specimens are to be plated in a biosafety cabinet.**

#### 1. Media and its location:

All routine media will be stored in the refrigerator. Microbiology media should be kept in the refrigerator until needed. Media should be allowed to warm to room temperature before use. A working supply (minimum amount) is left at room temperature for use.

#### 2. Loops, Swabs, and Pipettes:

- a. Loops - A wire loop is used for streaking specimens, with the exception of urines. A 0.001 mL calibrated loop must be used to inoculate urine. For sterile urines (cystoscopy, suprapubic aspirate, etc.) use both 0.001 and 0.01 mL calibrated loops.
- b. Swabs - Swabs are used in making the initial inoculation of plates, for preparing smears, and for inoculating specimens into broth media. If a specimen is submitted on a swab it must be submitted in a culturette containing holding medium to prevent drying out. A swab is convenient for inoculating certain specimens onto media, e.g., stool, sputum. Sterile swabs are available at the plating bench.
- c. Pipettes - A sterile pipette may be used to inoculate liquid specimens into broth media, such as thioglycolate, and any plated media. A sterile pipette should be used to inoculate CSF and other body fluids and environmental cultures of liquids. To inoculate thioglycolate with a pipette, insert pipette to bottom of tube and slowly evacuate sample as you withdraw the pipette.

#### 3. Preparation of Smears and Gram Stains:

- a. Smears - Write the accession number, specimen source, date and the patient's last name on the slide. Using a sterile loop or swab, make a smear about the size of a nickel near the center of the slide. Let the slide air dry, then heat fix. ALWAYS MAKE SMEAR AFTER INOCULATING MEDIA TO AVOID CONTAMINATING THE SPECIMEN.

#### 4. **Inoculation and Streaking of Media:**

- a. The first process in the cultural examination of clinical specimens is the selection of appropriate isolation media. The plating table AG.F191 lists the media suitable for the isolation of microorganisms most commonly recovered from various clinical specimens. It is desirable to inoculate more than one kind of isolation medium unless a single organism is targeted.
- b. The purpose of isolation is to obtain bacterial colonies representing progeny of a single cell and thus provide the source of a pure culture. The streaking of materials onto the surface of the medium provides such results. The streaking method must be such that, (a) part of the medium is inoculated with a large amount of material, and (b) subsequent streaking will allow for growth of isolated colonies. When done properly, the completed streaking should cover essentially the entire surface of the medium. When streaking plates, flame the loop between the first and second streak area to avoid overly heavy growth or use a disposable loop. An illustration and explanation of an acceptable method is provided in Addendum A, Figure 1. Addendum A, Figure 2 illustrates the proper streaking procedure for a urine colony count.

#### 5.2 **Incubation of Plates:**

Plates and broth media are incubated in a CO<sub>2</sub> incubator at 35 ±2°C except for chromogenic medium for MRSA which is incubated at 35-37°C in air (non-CO<sub>2</sub>).

#### 5.3 **Specimens and Special Requirements:**

##### 5.3.1 **Stool Cultures**

1. If *E. coli* O157 is ordered, use test code XECOL.
2. Stool in transport media is to be sent to Chantilly. Stool specimens will be plated in Chantilly.

##### 5.3.2 **IV Catheter Tips**

1. Perform all steps in a biological safety cabinet.
2. Using sterile forceps, remove catheter tip from transport tube.
3. Lay the catheter tip on a blood agar plate, and using sterile forceps, roll tip 4-5 times over entire plate. If the catheter tip is longer than 2 inches (5 cm), use sterile scissors or scalpel to cut the end closest to the top of the tube (proximal end) prior to rolling the distal end on the plate. The proximal end may be rolled on a second plate, if desired.
4. Leave the catheter tip on the plate, do not press it into the agar.

##### 5.3.3 **Cerebrospinal Fluid and Other Sterile Body Fluids**

1. Sterile body fluids must be processed immediately upon receipt.
2. Include a chocolate agar plate for all sterile body fluids submitted for culture.
3. Include thioglycolate broth for CSF and synovial/joint fluids

### 5.3.4 Environmental Specimens

1. Environmental samples are specimens not obtained from patients.
2. These specimens are primarily ordered by Infection Control and Prevention and the Laboratory. The most common samples received are listed below. For other requests contact a supervisor or Chantilly microbiology.
3. Laboratory
  - a. Water used for laboratory reagents
  - b. Accession, label and transport to Chantilly
  - c. Refrigerate the water and send to Chantilly to arrive on day shift.
4. Infection Control and Prevention
  - a. Most samples will be swabs from surfaces.
  - b. Order as XENVR and add the appropriate specimen description
  - c. Inoculate the swab onto the first quadrant of a Sheep blood agar plate (BAP)
  - d. Use a sterile loop to streak in 4 quadrants to obtain isolated colonies
  - e. For unusual requests, contact the Microbiology Director

### 6. RELATED DOCUMENTS

Plating Chart for Media by Source and Test Code (AG.F191)  
 Sources for Anaerobic Culture Table (AG.F326)

### 7. REFERENCES

N/A

### 8. REVISION HISTORY


Version	Date	Reason for Revision	Revised By	Approved By
		Supersedes SOP M006.007		
000	5/26/09	Addenda D: media change for MRSA screen	L. Barrett	R. Master
001	3/10/10	Section 5: Change stool culture, add IV cath tip	R. Master	R. Master
002	5/17/10	Section 5.1: Delete requirement for date on opened media 5.2: Change temperature to 35 ± 2°C	R. Master	R. Master
003	7/12/11	5.3.2 Specified catheter length	R. Master	R. Master
004	5/21/12	5.3.1 Deleted plated media for stool cultures	R. Master	R. Master
004	5/21/12	Figure 4: Updated stool cultures	R. Master	R. Master
005	4/9/13	5.2 Add exception for MRSA chromogenic medium 5.3.3 Add centrifugation of sterile body fluids 5.3.4 Add environmental cultures	R. Master	R. Master
006	4/27/15	5.3.3 Edited centrifugation of sterile body fluids Section 6: Moved plating chart from section 9, added Anaerobic culture table Footer: Version # leading zero's dropped due to new EDCS in use as of 10/7/13	R. Master	R. Master
7	3/21/17	Header: Add WAH 5.1.4.a: Changed location of plating chart. Clarified use of multiple media.	R. Master	R. Master

<b>Version</b>	<b>Date</b>	<b>Reason for Revision</b>	<b>Revised By</b>	<b>Approved By</b>
8	3/14/18	Section 9: Added Figure 4: Z streak technique	R. Master	R. Master
9	12/4/19	Header: Changed WAH to WOMC 5.3.4 Deleted dialysis and pharmacy cultures; Added sending lab water to Chantilly 5.3.3 only need thio for CSF and synovial fluid	R. Master	R. Master


**9. ADDENDA**

- A. Figure 1 – Acceptable Method of Plate Streaking
- B. Figure 2 – Proper Streaking for a Urine Culture and Colony Count
- C. Figure 3 – Proper Streaking for a Biplate
- D. Figure 4 – Z-Streak Technique


Figure 1 **Acceptable Method of Plate Streaking**



**Step 1:** Using a loop or a swab inoculate the specimen onto one edge of the plate, covering about one third of the plate.




**Step 2:** Flame loop and cool it by stabbing into the sterile agar. With cooled loop, streak at a right angle to the initial inoculum going back and forth




**Step 3:** Flame loop if specimen is likely to contain a lot of normal flora. Cool the loop. Rotate plate again, and entering only the isolation area, draw loop over the previously uninoculated portion of the plate. (Be careful to not streak over the initial inoculum.)

Figure 2 **Proper Streaking for a Urine Culture and Colony Count**




**Step 1:** Swirl urine to mix. Select sterile calibrated loop. Dip into the bottom of the urine sample and streak down the middle of the plate.



**Step 2:** Starting at the top, go back and forth numerous times over the initial streak line to facilitate the isolation of bacterial colonies.

Figure 3 **Proper Streaking of Bi-plates (non-urine specimens)**




Biplates are inoculated by initially streaking about 15-20% of the plate, then flaming before performing the downward streak.



Figure 4 Proper Z-streak Technique

Step 1  
Roll the swab over the plate in a “Z” pattern



Step2  
Using a sterile loop, streak back and forth numerous times over the initial inoculum

