



Current Status: <i>Pending</i>		PolicyStat ID: 8677574
Valley Laboratories	Origination:	N/A
	Effective:	<i>Upon Approval</i>
	Final Approved:	N/A
	Last Revised:	N/A
	Next Review:	<i>2 years after approval</i>
	Owner:	<i>Lindsey Westerbeck: Dir, Lab</i>
	Policy Area:	<i>Lab - Microbiology</i>
	References:	
	Applicability:	<i>Valley Laboratories</i>

Processing CSF Specimens for Culture in the Hospital Lab

PURPOSE

This procedure describes how to process CSF specimens for culture.

POLICY

- All work on CSF specimens is to be done in the biological safety cabinet.
- Order Code: **CSFC** (Culture, CSF)
 - Includes routine bacterial culture and STAT Gram stain
- Order Code: **GS** (Gram Stain)
 - STAT Gram stain is to be manually ordered upon receipt in laboratory when a CSF culture is requested, regardless of culture order priority.
- STAT Gram stain is to be setup after the culture has been inoculated to avoid contamination.
- When there is insufficient specimen, the Gram stain is to be omitted rather than the culture.
- If aliquoting of specimen is performed, no aliquots are to ever be returned to the original container to avoid potential cross-contamination.
- CSF culture and STAT Gram stain should be setup ideally within 1 hour of collection.

EQUIPMENT, REAGENTS AND SUPPLIES

- Biological Safety Cabinet
- Centrifuge
- Vortex
- Chocolate Agar Plate
 - Store at 2 to 8° C, protect from direct sunlight and UV light. May use until expiration date when kept refrigerated.
 - Alternatively, store at room temperature (no more than 25° C) for up to 24 hours. Plate **must** be discarded if not used within 24 hours.
- CO₂ Bio-Bags - BD Bio-Bag Type C
- Sterile Pipette
- Disposable Inoculating Loop

SPECIMEN REQUIREMENTS

- CSF is the only specimen type that will be plated on-site prior to sending to Sutter Health Shared Lab

(SHSL) Microbiology for culture work-up.

- All other culture specimen types will go direct to the SHSL Microbiology.
- To minimize contamination with normal skin flora, **culture** and **STAT gram stain** should be performed on **CSF tube 2 or 3** (unless otherwise specified in the order).

PROCEDURE

Organizing Specimens for Culture

Step	Action
1.	Print Sunquest labels for the specimen being set-up for culture.
2.	Obtain chocolate agar plate and make sure media is at room temperature prior to specimen inoculation.
3.	Place a barcode label on each specimen tube and plate.
4.	Prior to culture setup, it is important to check the following: <ul style="list-style-type: none"> • Patient label on specimen matches the Sunquest labels • Specimen received is acceptable for the test requested • Specimen source is appropriate for on-site culture setup
5.	IMPORTANT: Label plate with tech code and date/time of culture setup.
6.	Using a pencil, label a slide for STAT gram stain.

Processing Specimens for Culture

NOTE: Do not use a pipette to mix the sediment because the bacteria and cells may adhere to the sides of the tube resulting in false negative results.

Specimen Type	Process	
CSF	If Specimen is thin and volume is >1cc	Then <ul style="list-style-type: none"> • Centrifuge CSF specimen for 15 minutes at 1000 x g. • Decant supernatant into a labeled, sterile tube, leaving approximately 1ml at bottom of original tube. • Vortex the sediment for about 30 seconds to resuspend the pellet. • Using a sterile pipette place 1-2 drops of re-suspended sediment onto plate. • Make smear for STAT gram stain - refer to Valley Laboratories procedure, <i>Preparing a Gram Stain in the Hospital Laboratory</i>.
	Specimen is thin and volume is ≤1cc or Specimen is viscous or cloudy (<i>regardless of</i>	<ul style="list-style-type: none"> • Do not centrifuge. • Using a sterile pipette, place 1-2 drops of well-mixed CSF specimen directly onto plate. • Make smear for STAT gram stain - refer to Valley Laboratories procedure, <i>Preparing a Gram Stain in the</i>


volume)

Hospital Laboratory.

Streaking Plates for Isolation

Step	Action
1.	Inoculate plate by placing 1-2 drops of well mixed specimen to <u>one</u> quadrant with sterile pipette.
2.	Streak plate for isolation.
3.	Using the disposable loop spread the initial inoculum back and forth over a fourth of the plate, or the 1 st quadrant. <ul style="list-style-type: none">• NOTE: Avoid touching the sides of the Petri dish.
4.	Turn the plate 90°. Pass the loop through the edge of the 1 st quadrant approximately 2-3 times, while streaking into the 2 nd quadrant.
5.	Flip the loop and continue streaking in the 2 nd quadrant without going back to the 1 st quadrant.
6.	Rotate the plate another 90° and repeat the above procedure until two additional quadrants are streaked, for a total of 4 quadrants. <ul style="list-style-type: none">• NOTE: When streaking the additional quadrants, do not to go back into the 1st quadrant.

IMAGE: Streaking for Isolation - Total of 4 Quadrants



Preparing Samples for Transport

Step	Action
1.	Tape or parafilm the diameter of the plate to prevent any opening during transport.
2.	Place inoculated media into CO ₂ Bio-Bag. Seal and crush ampule to generate the CO ₂ environment.
3.	Place CO ₂ bag containing plates into designated pink transport bag for SHSL Microbiology.
4.	Transport at room temperature to SHSL Microbiology.
5.	Retain remaining CSF specimen on-site, clearly documenting on label the type of sample (i.e. supernatant).

RELATED DOCUMENTS

- Valley Laboratories procedure, *Preparing a Gram Stain in the Hospital Laboratory*

REFERENCES

- Centers for Disease Control and Prevention, *Laboratory Methods for the Diagnosis of Meningitis*, Chapter 6, rev April 15, 2016.
- *Clinical Microbiology Procedures Handbook*, 4th ed., Leber et al, 2016.

All revision dates:

Attachments

No Attachments

Approval Signatures

Step Description	Approver	Date
Lab Medical Directors	Kristen Vandewalker: MD	pending
Lab Medical Directors	Andrea Ong: MD	pending
Lab Medical Directors	Hannah Wong: MD	pending
Lab Medical Directors	Rowberry Ron: MD	pending
Lab Medical Directors	Jamie Cassity: MD	pending
Lab Medical Directors	Mary Keohane: MD	pending
Lab Medical Directors	Marian Butcher: MD	pending
	Lindsey Westerbeck: Dir, Lab	10/2/2020