

Aliquoting of Specimens

- Knows acceptable volumes for aliquoted urines
 - Urine CT/GC volume acceptable only between Black lines
 - Urine Culture acceptable only if filled to Black line on tube, **Do Not** open tube to add.
 - Urine Analysis (UA) – minimum 3mL or fwd to med-center lab.
- Knows how to review specimen requirements for pour-offs.
 - Reference LabNet for proper handling and transport
 - Review locations for all sendouts, evaluate for proper volumes to satisfy orders
 - Labels tubes, Aliquots appropriate volumes, adheres to appropriate transport temperature
 - Tests requiring Special Handling