

15.000 的复数电影 15.000 的复数自己的现在分词 15.000 的现在分词 15.000 的复数形式 15.000 的复数形式 15.000 的复数
DOCUMENT NUMBER: SCPMG-PPP-0123
DOCUMENT TITLE:  Tourniquet Application
DOCUMENT NOTES:
的。这种是一种,我们就是一种的,我们就是一个人的,我们就是一个人的,我们就是一个人的。""我们就是一个人的,我们就是一个人的,我们就是一个人的。""我们就是一

LOCATION: SCPMG-rel	VERSION: <sub>02</sub>	
DOC TYPE: Preanalytic	STATUS: Release	

NEXT REVIEW DATE: 23 Mar 2019	
EXPIRATION DATE:	

AUTHOR: K053260	PREVIOUS NUMBER:
OWNER: K053260	CHANGE NUMBER: SCPMG-CR-0309

Kaiser Permanente Medical Care Program California Division – South SCPMG Laboratory Systems Preanalytical Processing Policy

## **Tourniquet Application**

#### Introduction

The tourniquet is applied above the intended venipuncture site to restrict venous blood flow and make the veins more prominent.

#### Scope

This procedure is intended for any phlebotomist or health care worker who is certified to perform venipuncture or is within their scope of practice.

#### **Policy**

- The tourniquet will be placed 3 to 4 inches above the intended venipuncture site
- When used for drawing from a hand vein, the tourniquet will be placed proximal to the wrist bone.





- The tourniquet will be tight enough to restrict venous blood flow but not too tight which will restrict arterial flow.
- Never apply a tourniquet so tight that it pinches, hurts or causes the arm to turn red or purple.
- Never apply a tourniquet over an open sore.
- Do not allow the patient to pump their fist after a tourniquet is applied.
- Always use a new tourniquet with each patient draw.
- Never leave a tourniquet in place for longer than 1 [one] minute.
- The tourniquet will be inspected for rips or tears prior to use and disposed of in the trash if compromised.

# Non-Controlled Documents

McCall and Tankersley, Phlebotomy Essentials, 5<sup>th</sup> Edition, Lipponcott, Williams and Wilkins, 2012. Chapter 8.

#### Authors

- Rebecca Rosser
- Preanalytical Processing Work Group

### Signature Manifest

**Document Number: SCPMG-PPP-0123** 

Revision: 02

Title: Tourniquet Application

All dates and times are in Pacific Standard Time.

## New Lab Director - Preanalytic

### **Initial Approval**

Name/Signature	Title	Date	Meaning/Reason	
Janice Wolf (K119893)	Operations Director	17 Mar 2017, 08:37:13 AM	Approved	

## **Final Approval**

Name/Signature	Title	Date	Meaning/Reason
David Quam (P092597)	Rgnl Mg Admn-Pmg Executive	23 Mar 2017, 10:13:44 AM	Approved

#### **Set Effective Date**

Name/Signature	Title	Date	Meaning/Reason
Rebecca Rosser (K053260)	RRL ED CONSULTANT		
Matthew Jones (F754627)	Systems Consultant	23 Mar 2017, 12:36:19 PM	Approved

## **Tourniquet Application**

Reviewed and approved by (for Medical Center Area Approval Only):

SIGNATURE	DATE
Jeinny	5-10.18
Name: Janice M. Wolf	
Operations Director, Area Laboratory	
Mys	5/10/18
Name: Sony Wirio, MD	
CLIA Laboratory Director	