

SBAR COMMUNICATION

Requirements for COAG Specimens for Regional Lab.

SITUATION: The Coag specimen are being aliquot into an amber tube.

BACKGROUND: Once the Coag specimen that have been aliquot into an amber, these specimens are tracked and transfer to Sherman Way Regional Lab for testing.

ASSESSMENT: When specimen are received at Regional Lab, the testing site is unable to see the quantity of the specimen through the amber tube.

RECOMMENDATION: ONLY FOR COAG SPECIMENS All Coag specimens that must be tested at Regional Lab will now be aliquoted into a clear tube with red top. Please see LabNet for processing requirements.

Effective: November 09,2021