

ACCESSIONING IN KRMS CO FUNCTION

1. ENTER MR, ensure appropriate demographics populates.
2. Perform 2 tabs at the *ORDERING PROVIDER

```
COLLECTION LABEL OUTPATIENT ACCESSIONING SCREEN          MON - 03/27/23  13:43
                   *ACCESSIONING LOC: 020
MR# 00  0011544620
  -OR-NAME KPHC, TENMLZZTEST          &SEX 2      *DOB 0101900  AGE 123
    AKA:
SORT GRP:                *PRIORITY R0                *CALL N (Y/N)
*ORDERING LOCATION 020    *ORDERING MD: HSHFI  SHAFI, HEDYEH
SEND RPT TO:                PHYS LOCTN:
PHONE (ALT)                (DAY) 152  9991212 (EXT)    (EVE) 000  0000000
LAB CMT
-----
Please review this data before continuing.                SYSGEN(Y/N) Y
CMD:  |      RE=REINQ/F5      PV=PREV SCN/F6      RR=RETURN/F2      MM=MASTER/F3
      |      CA=COMP AC/F9      C1=COMP(W/LABEL)/F10      C2=COMP(WO/LABEL)/F11
```

3. Select Alpha Search (AS)

Reflection Workspace - [TPXPROD.rd3x]

TPXPROD.rd3x

COMMON ORDER SCREEN

PATIENT: KPHC, TENMLZZTEST MR# 00 001154-46-20 WED - 04/28/21 11:30
PAGE 1 OF 7

? HEMOGLOBIN A1C	? BILI TOTAL	? G-6 PD
? CBC W/O DIFF	? BILI DIRECT	? IRON/TIBC
? RETIC COUNT	? BUN	? LDH
? CBC W/DIFF	? CALCIUM	? LEAD
? SEMEN ANALYSIS	? CHOLESTEROL	? LIPID PROFILE
? T CELL PANEL	? CHOL/HDL RATIO	? MAGNESIUM
? SED RATE/ESR	? CPK ELECTROPHORESI	? PHOSPHORUS
? PREGNANCY/URINE	? CPK-TOTAL	? POTASSIUM
? PREGNANCY/SERUM	? CREATININE	? TOTAL PROTEIN-SERU
? URINALYSIS	? ELECTROLYTES	? TRIGLYCERIDES
? TOTAL HCG(INTACT/B	? FRUCTOSAMINE	? URIC ACID
CHEMISTRY	? GAMMA GT	**ELECTROPHORESIS*
? ALK. PHOS.	? GLUCOSE/FASTING	? IMMUNOGLOBULIN
? ALT/SGPT	? GLUCOSE/RANDOM	? HEMOGLOBIN
? AST/SGOT	? GLUCOSE TOLERANCE	? PROTEIN/SERUM
? AMYLASE	? GLUCOSE-2 HR PP	? PROTEIN/URINE

CMD/PAGE#(01-18): _ | PR=PROCESS/F5 HS=HELP SCREEN/F1
01 02 03 04 05 06 07 08 09
10 11 12 13 14 15 16 17 18
PF=PAGE FWD RR=RETURN/F2 AS=ALPHA SEARCH/F6
SRSEOMST PANGAN, AUGUST MM=MASTER/F3

TI 0 20,20 A

Connected to the host [tn3270.mvsprod.jp.org] CAP NUM SCRL 11:30 AM


4. Look for and select the appropriate test.

```
PROCEDURE SELECTION                                MON - 03/27/23  13:52
UNIT/RM-BD                                         PT NAME KPHC, TENMLZZTEST  MRN 00 001154-46-20
--PROBE THE DESIRED PROCEDURE--
-----
01=CA 125-SERUM                                   17=CARBAMAZEPINE
02=CA 15-3                                         18=CARBON MONOXIDE-WBLOOD
03=CA 19-9                                         19=CARCINOEMBRYONIC ANTIGEN
04=CADMIUM-U24HR                                  20=CARDIO CRP, HIGH SENS.
05=CADMIUM-WHOLE BLOOD                            21=CARDIOLIPIN ANTIBODY
06=CALCITONIN-SERUM                              22=CARNITINE-PLASMA
07=CALCIUM                                         23=CARNITINE-URINE
08=CALCIUM-IONIZED DIRECT                         24=CAROTENE-SERUM
09=CALCIUM-PRE DIAL                              25=CATECHOLAMINES-PLASMA
10=CALCIUM-URINE RANDOM                          26=CATH. TIPS CULTURE
11=CALCIUM-U24HR                                  27=CBC
12=CALPROTECTIN-STOOL                            28=CBC NO DIFF-PDIAL
13=CAMPYLO PYLORDIS CULTURE                      29=CBC W DIFF (ONCOLOGY)
14=CANA VAN DISEASE DNA                          30=CBC WITH DIFF
15=CANDIDA AURIS W/RFX FUNGAL CUL               31=CBC WITH DIFFERENTIAL
16=CANDIDA SP. ANTIBODIES                       32=CBC WITHOUT DIFFERENTIAL
CMD | -----
PB=PGBWD    PF=PGFWD                                FT=FREE TEXT
```

5. Select Process (PR)



6. Verify the tests



7. Press ENTER when tests are verified. PRINT LABEL.