



OSOM -SEKISUI MONO TESTING

KAISER LAMC (NEW DATE- NEW QC/each day of use)

DATE: _____ KIT LOT# _____ EXP. DATE _____ Kit paralleled(Y/N)? _____

QUALITY CONTROL - EXTERNAL			REFERENCE FOR INTERPRETAION OF RESULTS	
POS CONTROL RESULT Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No	NEG CONTROL RESULT Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No	CLS INITIAL	POSITIVE -The red line is indicating a valid internal control and the blue line is indicating a positive result.	
			NEGATIVE -The red line only is present, indicating a valid internal control.	

QC ACTION LOG AND COMMENTS:

PATIENT TESTING			
PATIENT	RESULT	COMMENTS	CLS INITIAL
	Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No		
	Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No		
	Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No		
	Internal control acceptable? <input type="checkbox"/> Yes <input type="checkbox"/> No		