LAB Dept MEETING – Huddles

Date of Meeting: October 8, 2020

Attendees: Myrna Ocab, Raquel Lecaro, Lourdes Maniago, Mey Phiri, Yetunde Kosoko, Bill Craig, Michelle Trammell, Mark Gomez, Marietes Gonzales, Juanita Fernandez, Vanessa Cardenas, Erica Torres, Patricia Chea, Marissa Calilung, Theda

Bryant

Topic	Details	Action Item, responsible person, date due, or informational only
KUDOS	KUDOS 100% ARUP and Quest logged in as SO. Thank you to Myrna for helping specimen processing area. KUDOS to all Phlebs, 3 rd highest STATs in all 15 facilities. Great work done. Thank you to all CLS and Validators for completing the 2020 competency Service awardees: (give engraved pen and letter) ✓ Raquel Lecaro – 10 years ✓ Theda Bryant- 10 years ✓ Janet Gerges – 5 years Congratulations and thank you for your years of service at KP!	Informational
SAFETY Manual Time Log	 Due to the Statewide Drill being cancelled because of the pandemic, each department is required to complete and validate the Annual Safety / Disaster Training review for all staff. We'll be doing this in the next few weeks, completing the Annual review checklist - work practices, Department resources and Department tour. Annual Environment Health & Safety and Chemical Hygiene Training is also coming up. This includes a KP Learn course training and a quiz. 	
Manual Time Log	 Need to reduce use of manual time log Some suggestion ✓ Let the automated voice finish ✓ Use speaker ✓ Use # ✓ Listen before walking away 	

24 hour Urine	Give all urine Cerner labels to CLS. CLS need the labels to identify if it is a Quest or RRL test.	
ARUP	Send transfer list that has Patient Name and MR number to ARUP. They cancel orders if paperwork does not have two patient identifiers	
Medicopia	 From Bruce Golovkin ✓ Drawing extra tube is non compliant. Only ED is allowed to do pre-emptive draws ✓ Pushing orders early is not a regional workflow. 	
STAT centrifuge	Reminder to use the STAT Spin centrifuge for all ED chemistries and the Lactic acid centrifuge for all lactic acid tests. DO NOT use the 15- minute big centrifuge for ED STATs.	ALL STAFF
STAT TAT	 We're still not meeting our STAT turn around time for ED CBC. The regional goal is 20 minutes from in lab to verify. We need to keep on working on this. Phlebotomists- bring all ED CBC's to testing area as soon as logged in. Do not wait for the Chemistry in the centrifuge. Announce the STAT CBCs. DO NOT batch CBC's, always bring to CLS as soon as logged in. Do not wait to log in the specimens that just got in before taking all at the back. CLSs- analyze promptly as soon as receive in the back. Release the Hgm part right away. Do not wait for the manual diff or smear review if there is. You can check the parameters that can be released and leave the ones that needs to be reviewed like the platelet count, etc. Some suggestions: Phleb to log in using Specimen Log In instead of ORV 	ALL STAFF
	CLS to log in ED specimens (CBCs only? but have the same acc # as some Chem and Coag)	
For non-HC send out orders (Quest, Mayo, Arup)	 Cerner Miscellaneous Outside Lab Order 1. Physician ordered test that is not orderable in HealthConnect 2. Fill out the requisition form 	Phlebotomist

	SCHOOL-TIGAL ANGER PERMANDENT : NON HEALTH CONCEST (POTTING LANGEST) WAS THE CONCEST (POTTING LANGEST) WAS T	
Courier Pathology Form	 Under Lab Staff name and NUID- Who fills up this up- the one who prepared the bin or the lab staff in processing when courier picks this up? Lab staff to sign when courier is here. Need WARM Hand Off. 	Phlebotomist
Phleb Workflow Reminders	 Lab Collect vs RN Collect. Lab Collect- Lab responsibility. RN Collect- RN responsibility. If Lactic acid reflexed which is a lab collect, while patient is still in ED, do we transfer to RN draw? Remember, lab collect orders is our responsibility so no. Call ED to collect and if they don't we have to draw. 	Phlebotomist
	 RN does not want the order anymore. What do you do? See policy MVPH4270. Do not transfer the order to RN collect. 	

If an RN is asking to change the time and/or date of collection, or to cancel an order; the lab assist should instruct the RN to discontinue the lab order in KP Health Connect. Instruct RN to re-order the lab order(s) as appropriate.

In order to update the collection list in the KPPI, the lab assistant must update the handheld (wireless) or dock (in non-wireless). If after 30 minutes the lab order is not discontinued, the lab assistant will proceed to collect the lab order(s).

- 3. You have a hard stick patient for routine; phleb tried 2x; nurse wants next shift to try. What do you do?

 See Venipuncture policy. You still have to follow CPOE guidelines (TC= 15 mins, Routine= 1.5 hours, Stat= Now or <30 mins)
- Draws are limited to 2 [two] attempts per phlebotomy for hard draws. A second person may attempt the draw with a limit of 2 [two] attempts.
 - o Ask patient permission for additional draws after each attempt is missed.
 - If unable to draw the blood, refer patient to a supervisor [outpatient] or nurse [inpatient] per local policy.
 - Inpatient after first phlebotomist misses, inform nurse of inability to draw specimen and ask for instructions.
 - 4. Can phleb modify the time bec it's a hardstick?
 - Follow the rules for modifying orders. If it does not fall on any of those then no. Must follow CPOE guideline.
 - Patient not in room.
 - · Patient receiving blood
 - Patient having procedure
 - Patient refused
 - Unable to draw 1st try
 - Unable to draw 2nd try
 - Can phleb cancel duplicate orders? Only if you are consolidating orders. E.g you have CBC and H&H order, you can consolidate by cancelling H&H. See policy MVPH4270


Phlebotomists are NOT ALLOWED to cancel any lab order(s) in any circumstances, except ONLY when he/she is consolidating laboratory tests into one specimen tube. Use cancellation reason code: "consolidate test".

Consolidation process is used primarily to avoid unnecessary repetition of tests and large blood draw volumes.

NOTE: Consolidate ONLY laboratory tests for the following scenarios:

- 1. Lab tests ordered with the same priority (stat, timed critical, routine) and collection time and date.
- Other scenarios: (see below)

CLS annual competency	 Thank you for submitting your annual competencies. In process of being reviewed and graded. Reminder to complete all assigned quizzes on Med Training. This is part of your annual competency. 	CLS
Urinalysis	 Heads up! UA area will be transferred to Immuno area to give way to additional bio hood. More on this to follow 	ALL STAFF
UBT	Do you have a good idea to improve workflow? Volunteer to be a UBT representative Reach out to Mitch or Suzanne Mitchell L. Marsile 8:41 PM MMarsile@seiu-uhw.org (626)644-2819 Suzanne M Bulf 8:41 PM	
	Susie Bulf 949-929-3017 cell 4 Meeting ended 18m 48s 8:46 PM	



KAISER PERMANENTE HEALTHEONNECT



Effective: September 25, 2020

Clinical Users

New Panel Order for Asymptomatic COVID-19 Testing

Based on recommendation from the SCAL COVID-19 Regional Command Center, we created a new order panel in support of COVID-19 testing expansion to asymptomatic patients.

The ASYMPTOMATIC COVID-19 SCREENING PANEL is under the Panels section during order search.

COVID Asx PNL -SCAL - EffSep2...

Open the attached PDF to review the entire Key Message.

Technical Bulletin

Laboratory Care Delivery System

CYSTATIN-C ADDED TO LAB FORMULARY

Effective immediately, CYSTATIN C WITH ESTIMATED GLOMERULAR FILTRATION RATE will be available for manual orders using the attached form (or <u>click here</u> to access the form on LabNet). Please complete the form and send it to your local laboratory or have the patient bring it with them. This form is also available in LabNet:

- Cystatin C, like creatinine, is filtered by glomeruli, but unlike creatinine, cystatin C is reabsorbed and metabolized by the proximal renal tubules.
- Calculation of eGFR using cystatic-C may be preferable to using creatinine in particular situations including the
 following: extremes of muscle mass (low or high), extremes of protein intake (low or high), concern about the
 validity of race adjustment for eGFR creatinine, and the need for confirmation of borderline low eGFR without
 other risk factors.
- · Ordering is currently restricted to nephrologists and pediatric nephrologists.

TEST INFORMATION		
Quest Test Name [Code]	Cystatin C with Glomerular Filtration Rate, Estimated (eGFR) [94588]	
Test Location	Quest Diagnostics, Nichols Institute, San Juan Capistrano, CA	
Specimen Requirement	1 mL Serum, drawn in GLD6 and separated from cells	
Transport Temperature / Stability	Room Temperature / 28 days	

QUESTIONS?

Client Service Center: 1-888-4LAB NFO, or tie line 8-397-7077



Indicate Med Center Lab:

C Antelope Valley

43112 10° St. W. Lancaster, CA 90534 (951) 726-2229

C Blatchwin Park

1011 Bakhwin Park Blird., Bakhwin Park, CA 91706 (826) 851-5144

C Downey

9333 E. Imperial Highway, Downey, CA 90242 (992) 657-8652

9961 Sierre Avenue, Fonteres, CA 92335 (90% 902-8060

C Kern County

3501 Stockdala Hwy, Bakersfield, CA 93309 (961) 396-6013

4867 W. Burset Blvd, 1th R. Los Angelos, CA 90027 (323) 783-4961

C Ontano

2295 S. Viveyard Ave, Ontario, CA 91761 (909) 724-2970

C Crange County - Anahem

3440 E. La Polma Ave., Anaheim, CA 92806 (714) 544-6110

C Crange County - Irvine

6640 Alton Parkway, Irvina, CA 92618 (949) 932-6100

19662 Carters St. Panorama City, CA 91402 (818) 375-3800

10800 Magnolia Ave., Riverside, CA 92505 (951) 353-4155

C San Drego - San Drego Medical Center

9455 Claremont Mesa Blvd, San Diego, CA 92123 (858) 266-2310

C San Drego - Zon Medical Center

4647 Zion Ave, Sun Diego, CA 92120 (819) 528-5395

C South Bay

25825 B. Vermort, Harbor City, CA 90710 (310) 517-2620

C West Los Angeles

8041 Cadillas Avenue, Los Angeles, CA 90004 (323) 857-2792

C Woodland Hills

5901 DeSoto Ave., Woodland Hills, CA 91365 (818) 719-2625

□ Hegional Heterence Laboratories - N. Hollywood 1668 Sherren Way, N. Hollywood, CA 91605 (818) 503-7077

C Hegional Heterence Laboratories - Chino Hills 1300 Payton Drive, Chino Hills, CA 91709 (818) 503-7077

MR#	DOB:	
Sex: F M		
Ordering Location:		
Ordering Provider:		
Signature:		
NPI / UPIN:		

PROVIDER INSTRUCTIONS:

Complete this fillable pdf form and forward to patient (print or e-mail)

PATIENT INSTRUCTIONS:

NOTE: LAB APPOINTMENTS FOR BLOOD DRAW CAN BE SCHEDULED ONLINE FOR YOUR CONVENIENCE. GO TO KP.ORG OR THE KAISER PERMANENTE MOBILE APP.

- REPORT TO THE LAB WITH THIS FORM. THE LAB WILL PERFORM THE BLOOD
- SPECIAL PATIENT PREPARATION INSTRUCTIONS: NONE

LAB INSTRUCTIONS:

COLLECTED BY:


CYSTATIN C WITH GLOMERULAR FILTRATION RATE, ESTIMATED (EGFR)

ESTIMATED (EGFK)	ESTIMATED (EGFR)		
Test Information			
Alternative Names:	+CysX +CysC		
KRMS Procedure Code:	Order Manuality		
Health Connect Order Code:	Order Manually		
CPT Code:	82610		
Methodology:	Particle Enhanced Turbidimetric Immunoassay (PETIA)		
Performed At:	Quest Diagnostics, Inc. 94588 Via Medical Center		
Formulary Status:	*Formulary		
Special instructions:	Manual Order Form Cystatin C with EGFr		
Specimen	Collection information		
Specimen Container:	GLDS		
Specimen Requirements:	1 mL serum		
Minimum Test Volume:	0.5 mL serum		
Pre-Analytical Handling and Transport			
Transport Temperature:	Ambient		
Specimen Stability Ambient:	26 days		
Specimen Stability Refrigerated:	26 days		
Specimen Stability Frozen:	10 years		

DATE/TIME:

The Riverside Legacy: 30 Years of Pride

Own a piece of our history – commemorative book now for sale!



Kaiser Permanente ONE Riverside is excited and proud to offer for sale our new book: The Riverside Legacy: 30 Years of Pride. This beautiful book tells the history of the Kaiser Permanente Riverside Service Area with lively stories and hundreds of photos starting with the early days planning Riverside Medical Center's construction continuing documenting 30 years of expansion across Riverside County.

Get your copy for \$25 at the volunteer services gift shop.

For more information, contact Luz Soto, Volunteer/Events Supervisor at 951-353-5759 or Luz.M.Soto@kp.org



Technical Bulletin

Laboratory Care Delivery System

MODIFICATIONS TO RAPID ABBOTT ID NOW MOLECULAR COVID-19 TEST

Effective immediately, sample collection instructions for the rapid Abbott ID NOW SARS-CoV-2 assay have been updated to reflect a recent Technical Brief distributed by Abbott entitled "ID NOW COVID-19 Labelling Updates".

The changes can be summarized as below:

- Swab is to be placed into a sterile labeled tube without media, preservatives or anticoagulants. It is then to
 be capped, transported at ambient temperature (15 30°C), and tested within one (1) hour of collection.
 - <u>Do not</u> return swab to original plastic sleeve for transport.
 - Ensure the swab fits securely within the tube and the cap is tightly closed.
 - If greater than a 1-hour delay occurs, properly dispose of the sample and recollect.
 - o Avoid touching the swab to the outside of the original wrapper or transport tube.
- Previous instructions to test within 2 hours, or store refrigerated up to 24 hours is no longer acceptable.

TEST INFORMATION		
KPHC Order Display Name [Code]	COVID-19 (SARS-COV-2), (RESTRICTED) PCR KP LAB RAPID [87635P]	
KRMS Procedure Name [Code]	8763504	
Test Location	Medical Center Laboratories	
CPT Code(s)	87635	
Specimen Requirement	<u>Do not use transport media</u> . Please use the swab provided in the Abbott test kit, provided by your Local Command Center, for optimal test performance. Alternative swabs will not be accepted for testing.	
Transport Temperature	Ambient	
Stability	1 hour	

QUESTIONS?

- Client Service Center: 1-888-4LAB NFO, or tie line 8-397-7077
- Jonathan C. Gullett, MD; Physician Director of Microbiology, jonathan.c.gullett@kp.org
- Ken Van Horn, PhD, D(ABMM); Technical Director of Microbiology, ken.van-horn@kp.org
- Tam Van, PhD, D(ABMM); Technical Director of Microbiology, tam.t.van@kp.org

Technical Bulletins are archived on LABNET for your convenience. http://kpnet.kp.org:81/california/scpmg/labnet/index.htm



1

This concludes the Minutes of the _October 8, 2020 Lab Staff Meeting.

Prepared by: Patricia Chea, Marissa Calilung Date: 10/20/2020_