QUALITY ADVISORY

Stoplight Report Update: Blood Culture Volume Metric

Quality Management & Improvement

For questions regarding this Quality Advisory, please email scpmq-lab-sys-quality@kp.org and reference "QA25-03-02."

Beginning with March 2025 data for Blood Culture Volume, the following updates are in place:

- Target = 10 mL ± 2 mL (8-12 mL) to better align with the manufacturer recommended blood culture collection volume of 10 mL.
- Metric = % of proper fill will replace Average collection volume. The indicator reflects total blood culture collections including laboratory and nursing.

In the monthly Blood Culture Volume report, the collection category by each collector will be classified as below:

Category	Blood Culture Collection Date		
	Prior to 1/1/2024	1/1/2024 to 2/28/2025	3/1/2025 or later
Underfill	< 5 mL	< 8 mL	< 8 mL
Proper fill	5-10 mL	8-10 mL	8-12 mL
Overfill	> 10 mL	> 10 mL	> 12 mL

The Quality Stoplight Report legend will be updated as follows:

	Blood Culture Collection Date			
	Prior to 1/1/2024	1/1/2024 to 2/28/2025	3/1/2025 or later	
Red	Average collection volume: ≤ 3.9 mL	Average collection volume: < 5.0 mL	% Proper Fill: less than 50% of collections	
Yellow	Average collection volume: 4.0 – 4.9 mL	Average collection volume: 5.1 – 7.9 mL	% Proper Fill: between 50% and 80% of collections	
Green	Average collection volume: > 5 mL	Average collection volume: > 8.0 mL	% Proper Fill: 80% or greater of collections	