




1) Procedure:
Log in to WakeOne.



CCL-QRG-007 How to Enter Comments and Edit Results in WakeOne

Encounters Notes Letters Labs Imaging Cardiology Tests Media Surgeries Anesthesia Meds Episodes Other Orders

Filters Text Search Preview Refresh Select All Deselect All Review Selected Side-by-Side Master Report

4 records match filters, all records loaded

Filters: Default filter

Date	CSN	Type	Department	Specialty	Provider
05/18/2015	60000175635	Initial consult	FMPP3A	Fam Med	Higginson, Crystal M
05/05/2015	60000175177	Initial consult	FMPP3A	Fam Med	Higginson, Crystal M
05/04/2015	60000175120	Initial consult	FMPP3A	Fam Med	Higginson, Crystal M
12/09/2014	60000171309	Lab Visit	PHLPP	Lab	

6. Look down and left.
Click **More Activities**

More Activities

Newtest, Test MRN: 61... PCP: None Resident PCP: None Needs Code Status Order

Male, 24 y.o., 01/01/1991 Ref Prov Nam...

Chart Review (Last refresh: 3:52:58 PM)

Encounters Notes Letters Labs Imaging Cardiology Tests Media Surgeries Anesthesia

Filters Text Search Preview Refresh Select All Deselect All Review Selected Side-by-Side Master Report

1 record matches filters, all records loaded

Filters: Default filter


Date	CSN
02/10/2015	60000173102


A dropdown box automatically appears.

7. Select **Enter/Edit Results**

- Encounters
- Enter/Edit Results
- FYI
- Growth Chart
- Open Orders
- Patient Files
- Document List

CCL-QRG-007 How to Enter Comments and Edit Results in WakeOne

8. Select 




Enter/Edit Results


Filters: Outstanding Point of Care Testing ordered on or after 6/13/2015 (0)

No.	Test	Code Type	Order Date	Auth. MD	Order Status	Result Status	Result Date
-----	------	-----------	------------	----------	--------------	---------------	-------------

9. Select Point of Care Testing Docked Device from the Order Type menu.

10. Click All

11. Click 



Enter/Edit Results

Filters: Outstanding Point of Care Testing ordered on or after 3/27/2015 (0)

Filter Orders

Order Type

- OB
- OB Ultrasound
- Ophthalmology
- Order Panel
- OT
- Outpatient Referral
- Pathology and Cytology
- PFT
- PFT Ancillary
- Point of Care Testing
- Point of Care Testing-Docked Device
- PR Charge
- Precaution

Result Status

- Outstanding, including preliminary
- Resulted, including preliminary
- All

All Order Types/Statuses

Search

Date range: 3/27/2015 to

Test:


Order number:

Maximum number of orders: 100

Accept Cancel

A single click shows the test result.

12. Double click test to get edit screen.



Enter/Edit Results

Filters: Point of Care Testing-Docked Device ordered on or after 3/27/2015 (22)

No.	Test	Code Type	Order Date	Auth. MD	Order Status	Result Status	Result Date
600535335	BASIC METABOLIC PANEL, WB POC	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015
600537398	BASIC METABOLIC PANEL, WB POC	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015
600535334	POCT IONIZED CA DEVICE [POC219]	Custom	35/18/2015	HIGGINSON	CF Resulted	Final-Edited	05/18/2015
600535333	POCT GLUCOSE, WB DEVICE [POC	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015
600535332	POCT POTASSIUM, WB DEVICE [PC	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015
600535331	POCT IONIZED CA DEVICE [POC219]	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015
600535330	POCT GLUCOSE, WB DEVICE [POC	Custom	35/18/2015	HIGGINSON	CF Resulted	Final	05/18/2015

CCL-QRG-007 How to Enter Comments and Edit Results in WakeOne

Filter: Point of Care Testing-Do* ordered on or after 3/27/2015 (22)

No.	Test	Code Type	Order Date	Auth. MD	Order Status	Result Status	Result
600535335	BASIC METABOLIC PANEL, WB POC Custom	Custom	35/18/2015	HIGGINSON, CF	Resulted	Final	05/18/
600537398	BASIC METABOLIC PANEL, WB POC Custom	Custom	25/18/2015	HIGGINSON, CF	Resulted	Final	05/18/
600535334	POCT IONIZED CA DEVICE [POC219 Custom	Custom	35/18/2015	HIGGINSON, CF	Resulted	Final	05/18/
600535333	POCT GLUCOSE, WB DEVICE [POC Custom	Custom	35/18/2015	HIGGINSON, CF	Resulted	Final-Edited	05/18/
600535332	POCT POTASSIUM, WB DEVICE [PC Custom	Custom	35/18/2015	HIGGINSON, CF	Resulted	Final	05/18/
600535331	POCT IONIZED CA DEVICE [POC219 Custom	Custom	25/18/2015	HIGGINSON, CF	Resulted	Final	05/18/
600535330	POCT GLUCOSE, WB DEVICE [POC Custom	Custom	35/18/2015	HIGGINSON, CF	Resulted	Final	05/18/

Specimen Type: []
 Collected by: []
 Collection date: 5/18/2015
 Collection time: 03 41 PM

Resulting Lab: ABBOTT [39]
 Technician: []
 Providers: []
 Billing: []
 Resulting: []

Results Message: Recipient: [] Modifier: []

Component	Value	Flags	Low	High	Range	Units	Comment
POC IONIZED CALCIUM [15	1.76	Low	1.00	1.30		mmol/L	
ISTAT - SAMPLE [1582]							CAPILLARY

Result: Date: 5/1
 Time: 03
 Abnormal: Ab
 Status: Fin

Accept/Next
 Accept

The test box will be highlighted and changes can be made.

Comments can be entered in the Comment field. Remember: SunQuest English Text Codes do not translate to WakeOne – don't use them.

Click Accept when all corrections and/or comments have been entered.

To correct a sample misidentification:

1. Replace all values with XXX
2. If a comment is present, replace the comment with XXX also
3. In the COMMENT section, enter this text: **INVALID RESULT, PATIENT SAMPLE MISIDENTIFICATION, NOTIFIED** enter the full name of the provider that was notified with credentials and the date and time of notification **BY** the name of the person that notified the provider

To flag questionable iSTAT results:

1. Replace all questionable results with XXX
2. In the COMMENT section, enter this text: **QUESTIONABLE ISTAT RESULTS, NOTIFIED** enter the full name of the provider that was notified with credentials and the date and time of notification **BY** the name of the person that notified the provider

2) References:

3) Review/Revision/Implementation:

- Review Cycle: All procedures must be reviewed at least every 2 years.
- Office of Record: Department of Pathology, Critical Care Laboratory

4) Guide Created: 1/8/16

5) Previous Revision Date(s):

6) Attachments:

7) Revised/Reviewed Dates and Signatures:

Reviewed	Date: _____	Signature: _____
Reviewed/Revision Date:	_____	Signature: _____
Reviewed/Revision Date:	_____	Signature: _____
Reviewed/Revision Date:	_____	Signature: _____