

PATIENT NAME: Thomas, Justin MRN: 414 5784

FORWARD			WEAK D			REVERSE			Additional Testing, Ag typings			ABO/Rh INTERP	
A	B	D	Saline control	Weak D	Check cells	Check control cells	A1 cells	B cells	A1 lectin			ABO	Rh
-	-	3+					-	3+					

ANTIBODY SCREEN (MEDIA INC TIME)										DIRECT COOMBS TESTING				
CELLS	IS	37C	AHG	Check Cells	Interpretation	TESTS	IS	5'	Check Cells	Interpretation				
1					ABID ?	Polyspecific AHG								
2						IgG								
3						C3d								
Auto						Saline Control								

UNIT CROSSMATCHING							
Donor #	Ag neg for ?	ABO/Rh	IS	37° C	AHG	CC	Interpretation (Circle One)
							Compatible / Incompatible
							Compatible / Incompatible

PANOCCELL -10

Master List

IMMUCOR, INC. Norcross, GA 30071 USA
 US LICENSE NO: 886
 LOT NO: 07141
 EXPIRES: 2020/04/24

NAME: Dick Donald
 NO: 1519113
 INSTITUTION: _____
 BLOOD GROUP: _____
 ANTIBODY IDENTITY: _____
 TECH: LV DATE: 3/21/20

Plasma
 Pel
 30"

VIAL	Special Type	Donor	Rh - Hr		Kell		Duffy	Kidd	Lewis	P	MN	S	Lu ^a Lu ^b Xg ^a	PATIENT'S SERUM TEST RESULTS TEST METHODS
			D	C	E	C								
1		RZR1 A4865	+	+	+	+	+	+	+	+	+	+	+	5TIGCL
2	Bgl(a ⁺)	R1WR1 B9954	+	+	+	+	+	+	+	+	+	+	+	2H12H
3		R2R2 C4460	+	+	+	+	+	+	+	+	+	+	+	3H12H
4	V ₂ , VS ₂	Ror D943	+	+	+	+	+	+	+	+	+	+	+	4H12H
5		Rr E762	+	+	+	+	+	+	+	+	+	+	+	5H12H
6		Rr F919	+	+	+	+	+	+	+	+	+	+	+	6H12H
7		Rr G1815	+	+	+	+	+	+	+	+	+	+	+	7H12H
8		Rr H1554	+	+	+	+	+	+	+	+	+	+	+	8H12H
9	Bgl(a ⁺)	Rr H177	+	+	+	+	+	+	+	+	+	+	+	9H12H
10		Rr G1843	+	+	+	+	+	+	+	+	+	+	+	10H12H
TC	TC (Y1b ⁺)	R1R1 B3715	+	+	+	+	+	+	+	+	+	+	+	TCNH-2H

NOTES:
 An antigen designated with a 'w' represent a weakened expression of the antigen that may or may not react with all examples of the corresponding antibody.
 TC: In most individuals, the amount of antigen present on red blood cells gradually decreases during the first 18 months of age, while the amount of I antigen increases. In rare individuals, this transition does not occur and the phenotype persists on the adult red blood cells. Trace amounts of I antigen may be

PANOSCREEN

Master List

IMMUCOR, INC. Norcross, GA 30071 USA
 US LICENSE NO: 886
 LOT NO: 07139
 EXPIRES: 2020/04/24

VIAL	Donor	Rh - Hr		Kell		Duffy	Kidd	Lewis	P	MN	S	Lu ^a Lu ^b Xg ^a	PATIENT'S SERUM TEST RESULTS TEST METHODS
		D	C	E	C								
I	R1R1 B9171	+	+	+	+	+	+	+	+	+	+	+	37TIGCL
II	R2R2 C4281	+	+	+	+	+	+	+	+	+	+	+	4H12H
III	Rr H1691	+	+	+	+	+	+	+	+	+	+	+	4H12H

Plasma
 Pel
 30"

PANUCCELL -1U Master List

IMMUCOR, INC. Norcross, GA 30071 USA
 US LICENSE NO: 886
 LOT NO: 03096
 EXPIRES: 2020/03/27

NO: Mouse, Mickey
 INSTITUTION: 1234567
 BLOOD GROUP: Plasma
 ANTIBODY IDENTITY: Plasma
 TECH: LN DATE: 3/21/20

30"

VIAL	Special Type	Donor	Rh - Hr			Kell					Lewis			MN			Lutheran		PATIENT'S SERUM TEST METHODS								
			D	C	E	C ⁻	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	Lu ^a	Lu ^b	Xg ^a	Xg ^b	
1	Col(b+)	R1R1 B7367	+	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
2	Col(b+)	R1wR1 B6333	+	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
3		R2R2 C5956	+	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
4		Ror D1863	+	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
5		r'r E170	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
6		r'r F854	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
7		rr G1841	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
8		rr H1506	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
9		rr N2078	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
10		rr N2994	0	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
TC	Va, VS+, He+	Ror D742	+	+	+	+	0	0	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
		Patient's Cells																									

NOTES:
 An antigen designated with a 'w' represent a weakened expression of the antigen that may or may not react with all examples of the corresponding antibody.

TC: He (Henshaw) is a low incidence antigen of the MN system. The antigen is detected on the red cells of approximately 2-3% of African Americans but is not associated with blood samples from Caucasians. The antigen is carried on an unusual form of the sialoglycoprotein that also carries S and s antigens. Ref: Issitt PD. Applied Blood Group Serology. 3rd ed. Montgomery Scientific; 1985.

PANOSCREEN Master List

IMMUCOR, INC. Norcross, GA 30071 USA
 US LICENSE NO: 886
 LOT NO: 03094
 EXPIRES: 2020/03/27

PANOSCREEN

TYPING	Donor	Rh - Hr			Kell					Lewis			MN			Lutheran		Xg									
		D	C	E	C ⁻	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b		P ₁	M	N	S	Lu ^a	Lu ^b			
I	R1R1 B9695	+	+	0	0	+	0	+	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
II	R2R2 C1498	+	+	0	0	+	0	+	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+
III	rr N3983	0	0	+	0	+	0	+	0	+	0	+	0	0	+	0	+	+	+	+	+	+	+	+	+	+	+

Plasma
 PEG
 30"

Plasma
 7:47 am, 3/20/20
 30''

SELECTED PANEL
 WAKE FOREST UNIV BAPTIST MEDICAL CENTER
 Technologist: LW

Patient's Name: Moose Mickey
 Patient's Number: 12345671
 Date: Collected 3/2/20 Date: Tested 3/2/20

Supplier Lot #	Donor/RhHr-Vial	RhHr			Kell			Duf	Kid	Lew	P	MN	Lut	X	Additional Antigens	Patient's Plasma Test Results	
		D	C	E	C	w	f										V
1 IMM PAN 10 07141	G1815 IR #7	0	0	0	+	+	0	+	+	0	+	+	0	+	+	+	57 IgG
2 BIOTEST-1 8009011-00-R1WR1 #3	101943 R1WR1 #3	+	+	0	0	+	+	+	0	+	+	+	0	+	+	+	NH 3+
3 IMM PANOSCR 07139	B9171 R1R1 #1	+	+	0	0	+	+	+	0	+	+	+	0	+	+	+	NH 3+
4 IMM PAN 10 07141	B9954 R1WR1 #2	+	+	0	0	+	+	+	0	+	+	+	0	+	+	+	NH 3+
5 IMM PANOSCR 07139	H1691 IR #3	0	0	0	+	+	0	+	+	0	+	+	0	+	+	+	NH-2+

Student Name: _____ Date: _____

PATIENT NAME: Butler, Danielle MRN: 475 95 50																							
FORWARD			WEAK D			REVERSE			Additional Testing, Ag typings														
A	B	D	Saline control	Weak D cells	Check cells	Control	Check cells	A1 cells	B cells	A1 lectin	ABO Rh												
ANTIBODY SCREEN (MEDIA)						DIRECT COOMBS TESTING																	
CELLS	IS	37C	AHG	Check Cells	Interpretation	TESTS	IS	S'	Check Cells	Interpretation													
1						Polyspecific AHG	—	+															
2						IgG	—	+															
3						C3d	—	+															
Auto						Saline Control	—	+															
UNIT CROSSMATCHING																							
Donor #			Ag neg for ? ABO/Rh			IS			37° C			AHG			CC			Interpretation (Circle One)					
Compatible / Incompatible												Compatible / Incompatible											

PATIENT NAME: Wright, Kelly MRN: 339 74 32																							
FORWARD			WEAK D			REVERSE			Additional Testing, Ag typings														
A	B	D	Saline control	Weak D cells	Check cells	Control	Check cells	A1 cells	B cells	A1 lectin	ABO Rh												
ANTIBODY SCREEN (MEDIA)						DIRECT COOMBS TESTING																	
CELLS	IS	37C	AHG	Check Cells	Interpretation	TESTS	IS	S'	Check Cells	Interpretation													
1						Polyspecific AHG	—	+															
2						IgG	—	+															
3						C3d	—	+															
Auto						Saline Control	—	+															
UNIT CROSSMATCHING																							
Donor #			Ag neg for ? ABO/Rh			IS			37° C			AHG			CC			Interpretation (Circle One)					
Compatible / Incompatible												Compatible / Incompatible											