

**CP RL6s Feb 2023**

2/2/2023	ROM sent to micro, not scanned/received. When nurse inquired, sample found past stability for testing.
2/6/2023	Urine chemistry sent in conical tube. Placed on UA analyzer, unable to use for chemistry.
2/6/2023	Ammonia scanned in tube room, did not make it to chemistry within stability. Location of specimen unknown.
2/8/2023	Unspun gold top placed in white/spun rack on track. Sample tested unspun and found by chemistry. Had to be spun and results corrected.
2/10/2023	CSF samples found on specimen cart, left at room temperature overnight. Glucose and protein canceled as stability exceeded.
2/10/2023	Glucose tolerance test - wrong specimen labeled as 1-hour specimen in CP (fasting specimen was labeled as 1-hour specimen)
2/15/2023	Short CBC not marked as short, put on track. Delayed redraw of stat test.
2/17/2023	Pediatric patient tests not scanned on packing list. One gold top sent to share with two accessions. When error found, sample out of stability.
2/20/2023	Specimen mislabeled. XEHSG specimen came with one patient's req, another patient's tube. Mislabeled in CP.
2/28/2023	CBCD not scanned in tube room. Went to error lane, not brought back to CP, lost.

**CP RL6s March 2023**

3/1/2023	Unspun gold top put on track in centrifuge bypass (white) rack
3/5/2023	Vanc Peak and Vanc Trough added to same sample as add-ons
3/6/2023	ED nurse called lab to inquire about tube types for sendout tests (third shift) - was not helped
3/20/2023	pedi bilirubin specimen "lost" - later found in biohazard bag trash
3/21/2023	urine chemistry in conical tube - went to UA and run on analyzer. Unable to complete urine chemistry testing.

3/23/2023	Ammonia removed from ice between CP and chemistry. Returned back to CP 1 hour after received in lab. Required redraw.
3/23/2023	Pedi tube not scanned in lab. Found after patient already redrawn. Extra unnecessary stick for peds patient.
3/29/2023	Blood specimens for chemistry testing were left uncentrifuged overnight on 3/28/23. OP Phlebotomy notified Client Services and sent tests for new collection. Central Processing received original specimens on 3/29/23 and processed for testing. When error was found, results and charges removed from specimens out of stability.
3/30/2023	Thoracentesis processed under CSF Glucose order