

Beaumont

SoftLab End User Training

QUERY: Order Status Monitor

Topic Overview

- Function of Order Status Monitor
- How to access Order Status Monitor
- Order Status Monitor Screens
- Search Criteria:
 - Patient Status
 - Result Status
 - Specimen Status
 - TAT Monitor

Definitions/Abbreviations


- TAT: Turn Around Time

Function of Order Status Monitor


- A sub-section of Specimen Query
- Phlebotomists use: Orders to be collected
- Techs use: Specimen TAT

How to Access Order Status Monitor

- To access the Order Status Monitor
 - Click on Query on Main Menu
 - Click on Order Status Monitor




Order Status Monitor Screen



- To search for desired orders, patients, or results:
 - Click on **Search by:** at bottom left hand corner
 - Note the four tabs:
 - Patient Status
 - Result Status
 - Specimen Status
 - TAT Monitor

Patient Status

- To access:
 - Click Patient Status
 - Search by MRN
- Top of screen displays the following:
 - MRN
 - Name
 - Admission Date
 - Diagnosis
 - Discharge date (if any)
 - Deceased (if any)
 - Clinic #
 - Room #
 - Bed #



Result Status

#	Patient Name	Order #	Tests
1	DEMOSOFT, ROFIPB...	182500004	NA K CL CO2
			AGAP GLU BUN CREI
			GFRNA GFRAA CA CHOL
			TRIG HDLC LDLC NHDLC
			CHOLR L H I
		181200243	DSDNA IFA
			10.5
		181200189	BUN CREI GFRNA GFRAA
			12 0.10 @GFR1 @GFR1
			BCRRC
			120
		181100191	CO2
			20
		180900284	WBCIR RBC HGB HCT

Patient Status **Result Status** Specimen Status TAT Monitor

▲ - panic result ▲ (yellow) - abnormal result ▲ - absurd result □ (yellow) - turnaround time status

Search by: MRN = 4001652 Search F12 Time Frame: Last month

Result: User: AMARM Env: ROLIVE Term: A0014 NUM 05/26/2011 10:18

- Results display by Order number
- Colored arrows correspond with abnormal results as indicated by the legend in the middle of the screen

Specimen Status

#	MRN	Patient Name	Tests	CLINIC
1	4001652	DEMOSOFT, ROFIP...	<div style="display: flex; justify-content: space-between;"> BMP LIPID L H I </div>	R5SB
			DNA	
			BCRR	
			CO2	
			PDPCH	
			GBMGT	GMYP
			CORT	
			CLCSF	
			NA	
			PHFLD	
			PHFLD	NAKW
			GLU	R5SB
			GLU	
			CREAT	
			CREAT	
			CO2	

Patient Status Result Status Specimen Status TAT Monitor

resulted - collected - received

Search by: MRN = 4001652 Search F12 Time Frame:

- Test display by test ID
 - Color surrounding the test corresponds with status listed at legend in middle of screen
 - TAT status is also displayed

TAT Monitor

#	Clinic Name	Room/Bed	Priority	Patient Name	Received	TAT	Order #	Specimen	Tests
1	5 SOUTH RO B (P...	5317/	Routine	DEMOSOFT, ROF...	05/03/2011 09:17	33197 min	180300086	SYR	PHIFL
2	5 SOUTH RO B (P...	5317/	Routine	DEMOSOFT, ROF...	05/12/2011 12:56	20018 min	181200243	57T	IFA
3	5 SOUTH RO B (P...	5317/	Routine	DEMOSOFT, ROF...	05/25/2011 05:25	1749 min	182500004	57T	NA, K, CL, CO2, ...
4	5 SOUTH RO B (P...	5317/	Routine	DEMOSOFT, ROF...	05/04/2011 21:03	0 min	180400600	NQ1	OLDHX

- Displays Received date/time
- TAT in minutes
- Green symbol = within time frame; Yellow=approaching deadline Red=exceeds TAT parameters for that test

Summary

The following topics have been discussed in this module:

- Function of Order Status Monitor
- How to Access Order Status Monitor
- Order Status Monitor Screen Layout
- Search Criteria tabs

Credits

- Questions on content and/or input or critique on required modifications should be addressed to:
 - Laboratory LIS Soft Support
LaboratoryLISSoftSupport@beaumont.edu