SERVICE UPDATE



CAMPUS AFFECTED: PEKIN, PROCTOR AND METHODIST			
AFFECTS PATIENT CARE			
TEST NAME: TACROLIMUS (FK506)			DATE OF CHANGE: 4/20/17
REPLACING: FK506A			
SPECIMEN TYPE:	WHOLE BLOOD, EDTA LAVENDAR TOP TUBE		
ORDERING CODES:	EPIC:	LAB876	NEW:
	SQ:	FK506	NEW:
OTHER USEFUL	CPT 80197		
INFORMATION:			
	Prep: Pre dose (trough) levels should be drawn or 12 hours post dose.		
	Storage/transport/volume: Refrigerated, whole blood 1 mL (Min 0.25mL)		
	Stability: Room temperature: 7 days, Refrigerated = 7days, Frozen = up to 6		
	months		
	Reported: Daily (Sun – Sat)		
	Normal Range: 5.0 – 15.0 ng/mL		
	Critical Value: >15 ng/mL		

For further Information Refer to UnityPoint Methodist Laboratory User Manual

For Technical Information Contact:

For Scientific and Clinical Information Contact:

Proctor Laboratory 691-1060

Amy Mayes, MT (ASCP) Chemistry Lead 309-691-1060

AmpAMayo 4/13/17

Cindy Schroeder, MT (ASCP) Laboratory Manager 309-683-5151

any Achine 4/13/17

Richard Borge, Director of Laboratories 309-671-8218

Richard J. Burge

Adam Quinn, DO Medical Director 309-691-1060

100 mg mod 4/13/17

Lori Racsa, DO, Clinical Pathologist 309-671-2181

1~10 4/17