



Technology Learning Services

Epic Functionality Job Aid

Inpatient Beaker - Enter or Edit Results for Inpatient Documentation

Below are instructions for using Enter/Edit Results for lab documentation.

IMPORTANT: Always log in to the correct care department in order for automatic charges to drop appropriately.

1. Look at the Brain or Work List to find the Task for the specific lab. Click on the hyperlink for **Enter/Edit Results**.

Ambrosia, Lorelai-LAB

To Do's

- 1116 POC Gastric Occult Blood with Result To be completed in Enter Edit Completed in Enter Edit
- 1116 POC Occult Blood Stool with Result To be completed in Enter Edit Completed in Enter Edit

[Enter/Edit Results](#) [Order Details](#) [Edit](#)

Go to Enter/Edit Results

+ Task + Appt + Patient Event

Task ID	Description	Frequency	Priority	Action	Status
1116	POC Occult Blood Stool with Result To be completed in Enter Edit	One Time	Routine	Enter/Edit Results	Completed in Enter Edit
1116	POC Gastric Occult Blood with Result To be completed in Enter Edit	One Time	Routine		Completed in Enter Edit

2. You can also use the Enter/Edit Results Activity to find the lab. Highlight the lab and double Click to open that Resulting tool.

The process is the same from here as it is for utilizing the hyperlink.

Work List

Summary Chart Re... Historic... Results Work List MAR Flowsheets Intake/O... Notes ...

More Activities

- Allergies/Contraindications
- Annotated Images
- BPA Review
- Cancer Staging
- Care Teams
- Charge Capture(non provider)
- Chart Central
- Code
- Communications
- Cutover
- Delivery Summary
- Demographics
- Device Data
- More
- Document List
- Enter/Edit Results**
- Episodes of Care

Time	Type	Task	Frequency	Priority	Doc
0700		Administer clindamycin (CLEOCIN) IVPB 600 mg Dose 600 mg : 100 mL/hr IntraVENous : Every 8 hours	Q8H	Routine	Doc
0900		Administer furosemide (LASIX) tablet 40 mg Dose 40 mg : Oral : 2 times daily	BID	Routine	Doc
1116		POC Occult Blood Stool with Result To be completed in Enter Edit	One Time	Routine	Completed in Enter Edit
1116		POC Gastric Occult Blood with Result To be completed in Enter Edit	One Time	Routine	Completed in Enter Edit
1120		Reassess Peripheral IV 4/22/2024 Left Antecubital	Q Shift	None	Flowsheets

Enter/Edit Results

Preview Refresh (11:28 AM) Exit Result Edit

Filters Outstanding Orders Default Types Clear Filters

Date/Time	Test	Auth MD	Order Status	Result Status
Today at 11:15	POC Gastric Occult Blood with Result	Rounds, Steve, MD	Active	
Today at 11:15	POC Occult Blood Stool with Result	Rounds, Steve, MD	Active	

Double Click on row for Lab you are Resulting

3. If there is *only one* lab to be resultated, Epic will default to that lab. If there is *more than one*, you will need to choose the lab by highlighting the row of the lab you are documenting.

Enter/Edit Results

Preview Refresh (11:28 AM) Exit Result Edit

Filters Outstanding Orders Default Types Clear Filters

Date/Time	Test	Auth MD	Order Status	Result Status
Today at 11:15	POC Gastric Occult Blood with Result	Rounds, Steve, MD	Active	
Today at 11:15	POC Occult Blood Stool with Result	Rounds, Steve, MD	Active	

Double Click on row for Lab you are Resulting

4. Complete all the necessary information including any Hard Stops. The Status will normally be **Final result**. Click **Accept** once completed.

Enter/Edit Results

Preview Refresh (11:26 AM) Ext Result Edit

Filters Outstanding Orders Default Types Clear Filters

Date/Time	Test	Auth MD	Order Status	Result Status
Today				
Today at 11:15	POC Gastric Occult Blood with Result	Rounds, Steve, MD	Active	
Today at 11:15	POC Occult Blood Stool with Result	Rounds, Steve, MD	Active	

Searched through 1/25/2024 Search Further

Scan Media Manager

POC Occult Blood Stool with Result

Component	Value	Flags	Low	High	Ref Range	Units	Comment
FECAL OCCULT BLOOD SCREEN, ...		⚠			Negative		⚠

Result

Result date: 4/24/2024

Result time: 11:39 AM

Abnormal:

Status:

Accept Accept/Next Discard Changes

5. Go back to the Brain or Work List and click **Completed in Enter Edit**. This documents the Task as **Completed**.

Ambrosia, Lorelai-LAB

To Do's

- 1116 POC Gastric Occult Blood with Result To be completed in Enter Edit Completed in Enter Edit
- 1116 POC Occult Blood Stool with Result To be completed in Enter Edit **Completed in Enter Edit**

[Enter/Edit Results](#) [Order Details](#) [Edit](#)

+ Task [Appt](#) [Patient Event](#)

1100

1116 POC Occult Blood Stool with Result To be completed in Enter Edit One Time Routine [Enter/Edit Results](#) **Completed in Enter Edit**

POC Gastric Occult Blood

Ambrosia, Lorelai-LAB

To Do's

- 1116 POC Gastric Occult Blood with Result To be completed in Enter Edit Completed in Enter Edit

Completed

- 1116 POC Occult Blood Stool with Result To be completed in Enter Edit

+ Task [Appt](#) [Patient Event](#)

6. Documentation can be viewed on the Results activity.

2024 4/10/24 10:20	
MICROBIOLOGY	⊗ ⤴
MYCOBACTERIUM TUBERCULOSI...	
REFERENCE LAB (S...	⊗ ⤴
Cysticercosis Ab IgG	
IMMUNOLOGY (AFTER)	⊗ ⤴
HIV, Antigen and Antibodies	
MICROBIOLOGY	⊗ ⤴
MYCOBACTERIUM TUBERCULOSI...	
Others	⊗ ⤴
Bands Present (Ext)	
Human Immunodeficiency Virus (H...	
Human Immunodeficiency Virus (H...	
Mitogen minus Nil Result (Ext)	
Mycobacterium tuberculosis (Ext)	
Mycobacterium tuberculosis Sour...	
NIL Result (Ext)	
Occult Blood Result (POC)	Positive !

7. If a clinician Enters or Edits a result which subsequently causes an "abnormal value" flag to appear, that flag will stay in place even if you change the results to a value that is no longer abnormal. Only lab personnel have the ability to change or remove the "abnormal value" flag. This is a safety mechanism making it a security point to change flagging of values within Epic Beaker.

APPLICATION:	COPYRIGHT:	PUBLISHED DATE:	WORK AREA:	ROLE:
ClinDoc	© 2024 Epic Systems Corporation. All Rights Reserved.	9/9/2024	Inpatient	Nursing

