

Phlebotomy Notes



JUNE 2015
DAVID KENDRICK, PHLEBOTOMY
SUPERVISOR

JULIE WEST
LABORATORY TECHNICAL SPECIALIST

Review ORDER OF DRAW



- Sterile tubes (Blood cultures, AFB)
- Red Top Tubes (no additives)
- Citrate tube (blue)
- Gold or Speckled (Serum) tubes
- Heparin (green)
- EDTA (purple), PPT (pearl)
- Light Yellow (ACDA acid citrate dextrose)
- Fluoride (gray)
- All other tubes

Blood Culture Collection



- Bionerieux blood culture bottles (aerobic and anaerobic)
- DVD

Quantiferon Draws



- https://www.youtube.com/watch?feature=player_detailpage&v=TOXF6CzPJYA

1.57 – 6.28

- Training video for drawing quantiferon tests (Cellestis – Qiagen)

Collection and Mixing of Quantiferon tubes



- How many times do we mix by inverting? (10 times, up and down)
- Do we draw a purge tube first, if using a butterfly set? (yes)

Review SOPs - Discussion



- SOPs

(General Lab and Phlebotomy)

Troubleshooting



- Hemolysis
- Clots
- QNS
- Diluted by IV fluid
- Mislabeled/Mismatch/Unlabeled

The Joint Commission regulations: NPSG



Two forms identification

Handwashing

Questions?



Thank you