

**Document Name:**

**Un-linking Patients in Soft Bank**

**Approved By:**

Jennifer G. Daley Bernier, A/ Manager, Laboratory Services

**Status: APPROVED**

**PURPOSE:**

This procedure outlines how to un-link patients in Softbank who have had their Medical Record Number (MRN) linked in error.



**POLICY:**

When a patient visits any of the LIS linked Health Centres they are issued a MRN or Client ID for that health centre (Prefixes of: IN, ST, HR, FS). This creates multiple entries in the LIS for the same patient. Each MRN entry can have different patient histories attached to it. Linking the patient MRN's creates a consolidated history and requires the user to review the history of both patient entries.




**\*CAUTION: It is possible to link patients that are not the same person. Ensure spelling of name, Birthdates, and/or PHN match prior to linking patients in SoftBank. If patients are linked that should not be, please un-link patients using the procedure below.**

**PROCEDURE INSTRUCTIONS:**

Follow the steps in the table below to Un-link patients in SoftBank

Step	Action
1	Go to Patient > Maintenance 
2	Go to Link 


**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

<b>3</b>	Enter patients last name (LN) and first name (FN)	
<b>4</b>	Click: F12-Accept or Enter	
<b>5</b>	If required: Select the correct patient from list – highlight and press Enter or double-click on the patient	
<b>6</b>	<b>If:</b>	<b>Then:</b>
	Patient caution window appears	Click: Esc-Quit after reviewing (or F12 if history review is required)
<b>7</b>	Linked patients will show in window 	
<b>8</b>	Using the mouse, unmark the patient(s) to unlink 	
<b>9</b>	F12 - to Accept	
<b>10</b>	Confirmation message will appear:  Click Yes	

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.


11

ULINK exception comment required:



F12 – Accept

12




Click - YES


**EXPECTED RESULTS:**


Patients with history links will be unlinked

- 1** Patients who have a link display an L in the LK field in the patient selection window – the unlinked patient will no longer have an L in the LK field.



MRnumber	Name (last,first m.)	LK	DOB	St
IX00000032	LINKED, ONE	L	1979-01-01	Y
SX00000114	LINKED, ONE	L	1979-01-01	Y
SX00000115	LINKED, TWO		1980-01-01	Y
- 2** The patient caution window will also display that there are linked patients (if more than one patient was linked) – highlight and press Enter to view details of the Link


- 3** Window will list patients linked – patient that was unlinked (Linked,Two) will no longer appear in list



MRnumber	Name (last,first m.)	DOB	MPI/BB
IX00000032	LINKED, ONE	1979-01-01	B
- 4** Click F8-Sel pat to view linked patient (if applicable)

**RELATED DOCUMENTS:**

- TMM20300 Patient History Review
- TMM20400 Linking of Patients in Soft Bank
- TMM20401 Linking patients Job Aid
- TMM20501 Unlinking Patients Job Aid
- IO500 Identifiers – The use of 2 Patient identifiers

**REFERENCES:**

1. SCC Soft Computer. (n.d.). *SoftBank II v25.1 System Design Training Manual*.  
SCC Soft Computer.

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	09-Nov-17	Initial Release	A. Richardson