

## **PURPOSE:**

To allow user to order fractionated products to a patient in the LIS in Softbank module POLICY:

When an order is received for a fractionated product it must be ordered in the LIS to the patient in order to issue it to the patient. This is to ensure a record of transfusion is maintained and accurate for the patient in the laboratory.

## PROCEDURE INSTRUCTIONS:

Follow the steps in the table below to add Fractionated orders to a patient in the LIS

Step	Action			
1	In SoftBank click Patients > Orders  SOFTBANK  1-PATIENT 2-DONOR 3-II  O Orders			
2	Go to Modify			
	■ SOFTBANK - PATIENT Orders			
	1 New_add 1 Print 1 Cancel 1 New_add 1 Print 1 Cancel 1 New 1 Display 1 re⊻iew 1 Actions 1 Results 1 Xmatch 1 Units 1 Transf X Quit			
3	Enter patients last name (LN) and first name (FN) and/or Medical Record Number OR			
	Click: F3-By order and enter the order number			
4	Click: F12-Accept or Enter			
5	Select the correct patient from list if applicable			
6	Select the correct stay from list if applicable			
7	If: Then:			
	Patient caution window appears Click: Esc-Quit after reviewing			
8	Click: F12-Accept			


**Document Name:** 

**Ordering Fractionated Products to a Patient in Soft Bank** 

**Document Number: TMM21600** 

Version No: 1.0 Page: 2 of 3

Effective: 09 November, 2017



# **METHOD PERFORMANCE EXPECTATIONS:**

• LIS will have a record of blood products ordered to patient.

#### **RELATED DOCUMENTS:**

SBM10100 – Order Entry in GUI Lab

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: TMM21600.1OrderingFractionatedProductstoaPatientPRO.docx

Print Date: 09/11/2017 4:01:00 PM

Document Number: TMM21600

Document Name:

Ordering Fractionated Products to a Patient in Soft Bank

Document Number: TMM21600

Version No: 1.0 Page: 3 of 3

Effective: 09 November, 2017

# **REFERENCES:**

SCC Soft Computer. (n.d.). SoftBank II v25 Super User Training Manual.

# **REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	09-Nov-17	Initial Release	A. Richardson
	_		

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: TMM21600.1OrderingFractionatedProductstoaPatientPRO.docx Print Date: 09/11/2017 4:01:00 PM