



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Maintenance Record – BACTEC FX

Distribution:
 Microbiology Instrumentation Manual

Document Number:
 MIC70500

Date Issued:
 06 November, 2017

Maintenance Record – BACTEC FX

Month: _____ Year: _____

Date	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Left Side																															
Temp: 33.5-36.5																															
Internal temp																															
Green LEDs																															
Red LEDs																															
Right Side																															
Temp: 33.5-36.5																															
Internal temp																															
Green LEDs																															
Red LEDs																															
External																															
Yellow lights																															
Red lights																															
Green lights																															
Audible alarm																															
Change CD																															
Current inventory																															
Initials																															

Corrective action:

Change air filters	✓ Performed:	Date:	Initials:
Black filter			
Silver filter			

Reviewed by:

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.