

SCC-LIS_TIP OF THE WEEK

MODULE: SoftBank	SOFTWARE VERSION: 25.1.3.0.16
TOPIC: Hollister Number Not Unique	ISSUE DATE: 17May2018
DISTRIBUTION: Technologists_Bank	PAGE: 1 of 1
ISSUED BY: M-L Dufresne, Territorial LIS Administrator	

TOPIC/QUESTION:

In Order Entry, why am I getting error "Hollister number is not unique"?

Specimen Hollister Number Enter a Ho Specim Order Entry A Hollister number is not unique QK V OK X Cancel

ANSWER/TIP:

"Hollister" or "BBID" number is the unique blood bank bracelet number. If you are getting this error at OE there is one of two possibilities:

- You are mis-entering the BBID on the current patient
 →Double check and try again
- 2) Someone has mis-entered the BBID on a previous patient
 - → Report the incident per established lab protocols
 - →Backtrack which patient is assigned to the BBID number via:

Softbank>5-Management><u>E</u> Report Print>BBspecimens (Double Click)>Date (mmddyyyy-"Enter") > Function S-Search > Type in BBID number and "search next"

EXAMPLE:

Who has BBID# B1111? → It's APPLE, TART!



