


SCC-LIS_TIP OF THE WEEK

MODULE: SoftLab/SoftMic	SOFTWARE VERSION: 4.0.8
TOPIC: Reviewing System Notifications	ISSUE DATE: 08Jun2018
DISTRIBUTION: Technologists	PAGE: 1 of 1
ISSUED BY: M-L Dufresne, Territorial LIS Administrator	

TOPIC/QUESTION: A system notification disappeared before I got to see it!! Can I find it again?



ANSWER/TIP:

All "Received Messages" are saved to your SoftComm until you log out. If you miss one go to SoftComm>View>Messages and you will see the list of everything that has fired since you have been

logged on.

EXAMPLE:

